



| | | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 1 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |
| Prepared by QA Committee | | |
| Issued by: Laboratory Manager | Revision Date: 3/11/2022 | |
| Approved by Laboratory Director: Microbiologist-in-Chief | Next Review Date: 3/11/2024 | |

Uncontrolled When Printed

Introduction 3

Collection and Transport & Storage..... 3

Materials, Equipments and Facilities 3

Specimen Processing 4

Procedure..... 4

 General Precautions..... 4

 Extraction Set up..... 5

 Dispense reagents (regular procedure) 9

 Dispense reagents (massive procedure)..... 17

 Off board preparation..... 19

PCR Set-up..... 22

 Bio-Rad Detection Area:..... 25

 Build Worklist..... 29

 Import Worklist..... 31

 Analysis 35

Interpretation 41

 Quality Control (QC) Samples 41

 Test Samples 42


 General BGI Reporting Process 43

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 2 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

BGI Reporting Process of Repeat Samples..... 45

Reporting 47

1. PRELIMINARY POSITIVE RESULTS (for UHN/MSH/STAT patients only)..... 47

2. FINAL RESULTS 48

2.1 NEGATIVE RESULTS 48

2.2 POSITIVE RESULTS 49

2.3 INVALID RESULTS 50

3. FINAL RESULT for REPEAT of LOW LEVEL DETECTION 51

3.1 CONFIRMED POSITIVE RESULTLS 51

3.2 INDETERMINATE RESULTLS 51

3.3 NSQ to Repeat RESULTLS 52

Cleaning 53

Quality Control..... 54

Troubleshooting..... 57

1. Technical Support 57

2. Common Troubleshooting Tips..... 57

c. Contaminated run..... 60

d. Fanning curves (uneven distribution of Mastermix)..... 61

4. MGI workflow checklist..... 62

5. MGI Reagent Preparation Log..... 64


Related Documents..... 65

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 3 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

| | |
|------------------------------------------------------------------------------------------------------------|----|
| References | 65 |
| Record of Edited Revisions | 66 |
| MGI Reagent Preparation Log (form) added in Troubleshooting section | 66 |
| Added steps to follow when there is no BioRad thermocycler available in Off board preparation section..... | 67 |

Introduction

The 2019-20 coronavirus outbreak is an ongoing public health emergency of international concern involving COVID. It is caused by the SARS-CoV-2, first identified in Wuhan, Hubei, China. SARS-CoV-2 is closely related to the original SARS-CoV. It is thought to have a zoonotic origin. Genetic analysis has revealed that the coronavirus genetically clusters with the genus Betacoronavirus, in lineage B or the subgenus Sarbecovirus together with the two bat-derived strains. It is 96% identical at the whole genome level to other bat coronavirus samples (BatCov RaTG13). In February 2020, Chinese researchers found that there is only one amino acid difference in certain genome sequences between the viruses found in pangolins and those from human patients, implying that pangolins may have been an intermediate host.

The BGI 2019-nCoV PCR assay is used with the MGI extraction system as a real-time RT PCR test able to detect the COVID gene (ORF1ab) and human beta actin gene.

Collection and Transport & Storage

Nasopharyngeal or throat swabs collected in viral transport media

- Store collected specimens at 4°C and process specimens as soon as possible.
- If delay of more than one day is expected, aliquot specimens and store at ≤-20°C.

Materials, Equipments and Facilities

Clean Room: Biosafety Cabinet (MIBCT3), freezer (MIFTG)



- Specimen Preparation area: Biosafety Cabinet (MIBCT7 or MIBCT8)
- BIO-RAD CFX96™ Real-Time System

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 4 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- BIO-RAD Hard-Shell PCR Plates 96-Well WHT/CLR
- BIO-RAD Microseal ®‘B’ seal Seals
- 96-Well Loading Block (pre-cooled to -20°C)
- Variable volume Rainin pipettes: 1 to 20 µL, 10 to 200 µL, 100 to 1000 µL
- Variable volume multichannel Rainin pipettes: 1 to 20 µL, 10 to 200 µL, 100 to 1200 µL
- Reagents: MGI Easy DNA/RNA, MGI Easy Magnetic Beads Virus DNA/RNA Extraction Kit, BGI Real time fluorescent RT-PCR kit for detecting 2019-nCoV,
- Positive and negative external Controls

Specimen Processing

Minimum volume for testing is 180 µL.

Heat fix 300 µL sample aliquots at 65°C for 30 minutes to inactivate virus.

***NOTE:** For dry swabs, refer the details to COVID 19 Pre-Analytic Processing Procedure from Microbiology Internal Manual*

Procedure

General Precautions


- There must be separate PCR work areas:
 - Clean room
 - Specimen preparation room
- Powder-free Gloves should only be in use in PCR areas.
- Change gloves frequently and keep tubes closed whenever possible.
- Prepare Working 1% sodium hypochloride daily.
- Specimen Preparation Supplies and equipment must be dedicated to Specimen Prep Area and not used for other activities and never used in Clean Room.
- Change lab coats and gloves between work areas.
- Use only Aerosol Resistant Tips (ART)
- Use only sterile RNase-free, DNase-free microtubes
- Thaw components completely at room temperature.
- PCR work areas (Clean Room and Specimen Preparation Area) bench tops and equipment after each shift.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 5 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Extraction Set up


1. Turn on the computer. Turn on the MGI platform.
2. Double click the **MGI SP 960 icon**.




3. Select **Real** in the 'Select the mode' window.




4. Click **Create** button.





| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 6 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


5. Authentication window pops up. Click **User Entry** button.




6. MGI SP960 is in the Initializing process.





| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 7 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |





- Click **Home** bar to go back the main menu
Select **Clean**.




| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|--------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 8 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

9. Select **Pre Clean**.
Click **Start** button.




| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|--------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 9 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

10. Confirm the operation listed in the confirmation window
Click **Continue** button.



11. MGI SP 960 will start the clean process. It takes 20 minutes to finish.





UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***


Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\


| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 10 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |



3. Choose **SARS CoV2 Reagent Setup** in the Solution drop-down menu



4. Choose **RNase Free Water.py** from **Script** drop-down menu to dispense water.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*


Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 11 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

5. Load tips, deep well plates and the tip lid according to the position on the screen.
6. Measure and transfer 45ml RNase Free Water into the tip lid at position 2.

Note: Make sure not to spill onto the surface of the instrument.





7. Press **Start** button.



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**. Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 12 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

8. Choose number of extraction plate (s) according to the table below
 Click **Continue** button

| Number of the Extraction Plate(s) | Total number of samples (patient samples + controls) |
|-----------------------------------|------------------------------------------------------|
| 1 | ≤96 |
| 2 | >96 and ≤192 |

 Require input


Comment Value

Number of Extraction Plates 1

| Name | Value | Comment |
|------|-------|---------|
| | | |

Continue Stop

Click **Continue** button

 Require input

Comment Value



Number of loading plates:

| Name | Value | Comment |
|------|-------|---------|
| | 1 | |
| | 2 | |
| | 3 | |


Continue Stop

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.



Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 13 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

10. After RNase Free Water is dispensed, take out the unused water and pour it back into the original container.
11. Seal the deep well plates with RNase Free Water and leave them on the bench at room temperature.
12. Dispense the next reagent-MW1.
 Choose **Buffer MW1.py** from **Script** drop-down menu to dispense MW1



13. Choose number of extraction plate (s) according to the table below
 Click **Continue** button

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 14 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

| Number of the Extraction Plate(s) | Total number of samples (patient samples + controls) |
|-----------------------------------|------------------------------------------------------|
| 1 | ≤96 |
| 2 | >96 and ≤192 |

Comment Value

Number of Extraction Plates

Name Value Comment


- Choose number of loading plate(s) (from 1 to 3) for MW1.
Click **Continue** button

Comment Value


Number of loading plate:

Name Value Comment


- Load tips, deep well plates and the tip lid according to the position on the screen.
- Measure and transfer enough volume of MW1 according to the screen into the tip lid.
Note: Make sure not to spill on the surface of the instrument.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 15 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

17. Label the deep well plates with MW1
18. Press **Start** button
19. After MW1 is dispensed, seal the deep well plates and leave them on the counter at Room temperature.
20. Pour the unused MW1 reagent back into the original container.
21. Dispense the next reagent MW2.



22. Choose number of extraction plate (s) according to the table below
Click **Continue** button

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 16 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

| Number of the Extraction Plate(s) | Total number of samples (patient samples + controls) |
|-----------------------------------|------------------------------------------------------|
| 1 | ≤96 |
| 2 | >96 and ≤192 |

23. Select the number of loading plate.

 Require input



24. For extraction of 96 samples choose 1, for 2 sets of 96 samples choose 2.

25. Dispense the next reagent MW2 buffer.



Set up tips, plates and tip lids according to the position on the screen. Label the deep well plates as MW2.

26. Measure the volume of the reagent according to the screen into the tip lid. Make sure not to spill on the surface of the instrument.

27. Press **Start**.


28. After the dispense seal the deep well plates and leave them on the counter at Room temperature.

29. Pour the unused reagent back into the original container.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|----------------------|
|   <p>Department of Microbiology</p> <p>Quality Manual</p> | <p>Policy # MI_MD_COVMGI</p> <p>Version: 3.12 CURRENT</p> | <p>Page 17 of 67</p> |
| <p>Section: Molecular Diagnostics Procedures</p> | <p>Subject Title: SARS-CoV-2 PCR by MGI</p> | |

Dispense reagents (massive procedure)

1. From the **Run Wizard** in the main menu, choose **0.SARS CoV2 BGI Reagent Setup (Massive)** in the **Solution** drop-down menu.
2. Choose **RNase Free Water (14 runs, 2-plate).py** from **Script** drop-down menu to dispense water.
3. Load tips, deep well plates and the tip covers according to the positions on the screen.
4. Measure and transfer 45ml RNase Free Water into the tip cover at position 3 and 4.
5. Press **Start** button to dispense RNase free water.
6. After RNase free water is dispensed, pour the unused RNase Free Water back into the original container.
7. Seal the deep well plates with RNase Free Water and leave them on the bench at room temperature.



The screenshot shows the 'Run Wizard' interface for a '0.SARS-CoV-2 BGI Reagent Setup' using the 'RNase-Free Water (14 runs, 2-plate).py' script. The interface is divided into two main sections: 'Operation Deck' and 'Records'.

Operation Deck: A grid of 24 positions (POS1 to POS24) is shown. POS1 contains a 'Full Tip Box 96 Tip Access' (TipGEBAF250A). POS2, POS3, POS4, POS6, POS7, POS8, POS10, POS11, POS12, POS13, POS14, POS15, POS16, POS17, POS18, POS19, POS20, POS21, POS22, and POS23 contain 'Deepwell Plate 1' through 'Deepwell Plate 14' (DeepwellPlateD77350504). POS5, POS9, POS24, and POS25 contain 'TipGEBCover'. POS18 and POS19 are labeled 'MagRack', and POS20 is labeled 'Shaker'. The 'Temp_Module' is also indicated.


Records: A table showing the sequence of activities and their completion status. The table has columns for 'No.', 'Time', 'Activity', 'Status', and 'Parameter'. The records show a sequence of actions from 11-05 08:28:43 to 11-05 08:30:22, including 'Ask', 'Reload deck', 'Workflow', 'Home', 'Message', 'Load Tips', 'Aspirate', 'Empty', 'Mix', 'Unload Tips', and 'Load Tips'.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 18 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


8. Similarly, choose **buffer MW1 (14 runs, 2-plate).py** from **Script** drop-down menu to dispense MW1 buffer.




The screenshot shows the 'Run Wizard' interface for a laboratory instrument. The top bar indicates the solution is '0.SARS-CoV-2 BGI Reagent Setup (Ma)' and the script is 'Buffer MW1 (14 runs, 2-plate).py'. The device status is 'Idle' and the completion status is 'Completed'. The interface is divided into an 'Operation Deck' and a 'Records' table.

Operation Deck: A 6x6 grid of 24 positions (POS1-POS24). Each position shows a thumbnail image of the instrument's deck and a label. POS1-POS4 show 'TipGEBCover' or 'Full Tip Box 96 Tip Access'. POS5-POS8 show 'Deepwell Plate 1-4' or 'Deepwell Plate 5-8'. POS9-POS12 show 'Deepwell Plate 1-4'. POS13-POS16 show 'Deepwell Plate 5-8'. POS17-POS20 show 'Deepwell Plate 9-12'. POS21-POS24 show 'Deepwell Plate 13-16'. POS18 and POS19 are labeled 'MagRack'. POS20 is labeled 'Shaker'.

Records Table: A table with columns: No., Time, Activity, Status, Parameter. It lists 19 records of activities such as 'Ask', 'Reload deck', 'Workflow', 'Home', 'Message', 'Load Tips', 'Aspirate', 'Empty', 'Mix', 'Unload Tips', and 'Aspirate'.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 19 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- Similarly, choose **Buffer MW2 (14 runs, 2-plate).py** from **Script** drop-down menu to dispense MW2 buffer.



Off board preparation

- Using the adjustable multichannel pipette, aliquot/pipet 200 μ L of samples into a labeled deep well plate.
- The negative control must be placed in the first well (A1) and the positive control in the last well (H8)

MGI # - Patient's samples from the space rack should be labeled in the same way as deep well plates.

e.g. MGI 1-1, MGI 1-2

MGI 2-1, MGI 2-2



MGI 3-1, MGI 3-2

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.


Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\


| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 20 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- Uncap each column of samples one by one and set aside the caps. Adjust the pipette to draw the volume into the pipette. Adjust the pipette again and deliver samples into the deep well plate. Cap each tube and put them on a different space rack at the same locations as on the original sample space rack.
- Put a gauze wet with disinfectant on the counter to catch any drips, also each time you pipetted a line touch the gauze to decontaminate your gloves. Repeat until one plate is done. Load the sample deep well plate into the MGISP according to the plate map. Do the second plate and load accordingly.
- Make the following Lysis buffer mixture into two of 50 mL tube if doing 2 extractions of 96, only 1 if doing 1 set of 96.

| REAGENT | VOLUME (1 PLATE) | PER TEST |
|------------------|------------------|----------|
| Buffer MLB | 16,000 µL | 160 µL |
| Absolute Alcohol | 20,000 µL | 200 µL |
| Proteinase K | 1,500 µL | 15 µL |
| Magnetic Beads | 1,500 µL | 15 µL |
| Enhancer | 100 µL | 1 µL |

- Pour into a reagent trough/ boat. Using a multichannel pipette dispense 360 µL into a clean deep well plate labeled as Lysis Buffer. Repeat until done. Load the lysis buffer plate according to the plate map. Repeat for another deep well plate if extracting 2 sets of 96 samples.
- For SARS-CoV-2 Extraction, go to Script and choose **Extraction 1 plate.py** if doing 1 plate of 96. Choose **Extraction 2 plate.Py** if doing 2 plates of 96.
- Check that all reagents/tips/ samples are loaded as the plate map. Label the Biorad PCR plate/s MGI machine#-Sample#
e.g. MGI1-1, MGI1-2
- Start
- Make your worklist according to the sample position into the deep well plate.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 21 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |




11. After extraction clean the platform, bag the left over sample deep well plates and secure the PCR plate/s.
12. For SARS-CoV-2 Amplification set up: From the Script drop down menu choose **1 Plate setup.py** if doing 1 plate set up. Choose **2 plate setup.py** if doing 2 plate set up.
13. Prepare PCR Master mix using reaction mix and enzyme mix in the clean room.
 - Option 1:** pipette 264 μ L of master mix into each well of the first column (A1-H1) of a deep well plate, MGI will dispense the master mix from well A1-H1 into a Biorad plate PCR plate on board.
 - Option 2:** Manually pipette 20 μ L of master mix into the PCR plate as detailed in **PCR Set-up** section below.
14. Load the PCR end product, tips and Master mix into the MGI SP

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 22 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

15. Press Start.



16. Seal the plate and load it into the Biorad Thermocycler.



17. If there are NO Biorad thermocycler available for

- ≤ 30 mins: put your plate in the clean fridge/reagent fridge near Virology Seniors Room;
- > 30 mins: put your plate in the -20°C freezer beside Virology Seniors Room

18. Leave the Worklist and a note (location of the plate) on the bench for the MLA and MLT to follow up.

PCR Set-up

In the Clean Room: (Change into dedicated clean room gown and gloves, work in Biosafety Cabinet)

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 23 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Remove the required boxes of BGI Real time fluorescent RT-PCR kit for detecting 2019-nCoV (50 reactions), for examples, 2 boxes if doing one extraction plate and 4 boxes if doing 2 extraction plates.

Storage at >4⁰C should not exceed 2 hours.

Prepare the number of Master Mix reactions required: Number of Test samples + 4 Controls + one extra test

Work quickly.

Mix gently. DO NOT VORTEX.

Make only enough master mix for the tests you are running.

| | |
|------------------------------|-----------|
| Number of Reactions | 1 |
| 2019-nCoV Reaction Mix | 18.5 |
| 2019 nCoV Enzyme Mix | 1.5 |
| Volume of Master Mix | 20 |
| Sample/Control Volume | 10 |

| No. of Test reactions | BGI nCoV PCR Mix | |
|------------------------------|-----------------------------|---------------------------|
| | 2019-nCoV Reaction Mix (µL) | 2019 nCoV Enzyme Mix (µL) |
| 1 | 18.5 | 1.5 |
| 50 (48rxn) | 925 | 75 |
| 100 (1 plate) | 1850 | 150 |
| 200 (2 plates) | 3700 | 300 |

To be loaded manually:


- Place 96-well plate onto pre-cooled block

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 24 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- Pipette **20** μ L prepared PCR Master Mix into each reaction well of the Biorad Hard Shell PCR plate using the multichannel pipette or repeater pipette.

On the MGI SP 960:

From the drop down menu

Solution: Choose **SARS-CoV-2 Amp Plate Set up**


Script: 1 plate setup (no aliquoting).py


If doing 1 extraction plate

From Script choose **2 plate setup (no aliquoting).py** if doing 2 extraction plates

Label your Mastermix according to run MGI number and run number.


Load according to the plate maps in MGISP 960



| | | |
|-------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|----------------------|
|  <p>Department of Microbiology</p> <p>Quality Manual</p> | <p>Policy # MI_MD_COVMGI</p> <p>Version: 3.12 CURRENT</p> | <p>Page 25 of 67</p> |
| <p>Section: Molecular Diagnostics Procedures</p> | <p>Subject Title: SARS-CoV-2 PCR by MGI</p> | |


Bio-Rad Detection Area:


- Double click **Bio-Rad CFX Maestro** icon
- Startup Wizard window pops up
- Make sure **CFX96** instrument is selected







- Run Set up window pops up, including three tabs-Protocol, Plate, and Start Run
- Click Open Lid button to open the lid

- Under Select run type, click User-defined button

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 26 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |





- Run Set up window pops up, including three tabs-Protocol, Plate, and Start Run
- Click Open Lid button  to open the lid
- Load sealed plate to the block
- Click Close Lid button  to close the lid
- WARNING! Do NOT manually close the motorized lid
- Under Protocol tab, select **BGI COVID.pcrd** from Express Load pull-down menu






Click Plate tab  on the top or click Next button  on the bottom right side to load the plate profile



- Select plate profile BGI COVID.pltd from the Express Load pull-down menu





| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 28 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


- Click Start Run tab  on the top or click Next button  on the bottom right side



- Confirm the following information:
 Protocol: **BGI COVID.pcrd**
 Plate: **BGI COVID.pltd**
 Scan Mode: **All Channels**
CFX Thermocycler instrument number: checked on
- If any information above needs to be edited, click Prev button 
- If all information is correct, click Start Run button 
- Save Optical Data File [CT014845] window pops up
- Change “admin” to MGI 1-Run no.
- Save the file in the designated folder:

T:Microbiology>Virology>Bio Rad COVID MGI 1-1>Save Run>

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 29 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


- Click Save button 

Build Worklist


- Open the worklist file according to the following path:


T:Microbiology>Virology>Bio Rad CFX96 PCR>Worklist Master copy

- Password window pops up, click **Read only** button



- Scan the samples' information

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 30 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |



- File > Save as CSV (Comma delimited) at


T:Microbiology>Virology>Bio Rad MGI 1-1>Import Worklist

File name: MGI 1-1yyyy.mm.dd.run No. eg.MGI 1-1_2020.05.20.1

Save as type: CSV (Comma delimited)

Note: Worklist can only be imported as CSV type.

- Click **Save** button
- Click **OK**
- The following window pops up, click **Yes** button




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.



Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 31 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- Close current excel file without saving.



Import Worklist

- Click **Realtime Status** tab
- Select **View/Edit Plate...** from **Plate Setup** pull-down menu



- Plate Editor window pops up

- Click Spreadsheet View/Importer button




- Plate Spreadsheet View window pops up

*NOTE: This document is **Uncontrolled When Printed.***


Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.


- Select **HEX** from **Fluors List** pull-down menu




- Click **Import** button 
- **Select Import File** window pops up

T:Microbiology>Virology>Bio Rad MGI 1-1 yyy-mo-date>Import Worklist



- Select appropriate worklist and click **Open**, all samples and controls information are imported
- Click **OK** button
- Select all wells by clicking top left side corner
- Type **beta actin** as target Name for **HEX** channel
- Hit **Enter** key on the key board
- Click **OK** button 
- Apply changes, click yes.





UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.


Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 35 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Click **Time Status** tab





from the Run Details window



Analysis

After the PCR is completed:

- Change the channel name from **Fluorophore** to **Target**
- Mark **beta actin** box under the graph area with Covid gene unmarked, click any point to the left of all the beta actin curves, drag the cursor to select all the graphs above the threshold. Examine and mark on the worksheet any sample and control with no beta actin or Ct of beta actin greater than 35.

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 36 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

For ALL samples from any client on the BGI on a 96 well plate:

- IF there is a widespread range of beta-actin curves
- IF >5 are invalid


Then:

- Adjust threshold of **beta-actin** to **100 RFU** line.
- Repeat remaining invalids on Seegene


NOTE: If there is NOT a widespread range of beta-actin curves and more than 5 invalids, DO NOT adjust curve, repeat invalids on Seegene.

❖ **Please consult microbiologists or/and seniors if you have any concerns about the adjustment of beta-actin threshold**


Example for a tight-spread range of beta-actin curves:




Example for a wide-spread of beta-actin curves:



- Mark **Covid gene** box and **log scale** box under the graph area with beta actin unmarked. Click positive control only, and drag the threshold line to the mid-point of the exponential amplification (linear portion of the sigmoidal amplification curve), and then unmark log scale and select all the samples. Examine and mark on the worksheet the result of all samples and controls for Covid gene.




- A positive PCR is observed if there is a rise or amplification in the Target channel e.g. FAM (Covid gene), & HEX (beta actin). The graph should be exponential and sigmoidal in shape. Conversely a negative PCR is a flat line or no signal.
- Click **Export**
- Select **Custom Export...**
- Click **Export** button



- “Save as” Window pops up
- T:\microbiology\Virology\Biorad CFX COVID MGI 1-1\COVID Excel results



- Click **Save** button
- The final results with Ct values are exported in Excel sheet
- Print the excel result sheet and attach it to the COVID Worksheet.
- Close excel.
- Click OK button.




ALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY



*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 39 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

- Click **Close** button.
- **LIS Interfacing:** exclude from analysis any sample that is not expected to cross over to LIS, e.g. low positive sample from Seegene repeated as negative in the current MGI run.
Click **Export > Custom Export**
Choose **Export Format** as **Text (*.txt)**
Click **Export** button and save the text file to T:\microbiology\Virology\MGI
Close the text file that is automatically opened.
- Close Biorad CFX96 Maestro Real-Time System. **Do you want to save the changes to MGI 1-1yyyy.mm.dd CT038327.pcrd?**
- Press **Yes**.
- Shut down computer
- Shut down Biorad CFX 96 Thermocycler.



| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 41 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Interpretation

Quality Control (QC) Samples

Assess QC samples for run validity prior to interpreting patient samples.


| QC Sample | VIC channel | FAM channel | Interpretation |
|---------------------|---------------------------------------------------|---------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No template control | No amplification | No amplification | Pass – Proceed to sample analysis |
| Positive control | Sigmoidal amplification curve and Ct value is <35 | Sigmoidal amplification curve and Ct value is <37 | Pass – Proceed to sample analysis |
| | | | |
| No template control | Sigmoidal amplification curve and Ct value is <35 | Sigmoidal amplification curve and Ct value is <37 | Fail - Repeat all positive and “repeat” samples from the run and report the negatives as negative. <i>Note: if curve is not sigmoidal, contact microbiologist-on-call for further instructions.</i> |
| Positive control | No amplification or Ct value is >35 | No amplification or Ct value is >37 | Fail - Do not report sample results. Repeat run |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.



Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 42 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Test Samples

Note: Interpretations below are based on **Sigmoidal** amplification curves.

| Sample | FAM (COVID gene-ORF1ab) | HEX (beta actin) | INTERPRETATION |
|--------|------------------------------------------------------------------------------------------------------------------|---------------------|---------------------------------------------------------------------------------|
| A | Ct value ≥ 37 | Ct value < 35 | SARS-CoV-2 RNA NOT Detected |
| B | Ct value < 35 | Ct value < 35 | SARS-CoV-2 RNA Detected |
| C | $35.0 \leq Ct < 37.0$ | Ct value < 35 | SARS-CoV-2 RNA Detected – Low level <i>Repeat on Seegene</i> |
| D | Ct value ≥ 37 | Ct value ≥ 35 | INVALID |
| E | Ct < 37.0 | Ct value ≥ 35 | Questionable Invalid <i>Repeat on Seegene</i> |
| | | | |
| Repeat | Follow General BGI Reporting Process and BGI Reporting Process of Repeat Samples | | FINAL RESULT for REPEAT of LOW LEVEL DETECTION |

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 43 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


General BGI Reporting Process

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\





UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is **Uncontrolled When Printed**.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 45 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


BGI Reporting Process of Repeat Samples


UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\



| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 47 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Reporting

1. PRELIMINARY POSITIVE RESULTS (for UHN/MSH/STAT patients only)

- Result in **the isolate window** as **COVID-19 virus**
- In the Isolate comment, choose **\DPCR & \indM (Low Prelim MGI)**
- Status as **PRELIMINARY**
- Inform all positive results according to **Isolate Notification and Freezing Table**

Example:

COVID-19 virus

DETECTED by real-time PCR.


* * * * *

(Low level detection)

Sample is being retested. Further report to follow.

Testing performed using the BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay.

NOTE: The BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 48 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

2. FINAL RESULTS

2.1 NEGATIVE RESULTS

- Result in the test window using COVID-19 Virus Not Detected test comment
- Choose from the keypad: }-NCB (Neg BGI)
- Status as FINAL
- When BGI .txt file for a run is exported to LIS, all negative samples should auto-finalize.


Example:

COVID-19 virus NOT detected by real-time PCR.

* * * * *

Testing performed using the BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay.

NOTE: The BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 49 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

2.2 POSITIVE RESULTS

- Result in **the isolate window** as **COVID-19 virus**
- In the Isolate comment, choose **\DPCR & \CVB+ (BGI)**
- Status as **FINAL**
- Inform all positive results according to **Isolate Notification and Freezing Table**

Example:



COVID-19 virus

DETECTED by real-time PCR

* * * * *

Testing performed using the BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay.

NOTE: The BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory.

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 50 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

2.3 INVALID RESULTS

- Result in the test window using COVID-19 virus PCR test unable to be completed test comment
- Choose from the keypad: }INCB (BGI INV)
- Status as **FINAL**
- Inform all invalid results according to **Isolate Notification and Freezing Table**



Example:

COVID-19 virus PCR test unable to be completed.

* * * * *

No gene targets were detected including the human beta-actin gene suggesting inhibition of the BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay PCR reaction or inadequate sampling. Please submit another sample for testing if clinically indicated.

NOTE: The BGI Real-Time Fluorescent RT-PCR 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory.

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 51 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

3. FINAL RESULT for REPEAT of LOW LEVEL DETECTION

Follow [BGI Reporting Process for Repeat Samples](#)

3.1 CONFIRMED POSITIVE RESULTLS

Follow [FINAL RESULTS for POSITIVE](#) reporting procedure.

- Finalize as COVID-19 virus Detected instead of “Low level detection confirmed”.

3.2 INDETERMINATE RESULTLS

- **Suppress** COVID-19 virus isolate
- Choose from the keypad: }**INDC (Indeterm)**
- Status as **FINAL**
- Inform all invalid results to the ward /clinic according to **Isolate Notification and Freezing Table**

Example:

INDETERMINATE for COVID-19 virus.

* * * * *



Results should be interpreted within the context of the clinical signs, symptoms, and history of the patient.

INDETERMINATE result may indicate the presence of low levels of virus, non-specific reactivity of the assay, or other unrecognized factors. Please submit a follow-up sample if clinically indicated.

Testing performed using the Seegene Allplex 2019-nCoV Assay.

NOTE: The Seegene Allplex 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory

* * * * *

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 52 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

3.3 NSQ to Repeat RESUTLS

- Initially COVID-19 virus **Detected – Low level but Non-Sufficient Quantity to repeat**
- Update Isolate Comment from the keypad: \DPCR + \liq
- Status as **FINAL**
- Inform all invalid results according to **Isolate Notification and Freezing Table**

Example:

COVID-19 virus

DETECTED by real-time PCR.


* * * * *

(Low level detection-insufficient quantity of specimen to confirm)

* * * * *

Testing performed using the Seegene Allplex 2019-nCoV Assay.

NOTE: The Seegene Allplex 2019-nCoV Assay has been approved by Health Canada for Emergency Use Access (EUA) and has been verified by the University Health Network/Sinai Health Microbiology Laboratory

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 53 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Cleaning

At shift end, perform cleaning protocol as outlined below:

Clean Room: Wipe down with RNAse Away or NucleoClean on paper towel, followed by distilled water, and then 70% alcohol

Biological Safety Cabinet

Pipettes

Bench tops

Specimen Preparation Area: Wipe down with Working 1% hypochloride (made daily), followed by distilled water, and then 70% alcohol

Biological Safety Cabinet (BSC), pipettes, centrifuge, and bench top.

Seal and discard BSC waste

Wash racks


Space Racks, and Multi-channel Pipette: Wipe down surfaces with RNAse Away or NucleoClean on KimWipe, followed by UltraPure water, and then 70% alcohol Wipe

Amplification Area: Wipe down surfaces with RNAse Away or NucleoClean on KimWipe, followed by UltraPure water, and then 70% alcohol Wipe

Seal & discard reaction microtubes into biohazard waste after each run.

Perform the cleaning procedure according the daily maintenance sheet.

MGI SP 960: wipe any spills with 70% Alcohol and do a post clean.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 54 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Quality Control

Reagent QCs:

An External Control (external to MGI BGI assay) is used to monitor the isolation, amplification and detection procedures. The result must correspond to expected value supplied by the manufacturer.


A positive and negative external control for Covid gene should be extracted on the MGI run



- with each new lot/new shipment of RNA/DNA MGI Easy Extraction kit and BGI *Real-Time Fluorescent RT-PCR 2019-nCoV* Assay and
- after major maintenance

Daily QCs:

Every run:



- Each patient specimen must show beta actin to monitor both extraction and PCR inhibition.
- Both a Positive and Negative control including target and beta actin will be included.
- For QC result acceptability refer to interpretation section.
- As a way of avoiding plates swap from extraction to reporting, please follow the graph below for addition controls:
Plate 1: add POS control in A2 well
Plate 2: add NEG control in A2 well

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 55 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.
 Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\


| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 56 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Report all failed QCs to senior/charge technologist.

Failed QC:

Test is invalid without satisfactory QC results.

- a. Do not release results pending resolution of QC failure.
- b. Inform charge/senior technologist.
- c. Record in Reagent Log Chart, Instrument Maintenance Log or Incident Report where appropriate.
- d. If the QC failure was due to a simple matter of position reversal or misplacement, the run can be released (positive QC material yielded positive result, negative yielded negative result).
- e. If negative QC material yielded positive result, it may be due to cross-contamination from adjacent positive sample within the run or carry-over contamination from previous runs via equipment or the environment. Review procedure and equipment to establish and eliminate potential sources of contamination.
- g. The extent and nature of contamination can also be evaluated by comparing the positive rate of the run with its expected positive rate.
- h. If the contamination is extensive, it is necessary to quarantine/discard potentially contaminated reagents and consumables and disinfect equipment and environment before repeating the run.
- i. If a carry-over contamination is suspected (e.g. two or more runs with negative QC being positive or patient samples have higher than expected positive rate and these samples are often non-repeatable positives), it is necessary to have a thorough environmental disinfection followed by swabbing to monitor.
- j. Successful ending to a carry-over contamination may be indicated by QC results and patient positivity rate falling back to the expected normal range and three negative environmental swabs.

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 57 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Troubleshooting

1. Technical Support

Collect information by taking pictures of reagents/supplies and run/area of issue, printing screen and exporting logs) and send them to MGI Canada Service Support (MGICSUPPORT@mgiamericas.com)


For details, please refer to [MGISP960 Troubleshooting Guide](#).

For urgent support, phone BGI technician support staff immediately.

2. Common Troubleshooting Tips

Non-Specific curves:

The non-specific curve (exemplified in the following picture) should be excluded from interfacing and reported as negative manually. If you still have concerns, please let the senior in your shift or microbiologist on call know.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 58 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


3. Failed run, possible Root cause and Graph
 a. Eluate plate use in PCR




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|  Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 59 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

b. Plate not sealed





UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.



Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 60 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |


c. Contaminated run





*NOTE: This document is **Uncontrolled When Printed**.*
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.
 Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 61 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

d. Fanning curves (uneven distribution of Mastermix)





| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 62 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

4. MGI workflow checklist

Many manual steps are needed in the MGI testing process. All the steps in the procedure should be followed to avoid potential sample mix-ups and plate swaps.

The following MGI workflow checklist is used to check and monitor all critical steps from aliquoting to reporting.

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 63 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

MGI WORK FLOW CHECKLIST

DATE:

| STEPS TO CHECK * | Initial § | Comments |
|-------------------------------------------------------|-----------|----------------------------------------------------------------|
| Aliquoting | | Please export STP worklist to STP folder if using MGI STP 7000 |
| MIANT 32-1/STP-1 | | |
| MIANT 32-2/STP-2 | | |
| Manual Aliquot | | |
| CONTROLS ADDED by | | |
| Plate 1 (2 Pos, 1 Neg) | | |
| Plate 2 (1 Pos, 2 Neg) | | |
| Lysis Buffer prep by | | |
| Reagents (H ₂ O/MW1/MW2) volume checked by | | |
| MGI extraction started by | | Check if all plates are in correct positions |
| Eluate/Mastermix volume checked by | | |
| Import Worklist and Time Stamp | | Convert STP/manual worklist to CSV |
| MGI 1-1 (Worklist created/time stamped by) | / | |
| MGI 1-2 (Worklist created/time stamped by) | / | |
| MGI 2-1 (Worklist created/time stamped by) | / | |
| MGI 2-2 (Worklist created/time stamped by) | / | |
| MGI 3-1 (Worklist created/time stamped by) | / | |
| MGI 3-2 (Worklist created/time stamped by) | / | |
| MGI 4-1 (Worklist created/time stamped by) | / | |
| MGI 4-2 (Worklist created/time stamped by) | / | |
| MGI 5-1 (Worklist created/time stamped by) | / | |
| MGI 5-2 (Worklist created/time stamped by) | / | |
| MGI 6-1 (Worklist created/time stamped by) | / | |
| MGI 6-2 (Worklist created/time stamped by) | / | |
| Biorad PCR/Reporting | | |
| Master Mix Prep | | |
| Load & Check PCR Plate Plates | | Make sure plate 1 & 2 at correct positions |
| Worklist Imported by | | |
| Analyzed by | | |
| Biorad checked by | | Open Biorads & check correct plates loaded before exporting |
| Exported by | | |
| Reported by | | |

*Please initial each step and pay attention to newly added critical check points as highlighted.



§ No need to initial in shaded area

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is *Uncontrolled When Printed*.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 64 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

5. MGI Reagent Preparation Log

Insufficient reagent volume during reagent preparation may cause run failure of various forms. Use of the following log sheet facilitates good documentation and tracking of MGI reagents.

Month/Year _____ Reviewed Monthly By Senior (initial) _____



| Day | MW2 Plates | MW1 Plates | H ₂ O Plates | Tech Initials/Shift | Lysis Buffer Vials | Tech Initials/Shift | NEGATIVE CONTROL Boxes | Tech Initials/Shift |
|-----|------------|------------|-------------------------|---------------------|--------------------|---------------------|------------------------|---------------------|
| 1 | | | | | | | | |
| 2 | | | | | | | | |
| 3 | | | | | | | | |
| 4 | | | | | | | | |
| 5 | | | | | | | | |
| 6 | | | | | | | | |
| 7 | | | | | | | | |
| 8 | | | | | | | | |
| 9 | | | | | | | | |
| 10 | | | | | | | | |
| 11 | | | | | | | | |
| 12 | | | | | | | | |
| 13 | | | | | | | | |
| 14 | | | | | | | | |
| 15 | | | | | | | | |
| 16 | | | | | | | | |
| 17 | | | | | | | | |
| 18 | | | | | | | | |
| 19 | | | | | | | | |
| 20 | | | | | | | | |
| 21 | | | | | | | | |
| 22 | | | | | | | | |
| 23 | | | | | | | | |
| 24 | | | | | | | | |
| 25 | | | | | | | | |
| 26 | | | | | | | | |
| 27 | | | | | | | | |
| 28 | | | | | | | | |
| 29 | | | | | | | | |
| 30 | | | | | | | | |
| 31 | | | | | | | | |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\



| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 65 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Related Documents

| |
|------------------------------------------------------------------------|
| Virology Accessioning Manual |
| Biorad Worksheet |
| Qualitative PCR External QC |
| Nucleic Acid Extraction – Biomerieux NucliSENS easyMAG |
| Pipetting by epMOTION Manual |
| Training Checklist |

References

MGI Easy DNA/RNA Extraction
 BGI Real-time fluorescent RT PCR kit for detecting nCoV

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 66 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

Record of Edited Revisions

Manual Section Name: SARS-CoV-2 PCR by MGI

| Page Number / Item | Date of Revision | Signature of Approval |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|-----------------------|
| Addition of neg and pos control positions on plate | July 17, 2020 | Dr. T. Mazzulli |
| Changed No template control interpretation when CT<37 from Repeat run to Repeat all positive and “repeat” samples and report the negative as negative. | September 14, 2020 | Dr. T. Mazzulli |

Full document review included in all updates. Bi-annual review conducted when no revision had been made within 2 years.



| Page Number / Item | Date of Revision | Edited by: |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|----------------------|
| Changed the analyzing method for the graphs of target gene and beta actin; Added the massive procedure for making extraction reagents; Added LIS interfacing procedure. | November 18, 2020 | Dorna Zareianjahromi |
| <ul style="list-style-type: none"> • Addition of General Seegene Reporting Process, BGI Reporting Process for Repeat Samples; • Updated BGI reporting rules, preliminary, final and repeat sample report statements in test window and isolate window • Addition of troubleshooting | January 08, 2021 | Dorna Zareianjahromi |
| Minor formatting change | April 11, 2021 | Jessica Bourke |
| MGI workflow checklist added in the trouble shooting session Graphs for additional controls added | May 17, 2021 | Oliver Li |
| MGI workflow checklist updated with volumes of reagents/eluates/mastermix to be checked | June 14, 2021 | Oliver Li |
| MGI Reagent Preparation Log (form) added in Troubleshooting section | July 22, 2021 | Oliver Li |
| Updated technical support by adding details as to how to collect information of runs of issue | September 20, 2021 | Oliver Li |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Molecular Diagnostics Procedures\

| | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|---------------|
|   Department of Microbiology Quality Manual | Policy # MI_MD_COVMGI | Page 67 of 67 |
| | Version: 3.12 CURRENT | |
| Section: Molecular Diagnostics Procedures | Subject Title: SARS-CoV-2 PCR by MGI | |

| Page Number / Item | Date of Revision | Edited by: |
|--------------------------------------------------------------------------------------------------------------|--------------------|------------|
| Added steps to follow when there is no BioRad thermocycler available in Off board preparation section | September 23, 2021 | Oliver Li |
| Updated BGI reporting flowchart for repeating samples | October 22, 2021 | Qin Liu |
| Added criteria for the adjustment of beta-actin threshold | October 27, 2021 | Qin Liu |
| Updated E gene only in “General BGI Reporting Process” | November 05, 2021 | Oliver Li |
| Minor formatting changes | November 08, 2021 | Oliver Li |
| Added trouble shooting for failed run: Possible root cause and Graph | March 11, 2022 | Qin Liu |
| | | |