

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/11/v03	Page 1 of 5
Section: Serology Manual	Subject Title: West Nile (Flavivirus) IgG	
Issued by: LABORATORY MANAGER	Original Date: April 25, 2001	
Approved by: Laboratory Director	Revision Date: May 27, 2003	

WEST NILE (FLAVIVIRUS) IgG

I. Introduction

The Focus Technologies Flavivirus (West Nile) ELISA IgG is intended for qualitatively detecting IgG antibodies to flaviviruses(West Nile) in human serum. In conjunction with the Focus Technologies Flavivirus (West Nile) ELISA IgM Capture ELISA, the test is indicated for testing persons having symptoms of arbovirus infection, as an aid in the presumptive diagnosis of flavivirus infection.

II. Specimen Collection and Processing

5 ml of blood is collected in a serum separator tube and separated by centrifugation. The serum is removed from the clot and refrigerated until testing. Specimens are stored at – 70⁰C after testing.

III. Procedure

- i) Reagent Preparation:
 - a. Wash Buffer (100 ml) 10X
Prepare working buffer (1X) by adding 100 ml of 10X Wash Buffer to 900 ml of Distilled H₂O. Mix well. Labeled and dated
 - b. Controls: 1 vial of Positive Control
1 vial of Negative Control
 - c. 1 vial of Cut-off Calibrator
 - d. IgG Conjugate
 - e. Substrate Reagent
2 vials of Tetramethylbenzidine (TMB) and Hydrogen Peroxide in buffer
 - f. Stop Reagent
1 vial of 1M sulfuric acid
- ii) Other Materials:
 - IgG Antigen Wells,96 wells
 - Sample Diluent

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/11/v03	Page 2 of 5
Serology Manual		

iii) Method:

1. Label a set of 12x75 mm glass tubes as follows: Blk, Pos,Neg, Cal,7,8,9,10.....,and External Control (when required, e.g. a new lot #).
2. Add 1 ml of Sample Diluent to all tubes.
3. Leave 1st tube blank; add 10 ul of Controls, Calibrator and patient samples to corresponding tubes.
4. Take out the required # of strips (Open bag containing strips when strips reaches room temperature, ~10-15 minutes). Return the unused strips to the bag, and refrigerate immediately.
5. Change bottle from D H2O to Wash Buffer.
 - a. Turn Washer 'ON'.
 - b. Press button under 'Prime' (F2).
 - c. Press button under 'washing' (F1).
 - d. 'Select Program': Press memory (F1), press Number (F1).
 - e. **For dispensing 350 ul of buffer use Program '2'.**
 - f. 'Program No : press '2',enter.
 - g. '002 Name 1'should appear.
 - h. Press 'Set' at this point, press 'Change' (F2), then press 'Enter' until 'Number of strips to wash' appears-check that # of strip appears is what is needed. If it is the same, press 'Escape'. If is not the same, enter the right # of strip, and press 'enter 'until ' 002 Name 01' appears, press 'Enter'.
 - g. Press' Run' (F1), and press 'No' (F2) under 'Prime'.
- h. The Washer will dispense the 350 ml of buffer into each well, and soak for 5 minutes.
6. After 5 minutes. Aspirate Wash Buffer from wells.

For aspiration of Wash Buffer, use Program ' 1 '.

 - a. 'Select Program': Press memory (F1), press Number (F1).
 - b. 'Program No: press 1,enter.
 - c. '001 Name 01'should appear.
 - d. Press 'Set' at this point, press 'Change' (F2), then press 'Enter' until 'Number of strips to wash' appears-check that # of strip appears is what is needed. If it is the same, press 'enter'. If is not the same, enter the right # of strip, and press 'enter 'until ' 001 Name 01' appears, press 'Enter'.
 - e. Press' Run' (F1), and press 'No'(F2) under 'Prime'.
 - f. The Washer will aspirate the buffer.
 - g. 'End washing' will appear when washing cycle is finished, remove plate.
 - h. 'Repeat Washing?' will appear, press 'No' (F2).

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/11/v03	Page 3 of 5
Serology Manual		

7. Add 100 ul of Blk to well # 1, 100 ul of Pos to well # 2, 100 ul of Neg to well # 3, 100 ul each of Cal to well # 4, 5, and 6. Add 100 ul of specimen to each well, starting from well # 7. Use last well as external Control well when needed.
8. Cover plate with sealing tape, and incubate for 60 minutes at Room Temperature.
9. Wash the wells three times using the Automatic Washer under Program '10'.
 - a. 'Select Program': Press memory (F1), press Number(F1).
 - b. For washing and soaking, use Program '10'.**
 - c. 'Program No : press 10,enter.
 - d. Program '10 Name IgG WN' should appear. Press 'Set' at this point, press 'Change' (F2), then press 'Enter' until 'number of strips to wash' appears-check that # of strip appears is what is needed. If it is the same, presses 'enter'. If is not the same, enter the right # of strip, and press 'enter 'until ' '11 Name IgG WN' appears, press 'Enter'.
 - e. Press 'Run' (F1), and press 'No' (F1) under 'Prime'. Will wash strip for 3 times, each time with working buffer.
 - f. 'End washing' will appear when washing cycle is finished.
 - g. 'Repeat Washing?' will appear, press 'No' (F2).
10. Remove plate from Plate Washer.
11. Add 100 ul of Conjugate into each well using an 8-channel pipette, cover plate with sealing tape and incubate for 30 minutes at room temperature.
12. Repeat wash steps -9a-g.
13. Add 100 ul of Substrate into each well using an 8-channel pipette, cover and incubate for 10 minutes at room temperature.
14. Add 100 ul of Stop Reagent into each well using an 8-channel pipette. In antibody-positive wells, color should change from blue to yellow.
15. Read the plate in spectrophotometer at dual wavelength 450 nm/630 nm within 1 hour of stopping the assay.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/11/v03	Page 4 of 5
Serology Manual		

IV. **Calculation:**

Index value= Optical Density of Sample/the mean Optical Density of the Cut-Off Calibrators

Interpretation of Index Values:

Positive: > 1.10

Negative: < 0.9

Equivocal: >=0.9 and <= 1.1

All positive and equivocal specimens should be repeated the next day .

V. **Validation:**

The mean value for the Cut-Off Calibrator must be within 0.100 to 0.700 O.D. units minus the blank well.

The Positive Control Index value should be between 1.5 and 3.5.

The Negative Control Index value should be less than 0.8.

VI. **Reporting**

Positive result: Sent to PHL for confirmation. Report PHL result in LIS

Negative result: West Nile Virus(Flavivirus) IgG by EIA: Negative.
This is a research test.

Equivocal result: repeat testing the next working day.

If result is negative,report as : West Nile Virus(Flavivirus) IgG
by EIA: Negative.
This is a research test.

If result is still equivocal /positive: Send to PHL for
confirmation. Enter PHL for
result in LIS.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/11/v03	Page 5 of 5
Serology Manual		

VII. Quality Control:

Each plate must include 1 Positive Control, 1 Negative Control, and 3 Cut-off Calibrators.

An External Positive Control (Accurun 165 West Nile virus IgM/IgG Positive Control Series 5000) must be included with each new lot. Result filed in Reagent Lot Binder.

Reference

Manufacture’s insert, Focus Technologies, Cypress, California 90630.U.S.A. 2003.