

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/17/04/v01	Page 1 of 5
Section: Serology Manual	Subject Title: Appendix IV The Thermo-Cycler Gene Amp PCR System 9600	
Issued by: LABORATORY MANAGER	Original Date: March 14, 2001	
Approved by: Laboratory Director	Revision Date: September 20, 2001	

APPENDIX IV

THE THERMO-CYCLER USED IN AMPLIFICATION IS GENE AMP PCR SYSTEM 9600

To use this machine for amplification, create the Method and the program for the experiment.

Creating a Method:

A method is a series of up to 17 programs linked together.
For example, you can link a HOLD program to a CYCL program and conclude you experiment. Access the main menu.

Select option 9600 <u>R</u> N-CREATE-EDIT-UTIL
--

1. Press the **OPTION** key to select CREATE, then press **ENTER**, the following display appear

Create program <u>H</u> OLD-CYCL-AUTO-METH

2. Press the **OPTION** key to select METH, then press **ENTER**. The following display appears.

Link progs: _ _ - - - - -

3. Type the program numbers of the programs you want linked. After you type each number, press **ENTER**. For example, to link programs 1,3,and 5, your display would look like this:

Link progs: 1- 2- 3- 5- - -

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/17/04/v01	Page 2 of 5
Serology Manual		

- Press **STEP** to move to the next display. If you are linking seven or more programs in a method, the following display appears:

- - - -
- - - -

This display let you link up to ten additional programs in the method.

NOTE: You cannot link a method.

- Type the program numbers of the programs you want linked. After you type each number, press **ENTER**.
- Press **STEP** to move to the next display.

Meth # ???	25.6 °C
<u>R</u> N-STORE-PRINT-HOME	

This display allows you to run, store, or print the method.

- Press **OPTION** key to move the cursor beneath **STORE**, then press **ENTER**. The following display appears:

Store
Enter program # <u>x</u>

X is the next available number.

- Press **ENTER** to store the method under the displayed number (1-150), the following display appears:

Store
Protect program? <u>N</u> O

Now your method stored, so you need to program your programs which linked in your method To do this follow the steps as following:

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/17/04/v01	Page 4 of 5
Serology Manual		

2. Press ENTER to store the method under the displayed number, OR
 To use a different number, type the number (1-150), then press ENTER.
 The following display appears:

Store Protect program? <u>NO</u>

3. The above display gives you the option of protecting the file under a user I.D.
 By doing so, the file cannot be overwritten or deleted unless the user number is
 entered. Select YES or NO by pressing the **OPTION** key, then press **ENTER**.
 [Always selects NO unless if you want to keep it forever and will be under your I.D.
 number].

CYCLE: The Cycle Program

A CYCL program contains the thermal ramp and hold segment patterns for PCR cycling.
 It contains up to nine set-points. Each set-point comprises a ramp segment and a hold
 segment for the specific target temperature. The pattern can be repeated up to 99 times.
 Typically there are two or three set-points in a CYCL program. A complete cycle
 typically required less than two to four minutes.

The Displaying in a CYCL Program: example for displays CYCL program

<u>3</u> Temperature PCR
Setpt #1 Ramp 0:0 <u>0</u> 94.0C Hold 0:30
Setpt #2 Ramp 0:0 <u>0</u> 55.0C Hold 0:30
Setpt #3 Ramp 0:0 <u>0</u> 72.0C Hold 0:30
Total cycles = 35 Pause during run? NO

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI/SER/17/04/v01	Page 5 of 5
Serology Manual		

Adjust your program temperature and number of cycles according your experiment protocol by following the table above. Store your programs as usual under the same number you create.