 Department of Microbiology Quality Manual	Policy # MI_EDU	Page
	Version: 1.1 CURRENT	1 of 56
<b>Section:</b> Bacteriology Procedures	Subject Title: <b>Education Manual</b>	
Prepared by QA Committee		
Issued by: Laboratory Manager	Revision Date: 9/14/2018	
Approved by Laboratory Director: Microbiologist-in-Chief	Annual Review Date: 5/1/2019	



**Uncontrolled When Printed**

Microbiology Rotation for Medical Laboratory Technology Students .....	2
Microbiology Rotation for Residents.....	7
Training Objectives for Adult Infectious Diseases: Medical Microbiology Rotation.....	11
<b>Resident Objectives at MSH</b> .....	12
Mycology - Training Manual for Residents.....	27
General Microbiology Review.....	30
APPENDIX I - Day One.....	31
APPENDIX II- Media.....	32
<b>General Overview of Weekly Bench Rotations</b> .....	39
Mycology - Training Manual for Students .....	44
Books in Microbiology Laboratory .....	51
Record of Edited Revisions .....	55

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 2 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Microbiology Rotation for Medical Laboratory Technology Students

### Welcome to Microbiology!

Microbiology is an interesting and dynamic science. Your time spent in our lab should be an enjoyable and successful learning experience.

The main objective of this training program is to produce a technologist who has the ability to apply theoretical knowledge to the practical aspects of the Microbiology laboratory, efficiently and with a caring attitude.

During your rotation you will be exposed to all aspects of the department. You will handle specimens from their receipt in the department through to the reporting of results. You will develop the ability to make a presumptive identification of common pathogens and be able to select the correct tests to fully identify them and perform antimicrobial sensitivity testing.

### Department Information

UNIVERSITY HEALTH NETWORK/Mount Sinai Hospital Department of Microbiology

14<sup>th</sup> Floor - Mount Sinai Hospital  
600 University Ave.  
Toronto, Ontario  
(416) 586-4432

### Regular Lab Hours

7:45am to 16:00pm **or** 8:00am to 16:15pm

Coffee and lunch schedules should be worked out with the bench tech you are working with.

If you are sick or cannot report to the lab please call the Senior Technologist's at (416)586-4800 ext 3822 and leave a message.



---

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 3 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Schedule for MLT Students Rotation

### 1. Orientation

- i. Introduction to Microbiologist-in-Chief Dr. Mazzulli and Microbiologist Dr. Poutanen. Introduction to Administrative Director - Katherine Wong, Supervisor - Pauline Lo and Charge Technologists - John Ng and Lilian Law.
- ii. Tour of lab
  - where to get clean lab coats, put dirty ones and wash up area
  - various benches, incubators, fridges, store room
  - bathroom and lunchroom facilities, including Food Hall on main floor
  - bulletin board with schedules, notices, lectures
  - conference room where journals and reference books are kept
  - where main lab manual, CLSI guidelines, Safety Manual, WHMIS Manual and other reference books are kept
- iii. Safety talk by Safety Officer
  - overview of lab safety procedures e.g. no open toe or open heel shoes, fire alarm pulls, stairway exits, fire extinguishers, eye wash / shower stations, WHMIS, HPTR etc



### 2. Competencies and Objectives

Please refer to objectives given to you by your individual teaching institutions. Objectives are also available with the senior teaching technologist.

**Note:** Competencies and Objectives are actually set by CSMLS.

### 3. Training

We have students from different institutions with variations in their length of rotation therefore, students will rotate differently on the benches. Some will be doing all and some will be doing specific benches. The benches include Specimen Processing, Urines, Respiratory, Enterics/Gynae., Wounds/Miscellaneous, Bloods, Infection Control, Quality Control, Virology/Serology and perhaps mycology. A schedule will be given to you on your first day.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 4 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

#### 4. Performance Assessments

Each student will be evaluated by the technologist working the bench that the student is assigned to that week. Students will also evaluate the benches on their rotation every week. Final evaluation will be done by the Senior Teaching Technologist at the end of rotation. Each student (except Mohawk) will be given "unknowns" throughout the microbiology rotation.

St Lawrence College students do weekly quizzes and final exam online. These students will also have a practical exam week 8.

Michener students do weekly quizzes online and submit competencies for sign-off via CompTracker.



UOIT students do weekly quizzes online. These students will also have a written and practical exam week 8. These students submit competencies for sign-off via CompTracker.

Mohawk students will have a quiz given by the teaching technologist at the end of Week 3.

#### 5. Meetings and Rounds

There are several tutorials/rounds you may attend:



- Departmental meetings and/or seminars are every Wednesday from 2:00-3:00 p.m. They are held in the 14<sup>th</sup> floor or 11<sup>th</sup> floor classroom.
- Other relevant seminars or infectious disease/micro rounds will be announced on an individual basis.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 5 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Bench Technologist – Teaching Guidelines

The bench technologist teaching the students should:

1. Discuss the nature of the specimens and the organisms (both pathogens and commensals) likely to be found in them. The student should be told how to identify these organisms and how the media is used for this purpose.
2. Make sure the organism(s) listed in each section are seen and worked up. If the organisms have not been seen in actual specimens, then the students should subculture the stock cultures and do the appropriate tests if time permits.
3. Demonstrate all procedures used in identification of pathogens on the bench, guide the students to perform all the procedures on their own. Due to the use of Vitek MS, students should be made familiar with classical methods as well.
4. Allow the students to work on the specimens independently for some (or all) of the time on the bench.
5. Read specimen gram smears. Correlate grams smear with culture results. The respiratory bench should show the students how to screen sputum smears by interpreting oropharyngeal contamination using epithelial cells in gram stains. The gynecological bench should show students bacterial vaginosis slides and also wet preps.
6. Complete bench evaluation on each student they work with.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 6 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Canadian Society for Medical Laboratory Science (CSMLS) - Code of Professional Conduct



1. Medical laboratory professionals are dedicated to serving the health-care needs of the public. The welfare of the patient and respect for the dignity of the individual shall be paramount at all times.
2. Medical laboratory professionals work with other health care professionals, to provide effective patient care.
3. Medical laboratory professionals shall promote the image and status of their profession by maintaining high standards in their professional practice and through active support of their professional bodies.
4. Medical laboratory professionals shall protect the confidentiality of all patient information.
5. Medical laboratory professionals shall take responsibility for their professional acts.
6. Medical laboratory professionals shall practice within the scope of their professional competence.
7. Medical laboratory professionals shall endeavour to maintain and improve their skills and knowledge and keep current with scientific advances. They will uphold academic integrity in all matters of professional certification and continuing education.
8. Medical laboratory professionals shall share their knowledge with colleagues and promote learning.
9. Medical laboratory professionals shall be aware of the laws and regulations governing medical laboratory technology and shall apply them in the practice of their profession.
10. Medical laboratory professionals shall practice safe work procedures at all times to ensure the safety of patients and co-workers and the protection of the environment.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 7 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Microbiology Rotation for Residents

### Welcome to Microbiology!

Microbiology is an interesting and dynamic science. Your time spent in the Microbiology Laboratory should be an enjoyable and successful learning experience.

During your rotation you will be exposed to all aspects of the department. You will handle specimens from their receipt in the department through to the reporting of results. You will develop the ability to make a presumptive identification of common pathogens and be able to select the correct tests to fully identify them.

In order to make your rotation a positive experience for both you and the other members of the department following these guidelines will be useful:

1. Your rotation in the clinical laboratory is first priority.
2. Don't change or delay your rotation without first discussing this with Dr. Poutanen.

### Department Information

UNIVERSITY HEALTH NETWORK/Mount Sinai Hospital Department of Microbiology

14<sup>th</sup> Floor - Mount Sinai Hospital  
600 University Ave.  
Toronto, Ontario  
(416) 586-4432

### Regular Lab Hours

7:45am to 16:00pm **or** 8:00am to 16:15pm

Coffee and lunch schedules should be worked out with the bench tech you are working with.



If you are sick or cannot report to the lab please email Dr. Poutanen or "[microbiologyspecialqueries@sinaihealthsystem.ca](mailto:microbiologyspecialqueries@sinaihealthsystem.ca)".

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

---

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 8 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Schedule for Resident's Microbiology Rotation

### Orientation



1. Introduction to Dr's Mazzulli, Poutanen, and McGeer. Introduction to Administrative Director Katherine Wong, Supervisor Pauline Lo and Charge Technologist Lilian Law & John Ng.
2. Tour of lab:
  - where to get clean lab coats, put dirty ones and wash up area
  - various benches, incubators, fridges, store room
  - bathroom and lunchroom facilities, including cafeteria location
  - bulletin board with schedules, notices, lectures
  - library where journals and reference books are kept
  - where main lab manual is kept and CLSI guidelines Safety Manual, WHMIS Manual and other reference books (QA bench)
3. Full safety talk by Safety Officer - Judith Cunningham,
  - short overview of safety rules e.g. no open toe or open heel shoes, fire alarm pulls, stairway exits, fire extinguishers, eye wash / shower stations, re: WHMIS, HPTR etc.
4. Arrange time for Quality Management talk by Quality Assurance technologist.
5. Arrange time for IC talk by an Infection Control Educator.

### Training

One to two weeks on the following benches – Specimen Processing, Urines, Respiratory, Enterics/Gynae., Wounds/Miscellaneous, Bloods, Infection Control, Quality Control and Virology/Serology. Extra weeks in certain areas if desired and time permitted.

A schedule will be given to you on your first day.



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 9 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### **Additional Training**

One month with Infection Control Department.  
 One week at Public Health Lab (Optional)

### **Performance Assessments**

1. Bench evaluations will be done by both the bench technologist and the student weekly.
2. Unknowns: series of unknown cultures to be given on the second week of the bench rotation.
3. Weekly Teaching Sessions / Questions
  - one per week that will cover the area of the laboratory that you have been exposed to.
  - Consists of review and discussion with staff microbiologists.
4. Examinations (For Microbiology Residents Only)
  - One practical exam will be given consisting of specimens to work up, and smears to read.
  - An oral exam will also be given by the Microbiologists



### **Meetings and Rounds**

#### **1. Infectious Diseases/Microbiology Rounds**

- 1300 - 1400 every Tuesday, at TGH Radiology Conference Room, East South, 1-452.
- see bulletin board for list of topics

#### **2. Departmental Seminars**

- 1400 - 1500 every Wednesday
- you may be required to present a topic and/or fascionoma rounds.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 10 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### 3. Infection Control Meetings

- 1430 - 1530 every Thursday

### 4. Lab Management Meetings

- every Wednesday

### 5. Quality Team Meetings



- every Monday

### 6. Other Rounds / Seminars

- Please refer to Academic Events Schedule

### Projects

- you may be given a project to work on with guidance from the Microbiologist
- plans for the project should be made halfway through the diagnostic rotation and media etc. ordered. (Please check with Pauline)

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 11 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

**University of Toronto Training Objectives**  
**Postgraduate Training Program in Infectious Diseases**

**Training Objectives for Adult Infectious Diseases: Medical Microbiology Rotation**

The department provides a six month rotation for infectious diseases trainees. The objectives of this rotation are:



1. To familiarize the resident with all microbiology services provided by the department including routine bacteriology, mycology, parasitology, virology, serology, and quality assurance.
2. To familiarize the resident with both the laboratory and clinical aspects of infection by having them rotate on the infection control bench as well as actively participating on the Infection Control Team, assisting with the Surveillance Program, taking part in the investigation of any relevant problems that may arise during their rotation, attend the weekly infection control team meetings and the monthly infection control team meetings.
3. To expose the resident to clinical and applied research by having them carry out a research project over the six months of their rotation. The candidate is expected to submit an abstract to a major national or international meeting and prepare a manuscript for submission for publication.
4. To allow the resident to attend ongoing microbiology and infectious diseases rounds/seminars as well as participate in, and present at, the departmental laboratory rounds and journal club.
5. To familiarize the resident with common clinical problems presented to the medical microbiologist by participating in a discussion of a pre-assigned weekly "question-of-the-week" regarding clinical microbiology and/or laboratory management issues.
6. To familiarize the resident with aspects of laboratory management and quality assurance by having the resident participate in management meetings and as member of the quality committee.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 12 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Resident Objectives at MSH

### SPECIMENS

**Reference:** Cumitech 9. Collection and Processing of Bacteriological Specimens

The resident shall:

- 1.0 Identify properly collected and transported specimens
- 2.0 Recognize and state problems in poor or unsuitable specimens
- 3.0 Demonstrate proper and safe handling of all specimens
- 4.0 Demonstrate proper use of the biological safety cabinet and how it works

### CULTURAL TECHNIQUES

The resident shall:



- 1.0 Select suitable primary media for all types of specimens
- 2.0 Recognize specimens that should be set up for anaerobes and understand the importance of handling anaerobic specimens promptly and planting them on pre-reduced media
- 3.0 Demonstrate proficiency inoculating and streaking a specimen on agar media
- 4.0 Be able to set up cultures in an anaerobic atmosphere and a 5-10% CO<sub>2</sub> atmosphere
- 5.0 Perform quantitative cultures and report the colony count
- 6.0 Prepare and read a simple wet mount for *Trichomonas vaginalis* specimen
- 7.0 Examine and report on the suitability of a sputum for culture (screening of a sputum)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 13 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## STERILIZATION

**Reference:** Disinfection, Sterilization, and Preservation, Block, Fourth Edition.

The resident shall:

- 1.0 Demonstrate aseptic technique
- 2.0 Select and perform correct sterilization for all types of media and objects
- 3.0 Demonstrate an understanding of these methods of sterilization and how to quality control them
  - 3.1 Dry heat
  - 3.2 Moist heat
  - 3.3 Filtration
  - 3.4 Ultra violet irradiation
  - 3.5 Ethylene oxide
  - 3.6 Gamma radiation
- 4.0 Demonstrate an understanding of disinfection
  - 4.1 Select proper disinfection for cleaning lab equipment and spills

## MEDIA

**Reference:** Media for Isolation - Cultivation - Identification Maintenance of Medical Bacteriology, McFadden

For media used in the laboratory the resident shall:



- 1.0 List the main ingredients and describe their functions
- 2.0 State the classification of the media
- 3.0 Quality control the medium by
  - 3.1 Checking sterility
  - 3.2 Checking pH and appearance
  - 3.3 Checking biochemical and growth of control organisms
  - 3.4 Checking media storage is at proper temperature and length of time

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

---

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 14 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## GRAM STAIN

For the Gram stain, the resident shall:

- 1.0 State principle of Gram's method
- 2.0 Perform a Gram stain
- 3.0 Recognize and state problems in poor or incorrectly stained slides
- 4.0 Recognize and report on Gram stained slides of specimens and cultures
- 5.0 Quantitatively report on direct smears

## MICROSCOPY

The resident shall:

- 1.0 State the principles of operation of light, phase contrast and UV microscope
- 2.0 Demonstrate proper use and care of the microscopes
- 3.0 Be able to perform Kohler illumination

## SAFETY / WHIMIS

The resident shall:



- 1.0 State the safety rules for a clinical laboratory
- 2.0 State procedure to be used for cleanup of various spills:
  - 2.1 Chemical
  - 2.2 Radioactive
  - 2.3 Blood
  - 2.4 Specimens
- 3.0 Demonstrate proper safety techniques for handling specimens

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 15 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- 4.0 State the principles of WHMIS
- 5.0 Demonstrate the proper handling of chemicals while using MSDS
- 6.0 Perform all procedures using chemicals using safety techniques as stated in the MSDS

**DIAGNOSTIC OBJECTIVES** (urine, genital, enteric, respiratory, wound and fluid benches).

- Reference:** Cumitechs:
- 4A. Laboratory Diagnosis of Gonorrhea
  - 14A. Laboratory Diagnosis of Central Nervous System Infections
  - 12A. Laboratory Diagnosis of Bacterial Diarrhea
12. Laboratory Diagnosis of Female Genital Tract Infections

Clinical Microbiology Procedures Handbook. Editor: Isenberg

For the following sites:

- a) respiratory tract (upper and lower)
- b) intestinal tract
- c) genitourinary tract
- d) skin
- e) urine
- f) blood
- g) CSF
- h) sterile fluids
- i) tissues

The resident shall:



- 1.0 State which specimens to be taken for the clinical picture of the patient.
- 2.0 State the specimens to be taken for proper lab work up.
- 3.0 State the proper collection methods including appropriate time for specimens to be taken, the amount and how to be transported to the laboratory.
- 4.0 State the organisms that commonly occur as pathogens, plus state what organisms may occur as commensal flora.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*



*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 16 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

5.0 Recognize and differentiate pathogens from commensal flora in the specimens.



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 17 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

For all the organisms listed below:

<i>Aerococcus urinae</i>	<i>Cryptococcus neoformans</i>
<i>Staphylococcus aureus</i>	<i>Haemophilus influenzae</i>
<i>Coagulase negative Staphylococcus</i>	<i>Haemophilus</i> species
<i>Staphylococcus saprophyticus</i>	<i>Campylobacter jejuni/coli</i>
<i>Streptococcus pyogenes</i>	<i>Campylobacter fetus</i>
<i>Streptococcus agalactiae</i>	<i>Listeria monocytogenes</i>
<i>Streptococcus pneumoniae</i>	<i>Corynebacterium pseudodiphtheriae</i>
<i>Viridans Streptococcus</i>	<i>Corynebacterium diphtheriae</i>
<i>Enterococcus</i> spp.	<i>Corynebacterium jeikium</i>
<i>Streptococcus anginosus</i> group	<i>Actinomyces israelii</i>
<i>Neisseria gonorrhoeae</i>	<i>Nocardia asteroides</i>
<i>Neisseria meningitidis</i>	<i>Clostridium difficile</i>
<i>Moraxella catarrhalis</i>	<i>Clostridium perfringens</i>
<i>Pseudomonas aeruginosa</i>	<i>Bacteroides fragilis</i> group
<i>Burkholderia cepacia</i>	<i>Prevotella melaninogenica</i>
<i>Acinetobacter</i> species	<i>Fusobacterium</i> species
<i>Pasteurella multocida</i>	<i>Propionibacterium</i> species
<i>Shigella</i> species	<i>Mycobacterium</i> species
<i>Salmonella typhi</i>	<i>Mycobacterium tuberculosis</i>
<i>Salmonella</i> spp.	<i>Candida albicans</i>
<i>Escherichia coli</i>	<i>Candida</i> species
<i>Yersinia enterocolitica</i>	
<i>Vibrio</i> spp.	
and all <i>Enterobacteriaceae</i>	

The resident shall:



- 1.0 Identify and record colonial appearance
- 2.0 Select optimal temperature, gaseous requirement and suitable media for isolation.  
Describe any special media that may be required, their main ingredients and functions.
- 3.0 Describe as applicable the typical cellular morphology, Gram stain reaction or wet mount appearance in direct smears and in culture.
- 4.0 Select and perform differential tests for identification.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 18 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- 5.0 Interpret all tests performed and identify the organisms.
- 6.0 State the principles, reagents and controls involved in all the differential tests for identification.
- 7.0 State the optimal method of susceptibility testing and the typical susceptibility pattern.
- 8.0 State the epidemiology, clinical significance and infection control implications.
- 9.0 State the problems of laboratory detection, identification and susceptibility testing.

## **BLOOD CULTURES**

**Reference:** Cumitech 1A. Blood cultures II

The resident shall:



- 1.0 State in what clinical situations blood cultures are to be taken.
- 2.0 State the proper collection methods.
- 3.0 For the blood cultures bottles
  - 3.1 list the major ingredients, and understand the purpose of each.
  - 3.2 state the recommended incubation time and temperature for all types of blood cultures.
  - 3.3 explain how and why blood cultures are subcultured.
- 4.0 For the Blood Culture machine in the laboratory
  - 4.1 state the principle of operation
  - 4.2 demonstrate proficiency in operating the machine
  - 4.3 perform regular routine maintenance
  - 4.4 recognize and solve minor problems
- 5.0 Demonstrate the proper procedures for working up positive blood cultures.
- 6.0 Describe the operation of the following blood culture system
  - 6.1 conventional

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY



---

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 19 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- 6.2 lysis centrifugation
- 6.3 infrared
- 6.4 radiometric
- 6.5 colorimetric detection (Virtuo)

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 20 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## ANAEROBES

**Reference:** Cumitech 5A. Practical Anaerobic Bacteriology



The resident shall:

- 1.0 Describe the principles of the anaerobic jar and anaerobic chamber, PRAS and rolled tube systems.
- 2.0 Operate the anaerobic jar utilizing proper techniques for achieving and monitoring anaerobiosis.
- 3.0 Differentiate organisms by their oxygen requirements.
- 4.0 State which specimens are suitable for anaerobic culture.
- 5.0 Describe the proper collection technique and transport for anaerobic specimens.
- 6.0 Discuss the type of media to be used for anaerobes.

## GENTALS

The resident shall:

- 1.0 Define bacterial vaginosis.
- 2.0 Differentiate between vaginosis, vaginitis or normal vaginal flora by Gram stain and wet preparations.
- 3.0 For *Trichomonas vaginalis* the resident shall:
  - 3.1 describe how to examine a wet mount.
  - 3.2 state it's main differentiating characteristics
  - 3.3 explain clinical significance

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 21 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## ANTIMICROBIAL SUSCEPTIBILITY TESTING

- Reference:** Cumitechs 25. Current Concepts and Approaches to Antimicrobial Agent Susceptibility Testing
- 6A. New Developments in Antimicrobial Agent Susceptibility Testing: A Practical Guide

### CLSI Guidelines



The resident shall:

- 1.0 Discuss all the antimicrobial agents used in the laboratory by:
  - 1.1 describing their antimicrobial spectrum
  - 1.2 stating the kind of activity they have (bactericidal, bacteriostatic)
  - 1.3 identifying the agents by their class
  
- 2.0 For each of the following antimicrobial susceptibility testing methods:
  - Kirby Bauer
  - MIC / MBC broth dilution
  - MIC agar dilution
  - MIC micro broth dilution
  - 2.1 state the principle
  - 2.2 perform and report results of the methods used in the laboratory
  - 2.3 select proper quality control procedures and indicate what troubleshooting would be done if the QC did not work
  - 2.4 explain limitations of each method
  
- 3.0 Describe the Beta Lactamase Test by:
  - 3.1 explaining the principle of the nitrocefin, acidometric, and iodometric methods
  - 3.2 performing and reporting the results
  - 3.3 stating the clinical significance of  $\beta$ -lactamase production

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 22 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## ANTIMICROBIAL SUSCEPTIBILITY TESTING

The resident shall:

- 4.0 Identify Methicillin Resistant *Staphylococcus aureus* (MRSA) by:
  - 4.1 explaining how MRSA's are detected
  - 4.2 stating the clinical significance of such an isolate
  - 4.3 discussing the problems in detecting MRSA
  
- 5.0 Identify Penicillin resistant Pneumococcus by:
  - 5.1 explaining the procedure used for screening and confirmation
  - 5.2 stating the criteria used to define resistance, relative resistance, and susceptible to Penicillin
  - 5.3 describing the clinical significance
  
- 6.0 Discuss the emerging new trends for susceptibility testing of *Enterococcus* by:
  - 6.1 explaining what synergy is and why it's important for *Enterococcus*
  - 6.2 determining synergy between Ampicillin and Gentamicin by the different methods available
  - 6.3 performing high level aminoglycoside testing and interpreting the results
  - 6.4 VRE
  
- 7.0 ESBLs and CRES

## MYCOLOGY

The student shall:



- 1.0 Identify properly collected and transported specimens for culture.
- 2.0 Perform and report on microscopic preparations from direct specimens by identifying fungal elements.
- 3.0 Select suitable primary media and atmospheric conditions for each specimen.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 23 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

Describe cultural requirements and differentiate the following microscopically and culturally:

*Microsporum*  
*Trichophyton*  
*Epidermophyton*  
*Aspergillus* species  
*Blastomyces dermatitidis*  
*Coccidioides immitis*  
*Histoplasma capsulatum*  
*Sporothrix schenckii*  
*Penicillium* species  
 Yeast (eg. Candida, Crypto)

5.0 Describe the appearance of *Malassezia furfur* in a direct mount of skin scales.

## SEROLOGY

**Reference:** Manual of Clinical Laboratory Immunology, Fourth Edition.

The resident shall:

- 1.0 Define antibody and antigen
- 2.0 State the principle, perform, read, report and interpret the results for the following tests:

Observe the following:

Architect i4000 reactive for serology tests

Latex Testing:           CMSE – CMV Antibody Capture Test  
                                   Monospot  
                                   Cryptococcal Antigen

EIA Tests:               Galactomannen



Vidas machine for VZ

DFA:                      PCP

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 24 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## VIROLOGY

PCR Tests: Chlamydia / NG (Roche 4800)  
 West Nile – NAT & Donor NAT testing (MPX)  
 Enterovirus/EBV/HSV/CMV/Parvovirus/VZ/  
 Adenovirus//HPV  
 Influenza A, B & RSV

### Reference:

- Schmidt N., Emmons R. Diagnostic Procedures for Viral, Rickettsial and Chlamydial Infections. 1989: American Public Health Association, New York. Sixth Edition.
- Wiedbrauk D., Johnston S. Manual of Clinical Virology. 1993: Raven Press, New York.
- Lennette E. (Ed.) Laboratory Diagnosis of Viral Infections. 1992: Marcel Dekker, New York. Second Edition.

The resident shall:

- Describe how a virus differs from bacteria.
- For the following virus families:
 

a) Herpesviridae	e) Retroviridae
b) Orthomyxoviridae	f) Togaviridae
c) Paramyxoviridae	g) Picornaviridae
d) Parvoviridae	h) Reoviridae



  - Classify according to:
    - RNA or DNA
    - Single or double stranded
    - Nucleocapsid symmetry
    - Positive or negative sense
    - Presence / absence of envelope
    - Size range
  - State the important human pathogens within these families.
- Explain what is meant by the following terms and give 2 examples for each:

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 25 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- a) Primary cell line
- b) Diploid cell line
- c) Continuous cell line

4.0 List the major components of:

- a) Viral Transport Media
- b) Growth Media
- c) Maintenance Media

5.0 For the following culture methods:

- a) Shell vial technique
- b) Tube culture technique

5.1 Briefly describe the method

5.2 List advantages and disadvantages for each method.

6.0. For the following viral illnesses:

- a) Herpes zoster
- b) Viral meningitis
- c) Systemic CMV disease
- d) Viral respiratory illness

6.1 List the appropriate specimens to be collected.

6.2 Describe appropriate specimen collection, transport and storage techniques.

6.3 Describe suitable inoculation and incubation techniques including selection of appropriate cell lines, culture methods and incubation.

6.4 State the virus(es) which may be expected to be detected.

6.5 Describe any rapid tests which may be used to aid in diagnosis.

7.0 For the viruses listed below:



- a) Herpes simplex type 1 and 2
- b) Cytomegalovirus
- c) Varicella Zoster
- d) Influenza A and B
- e) Parainfluenza 1, 2, and 3
- f) Respiratory syncytial virus
- g) Adenovirus

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 26 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- h) Enterovirus
- i) Coxsackievirus
- j) Poliovirus



- 7.1 State the optimal cell line(s) for isolation.
- 7.2 State the virus' effect on that cell line including a description of any cytopathic effect produced.
- 7.3 Describe the test(s) used for definitive identification once isolated, including any controls needed. (e.g. IFA stains)

8.0 Describe one method of anti-viral susceptibility testing.

### CMV QUANTITATIVE PCR

The resident shall:

- 1.0 Explain the significance of determining the presence of CMV antigen in the blood.
- 2.0 Show proficiency in CMV PCR testing by:
  - 2.1 Discussing proper specimen collection.
  - 2.2 Stating the principle of the procedure.
  - 2.3 Performing the procedure with quality control.
  - 2.4 Reporting results including quantitation of positive samples.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 27 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## **Mycology - Training Manual for Residents**

Teaching and Supervising Technologist: Mycology Technologist

Reference:

1. Davise H. Larone: Medically Important Fungi, A guide to identification, 3<sup>rd</sup>. Edition, ASM Press, 1995
2. Guy St – Germain, Richard Summerbell; Identifying Filamentous Fungi, A Clinical Laboratory Handbook, Star Publisher, 1996
3. Martha E Kern: Medical Mycology; Self Instructional Text, F A Davis, 1985 (1997)
4. Selected Reference Articles

**The Trainee at all levels of training shall know:**

### 1.0 *Processing of Specimens*

- 1.1 Know about proper mycology specimen collection, transportation and planting and the use of selective and/ or non-selective media
- 1.2 State the appropriate incubation temperature and length of incubation for recovery of various pathogens

### 2.0 *Direct Microscopy*

- 2.1 Demonstrate the proper of use of microscopy in the detection of fungal elements in clinical specimens using: Fungi-Fluor (FS), 10% KOH and Calcofluor White (CW)
- 2.2 Demonstrate and interpret fungal elements examined microscopically such as yeast, pseudo hyphae, septate or non-septate hyphae, narrow or broad base attachment etc. and/ or suspect or identify certain pathogens presumptively based on structure
- 2.3 Explain the different uses, advantages and disadvantages of the FS, CW and 10% KOH



### 3.0 *Procedure for identification of Yeast*

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 28 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	Subject Title: <b>Education Manual</b>	



- 3.1 Demonstrate the ability to identify various yeasts using Germ Tube, Oxgall, Cornmeal, Urea, EBM, Rapid Yeast ID System, API 20C procedures
- 3.2 Explain the principles, uses, advantages and disadvantages of the above methods
- 3.3 Ability to recognize organisms looking like yeast such as Prototheca based on structure
- 4.0 ***Protocol for reporting results as per manual policy***
  - 4.1 Reporting positive FS results in SOFT and phoning critical results to ward
  - 4.2 Reporting positive culture results in SOFT and phoning ward where applicable
- 5.0 ***Procedures for identification of moulds***
  - 5.1 Demonstrate the following methods used: Macroscopic Examination of Cultures, Microscopic Examination of Cultures, Scotch Tape Preparation, Tease Mount Preparation, Slide Culture, Sub-Culture, Phase Conversion, Permanent Mounts and Stock Preparation
  - 5.2 Explain the purpose of each of the above
- 6.0 ***Identification of the following fungi to the genus and/ or species level where indicated:***

Rapid Grower:

*Aspergillus fumigatus, niger, flavus and terreus*  
*Penicillium species and Paecilomyces species*  
*Fusarium species and Acremonium species*  
*Scopulariopsis species and Trichoderma species*  
*Scedosporium species (apiospermum and prolificans) and P. boydii*

Zygomycetes:

*Rhizopus species*  
*Mucor species*

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 29 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

*Absidia species*

Dematiaceous Moulds:

*Cladosporium species*

*Curvularia species* and *Bipolaris species*

*Alternaria species* and *Exophiala species*

Superficial Mycotic Agents:

*Dermatophytes* and *Scytalidium species*

*Trichosporon beigeli*

*Malassezia furfur* (Pityriasis versicolor – microscopic only, no culture necessary)

Dimorphic Fungi:

*Histoplasma capsulatum* and *Blastomyces dermatitidis*

*Sporothrix schenckii* and *Penicillium marneffe*



(*C. immitis* – Demonstration only)

7.0 ***Identification of Actinomycetes and other Filamentous Organisms***

*Nocardia species*

*Streptomyces species*



(Compare to Rapid Grower *Mycobacteria*)

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 30 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### General Microbiology Review

1. Review of common organisms should occur in the first week. This review will include:
  - #1 *Staphylococcus aureus*
  - #2 CNST
  - #3 *Enterococcus faecalis*
  - #4 *Escherichia coli*
  - #5 *Pseudomonas aeruginosa*
  - #6 Group B Streptococcus

(See Appendix I)
  
2. Kohler illumination (See Technical Section in Lab Manual)

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 31 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### **APPENDIX I - Day One**

1. On Day 1, subculture *S. aureus*, CNST, *Enterococcus* spp., *Group B Streptococcus*, *E. coli* and *Pseudomonas aeruginosa* onto MacConkey and Blood agar plates for isolated colonies. (See Appendix II).
2. On the second day gram stain each organism from any one plate. (See Appendix III). Also describe colony morphology for each organism on each of the plates.
3. On the second day also do the following tests and record results: (See Appendix IV).



Catalase - For either *Staph.* spp and *Enterococcus* spp.

Coagulase - For *S. aureus* and CNST

Bile esculin - For *Enterococcus* spp. and *Group B Streptococcus*

MUG - For *E. coli*

Oxidase - For *Pseudomonas aeruginosa* and *E. coli*

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 32 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## APPENDIX II- Media

1. MacConkey Agar
  - Is a selective and differential medium to isolate Gram negative bacilli
  - Lactose fermentors are pink
  - Non-lactose fermenters are pale
  - Peptone base with lactose and neutral red as indicator
2. Blood Agar
  - Typticase soy agar, Brucella agar or beef heart infusion base and 5% sheep blood
  - General purpose medium
  - Used to cultivate fastidious microorganisms
  - To determine haemolytic reactions
3. Chocolate Agar
  - Peptone base enriched with 2% haemoglobin
  - To cultivate *Haemophilus* and *Neisseria* species
  - Also a general purpose medium
4. Haemophilus Agar
  - Selective and differential medium
  - Chocolate base containing bacitracin to inhibit respiratory flora
  - To isolate *Haemophilus*
5. *Burkholderia cepacia* Agar
  - Selective agar for *B. cepacia*
6. CNA Agar
  - Columbia colistin nalidixic acid agar with 5% sheep blood
  - Selective and differential medium
  - Used to isolate gram positive organism from mixed cultures
7. Brucella Agar
  - With hemin and vitamin K
  - For the isolation and cultivation of anaerobic bacteria
8. LKV Agar
  - Kanamycin-vancomycin laked sheep blood



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***


*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 33 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

- Enriched selective and differential medium for the isolation and cultivation of anaerobic bacteria e.g. *Bacteroides* spp. and *Prevotella* spp.
9. Fastidious Anaerobic Broth
    - Broth culture for anaerobes
  10. Hektoen Enteric Agar
    - Selective and differential medium
    - Bile salts and indicator dyes select and inhibit gram positive organisms
    - Lactose, sucrose and salicin present to differentiate between fermenters and non-fermenters
    - Sodium thiosulfate and ferric ammoniumcitrate detect H<sub>2</sub>S production
    - *E. coli* are yellow-pink in colour.
    - *Shigella* are green or transparent Salmonella green or transparent with black centers.
  11. SMAC
    - MacConkey agar – with sorbitol
    - Selective and differential media used to isolate sorbitol – negative *E. coli*
    - Non-sorbitol fermenters appear colorless
  12. Selenite Broth
    - Enrichment broth for isolation of *Salmonella* and *Shigella*
    - Inhibits enterococci and coliforms
  13. CIN Medium
    - Cefsulodin-irgasan-novobiocin medium
    - Selective and differential for isolation of *Yersina enterocolitica*
  14. Campylobacter Blood Agar
    - Enriched selective blood agar
    - Trimethoprim, Vancomycin, Amphotericin B, Polymyxin B and Cephalothin in medium
  15. Martin-Lewis Agar
    - Enriched and selective medium for *N. gonorrhoeae*
    - Vancomycin, anisomycin, trimethoprim and colistin in medium



 Department of Microbiology Quality Manual	Policy # MI_EDU	Page 34 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

16. Strep B Carrot Broth
  - Selective enrichment and identification of beta-haemolytic Group B Streptococci
17. Denim Blue Oxoid Chromogenic Agar for MRSA
18. Brilliance VRE Chromogenic Agar

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY



*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 35 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	


### **APPENDIX III - Gram Stain Reactions**



- |                        |   |
|------------------------|---|
| Gram positive cocci    | <ul style="list-style-type: none"> <li>- <i>Staphylococcus aureus</i></li> <li>- Coagulase negative <i>Staphylococcus</i></li> <li>- <i>Enterococcus faecalis</i></li> <li>- Viridans streptococcus</li> <li>- Beta-haemolytic group B streptococcus</li> <li>- Peptostreptococcus species</li> </ul> |
| Gram positive bacilli  | <ul style="list-style-type: none"> <li>- <i>Clostridium novyi</i></li> <li>- <i>Bacillus</i> species</li> <li>- <i>Corynebacterium</i> species</li> </ul>   |
| Gram negative coccus   | <ul style="list-style-type: none"> <li>- <i>Moraxella catarhalis</i></li> <li>- <i>Neisseria</i> species</li> <li>- <i>Veillonella</i> species</li> </ul>   |
| Gram negative bacillus | <ul style="list-style-type: none"> <li>- <i>Escherichia coli</i></li> <li>- <i>Proteus mirabilis</i></li> <li>- <i>Pseudomonas aeruginosa</i></li> <li>- <i>Bacteroides</i> species</li> </ul>  |

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 36 of 56
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Education Manual</b>	

### APPENDIX IV - Tests

1. Catalase
  - bacteria that contain cytochrome enzymes are catalase positive and those that don't are catalase negative
  - *Staphylococcus* +
  - *Streptococcus* and *Enterococcus* -
  - flood bacteria with 3% hydrogen peroxide and observe for bubbles catalase
  - a)  $\text{H}_2\text{O}_2$  (3%)  $\rightarrow$  Catalase peroxide +  $\text{H}_2\text{O}$
  - b) Catalase peroxide +  $\text{H}_2\text{O}_2 \rightarrow$  Catalase  $\text{H}_2\text{O} + \text{O}_2$
  
2. Coagulase
  - the ability to clot plasma
  - two different coagulase tests can be performed, a tube test for free and bound coagulase and a slide test for bound coagulase, or clumping factor
  - *Staphylococcus aureus* +
  - Coagulase negative *Staphylococcus* -
  
3. Bile esculin
  - to detect beta glucoside which breaks down esculin to form a black precipitate due to the presence of ferric ions
  - *Enterococcus faecalis* +
  - Beta-haemolytic group B streptococcus -



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 37 of 56
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Education Manual</b>	

4. MUG
- this test is mainly done on the urine bench
  - to detect beta-glucuronidase
  - ortho-nitrophenyl-beta-glucuronide substrate is broken down by enzyme glucuronidase to produce ortho-nitrophenyl which produces a yellow color.
  - Indole production
- a) tryptophanase (produced by bacteria)  
 Tryptophan → ammonia + pyruvic acid + indole  
 ↓  
 reacts with p-dimethylamino benzaldehyde (Kovac's reagent)  
 ↓  
 quinoidal red violet
- b) Indole spot test uses Erlich's reagent 1% dimethylaminocinnamaldehyde  
 ↓  
 blue colour



*Escherichia coli* is MUG + and indole +

5. Oxidase
- to test for the production of oxidase
  - spot inoculate organism on to a filter paper soaked with 1% tetramethylphenylene diamine dihydrochloride - positive is purple, negative is yellow
  - *Pseudomonas aeruginosa* +
  - *Escherichia coli* -
- Oxidizing reaction  
 Reagent 1% Dimethyl or Tetramethyl para-phenylenediamine  
 ↓ on colonies  
 ↓  
 Indophenoloxidase (produced by bacteria)  
 ↓  
 ↓  
 Indophenol = black colonies (with dimethyl)  
 = magenta colonies (with tetramethyl)
6. Urease
- Detects urease production
  - Peptones in media utilized producing an alkalinity. Phenol red indicator
- Urea →→ ammonia

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY



*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 38 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

7. TSI
- Triple sugar iron agar. 0.1% glucose, 1% sucrose and lactose.
  - If the glucose is fermented, only a small amount of acid will be produced, which will be neutralized by alkali from peptone metabolism along surface of slant.
  - Oxidation of peptone cannot take place in the anaerobic conditions in the depth of the medium.
  - Therefore when glucose only is fermented, the butt of the medium becomes yellow and the slant remains red.
  - If lactose or sucrose is fermented, the amount of acid produced is large enough to offset alkali production and a yellow slant is produced.
  - Production of hydrogen sulphide is shown by formation of iron sulphide from the ferrous sulphate.
8. ONPG
- To detect enzyme  $\beta$ -D-galactosidase in lactose fermenting organisms. O-nitrophenol-B-D-galactopyranoside  $\rightarrow\rightarrow$  O-nitrophenol (yellow)
9. PPA
- To detect Phenylalanine deaminase production. Phenylalanine  $\rightarrow\rightarrow$  Phenylpyruvic acid +  $\text{FeCl}_3$  (ferric ions) = blue green
10. PYR
- To detect pyrrolidonyl peptidase
  - L-pyrrolidonyl- $\beta$ -naphthylamide  $\rightarrow$   $\beta$ -naphthylamine  $\rightarrow$  add cinnamaldehyde reagent = red colour

\*\*More extensive details for all tests performed in lab can be found in the Technical Manual.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 39 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## General Overview of Weekly Bench Rotations

### Specimen Management and Media

#### Learn types of media we use in lab and their purpose

Culture the following specimens:

- Urines
- Stools
- Cervical swabs
- Vaginal swabs
- Throat swabs
- Sputa (including screening)
- Bronchial washings
- Miscellaneous swabs
- Tissues
- Fluids, including CSF

Stain and read Gram smears.

Stain and read Kinyoun, Modified Kinyoun and Eosinophil smears..



Be familiar with operation of autoclave.

### Urine Bench

Read plates for colony counts and interpret significant colony counts.

Bacteria to have seen and identified:

- Proteus mirabilis*
- Escherichia coli*
- Klebsiella* species
- Staphylococcus saprophyticus*
- Enterococcus*
- Yeast
- Pseudomonas aeruginosa*

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 40 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### Stool Bench

Know biochemical reactions and specific media for:

*Salmonella* spp.  
*Shigella* spp.  
*Campylobacter* spp.  
*Yersinia enterocolitca*  
*Escherichia coli* 0157:H7  
*Vibrio* spp.

### Genital Bench

Read vaginal smears and wet preps. Should see lots of normal smears and vaginosis smears plus see yeast and *Trichomonas*.

Bacteria to have seen and identified:

*Neisseria gonorrhoeae*  
*Streptococcus agalactiae* (Group B strep)



### Respiratory Bench

#### Read sputum gram stains and know how to screen smears for oropharyngeal contamination.

Bacteria to have seen and identified:

*Streptococcus pyogenes* (Group A strep)  
*Pseudomonas aeruginosa*  
*Burkholderia cepacia*  
*Haemophilus influenzae*  
*Haemophilus parainfluenzae*  
*Streptococcus pneumoniae*  
*Moraxella catarrhalis*  
*Candida albicans*  
*Cryptococcus neoformans*



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 41 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## Miscellaneous & Sterile Fluids

### Read gram smears and correlate with cultures.

Bacteria to have seen and identified:

*Staphylococcus aureus* including MRSA  
*Serratia marcescens*  
*Enterococcus* species  
*Streptococcus anginosus* group  
*Streptococcus pyogenes*  
*Streptococcus pneumoniae*  
*Bacteroides* species  
*Prevotella* species  
*Clostridium perfringens*  
*Fusobacterium* species  
*Enterobacteriaceae*  
*Pseudomonas* species  
*Corynebacterium jeikeium*

## Blood Culture Bench

Operate the BacT Alert machine and assist with processing the specimens.  
Assist in working up the positives.

Bacteria to have seen and identified:


Staph. aureus  
*Neisseria meningitidis*  
*Listeria monocytogenes*  
*Corynebacterium* spp.  
*Viridans streptococcus*  
 Coagulase negative staphylococcus  
*Streptococcus pneumoniae*  
*Enterobacteriaceae*  
*Enterococcus* spp.  
 Anaerobes  
 Fastidious Gram negative bacilli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_EDU	Page 42 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	



Infection Control Bench

Learn about the different selective media used on the specific benches to identify MRSA, VRE and CPE.

Bacteria to have seen and identified:

- MRSA
- VRE
- ESBL
- CPE

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 43 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### Virology/Serology

Observe the following:

Accessioning and processing of specimens  
 Preparing smears for DFA (PCP)  
 Processing of FLU/RSV samples  
 Cryptococcal Latex Agglutination test

Architect Bench: theory , loading specimen, reflex tests, reporting

MPX Bench: Donor NAT testing, Monospot testing,  
 routine and STAT VZV Ab on Vidas

CMV PCR other PCR tests



Performing Chlamydia/GC/HPV PCR

### Mycology

Observe the following procedures:

Fungi Fluor stain for fungus

Isolation of various fungi and yeast, including knowledge of different media used.

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 44 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

## **Mycology - Training Manual for Students**

Teaching and Supervising Technologist: Mycology Technologist

Reference:

1. Davise H. Larone: Medically Important Fungi, A guide to identification, 3<sup>rd</sup>. Edition, ASM Press, 1995
2. Guy St – Germain, Richard Summerbell; Identifying Filamentous Fungi, A Clinical Laboratory Handbook, Star Publisher, 1996
3. Martha E Kern: Medical Mycology; Self Instructional Text, F A Davis, 1985 (1997)
4. Selected Reference Articles

**The student should know the following:**

### 1.0 *Processing of Specimens*

- 1.1 Know about proper mycology specimen collection, transportation and planting and the use of selective and/ or non-selective media
- 1.2 State the appropriate incubation temperature and length of incubation for recovery of various pathogens

### 2.0 *Direct Microscopy*

- 2.1 Demonstrate the proper of use of microscopy in the detection of fungal elements in clinical specimens using: Fungi-Fluor (FS), 10% KOH and Calcofluor White (CW)
- 2.2 Explain the different uses, advantages and disadvantages of the FS, CW and 10% KOH

### 3.0 *Procedure for identification of Yeast*



- 3.1 Demonstrate the ability to identify various yeasts using Germ Tube, Oxgall, Cornmeal, Urea, EBM, Rapid Yeast ID System, API 20C procedures
- 3.2 Explain the principles, uses, advantages and disadvantages of the above methods

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 45 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

#### 4.0 *Protocol for reporting results as per manual policy*

4.1 Reporting positive FS results in SOFT and phoning critical results to ward



4.2 Reporting positive culture results in SOFT and phoning ward where applicable

#### 5.0 *Procedures for identification of moulds*

Demonstrate the following methods used: Macroscopic Examination of Cultures, Microscopic Examination of Cultures, Scotch Tape Preparation, Tease Mount Preparation, Slide Culture and Sub-Culture. State the purpose of each step as listed above





  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 48 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

**6. Interpersonal relation:**

- |    |  |   |   |   |   |   |
|----|--|---|---|---|---|---|
| a. | Co-operation with others                         | 1 | 2 | 3 | 4 | 5 |
| b. | Ability to accept constructive criticism         | 1 | 2 | 3 | 4 | 5 |
| c. | Tactful in dealing with colleagues and personnel | 1 | 2 | 3 | 4 | 5 |



**7. Communication:**

- |    |   |   |   |   |   |   |
|----|---|---|---|---|---|---|
| a. | Ability to give clear and precise instructions and/ or explanations | 1 | 2 | 3 | 4 | 5 |
| b. | Good communication  | 1 | 2 | 3 | 4 | 5 |
| c. | Ability to ask for clarification in case of doubt                   | 1 | 2 | 3 | 4 | 5 |

**Comments:**

Supervising Technologist's signature



  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 49 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

**Bench Training Evaluation Form for Trainee - Mycology Section**

To be completed by the trainee at the end of training

Name of the trainee: \_\_\_\_\_



Dates of training: \_\_\_\_\_

Supervising Technologist:

- |  | Disagree<br>1 | 2 | Agree<br>3 | 4 | Strongly Agree<br>5 |
|--|---------------|---|------------|---|---------------------|
| 1. The objectives specified in the bench manual as applicable to your specific needs were covered during bench training.<br><b>If disagree please explain:</b>           |               |   |            |   | 1 2 3 4 5           |
| 2. The supervising technologist was patient and supportive (Approachable)  |               |   |            |   | 1 2 3 4 5           |
| 3. Frequent opportunities for questions and discussions were offered   |               |   |            |   | 1 2 3 4 5           |
| 4. The supervising technologist was able to give:  |               |   |            |   |                     |
| a. Clear and precise instructions  |               |   |            |   | 1 2 3 4 5           |
| b. Explanations  |               |   |            |   | 1 2 3 4 5           |
| 5. Feed back about your performance was given throughout the rotation  |               |   |            |   | 1 2 3 4 5           |
| 6. There were adequate hands-on participation on the bench   |               |   |            |   | 1 2 3 4 5           |
| 4. The supervising technologist's teaching abilities are:  |               |   |            |   |                     |
| <input type="checkbox"/> Poor <input type="checkbox"/> Fair <input type="checkbox"/> Average <input type="checkbox"/> Above Average <input type="checkbox"/> Outstanding |               |   |            |   |                     |

Comments:

Trainee's signature:

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 50 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

**Trainee Evaluation Form For Virology/Serology Section**

Name of the Trainee:

Dates of training:



	<b><u>Observed</u></b>
1. Virology Accessioning Bench:	<input type="checkbox"/>
a. Accessioning and processing of specimens	
b. Preparing smears for DFA (PCP)	
c. Run FLU/RSV samples	
d. Perform Cryptococcal Latex Agglutination test	
2. Architect Bench:	<input type="checkbox"/>
a. Theory of Architect testing	
b. Loading of specimens	
c. Reflex tests. (Patient vs. Donors)	
d. Printing of results and Reporting	
3. MPX Bench:	<input type="checkbox"/>
a. Donor NAT testing	
b. Perform Monospot testing	
c. Perform routine and STAT VZV Ab on Vidas	
4. PCR Benches :	<input type="checkbox"/>
a. CMV PCR	
b. other PCR tests	
5. Viral Load Bench:	<input type="checkbox"/>
a. Performing Chlamydia/GC/HPV PCR	

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 51 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	



### **Books in Microbiology Laboratory**

1. A guidebook to microscopical methods. A.V. Grimstone and R.J. Skaer. Cambridge University Press. 1972.
2. Anaerobe Laboratory Manual. 4<sup>th</sup> Edition. Lillian V. Holdeman et al. V.P.I. Anaerobe Lab. Blacksburg Virginia. 1977.
3. Antibiotics in Laboratory Medicine. 2<sup>nd</sup> Edition. Victor Lorian M.D. Williams and Wilkins. Baltimore U.S.A. 1986.
4. Approved Lists of Bacterial Names. V.B.D. Sherman et al. ASM Washington D.C. 1980.
5. Atlas of Human Parasitology. 2<sup>nd</sup> Edition. Lawrence R. Ash and Thomas C. Orihel. American Society of Clinical Pathologists Press, Chicago, U.S.A. 1984.
6. Bailey and Scott's Diagnostic Microbiology, 8<sup>th</sup> Edition, E.J. Baron and S.M. Finegold, C.V. Mosby Company, St. Louis, Missouri, 1990.
7. Basic Malaria Microscopy. Part 1. Learner's Guide. WHO, Geneva, 1991.
8. Bergey's Manual of Systematic Bacteriology. Vol. 1 Noel R. Krieg. 1984. Vol. 2 Peter H.A. Sneath 1986. Vol. 3 James T. Stanley 1989. Williams and Wilkins. Baltimore MD U.S.A.
9. Biochemical Tests for Identification of Medical Bacteria. 2<sup>nd</sup> Edition. Jean F. MacFaddin. Williams and Wilkins. Baltimore U.S.A. 1980.
10. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> Edition. Jean F. MacFaddin. Williams and Wilkins. Baltimore U.S.A. 2000.
11. Biosafety in Microbiological and Biochemical Laboratories. U.S. Dept. of Health and Human Services, OHS CDC. 2<sup>nd</sup> Edition. U.S. Gov. Washington, May 1988.
12. Clinical and Pathogenic Microbiology. Barbara J. Howard. Mosby Co. St. Louis Missouri. 1987.
13. Clinical Microbiology Procedures Handbook. Vol. 1 & 2. Isenberg ASM. 1992.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 52 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	



14. Clinics in Laboratory Medicine. Raymond D. Aller et al. Saunders Co. Philadelphia 1985.
15. Clinics in Laboratory Medicine. Blood Cultures, Vol. 14, No. 1. Michael L. Lulson et. al., W.B. Saunders Co. , March 1994.
16. Collection and Handling of Laboratory Specimens - A Practical Guide. Jean M. Slockbower et al. J.P. Lippincott Co. Philadelphia. 1983.
17. Color Atlas and Textbook of Diagnostic Microbiology. 2<sup>nd</sup> Edition. Elmer W. Koneman et al. J.B. Lippincott Co. Philadelphia Pennsylvania 1983.
18. Color Atlas and Textbook of Diagnostic Microbiology. 3<sup>rd</sup> Edition. Elmer W. Koneman et al. Lippincott. 1988.
19. Cowan and Steel's Manual for the Identification of Medical Bacteria. 3<sup>rd</sup> Edition. G.I. Barrow and R.K. A Feltham. Cambridge Univ. Press, Cambridge, G.B., 1993.
20. Disinfection Sterilization and Preservation. 4<sup>th</sup> Edition. Seymour S. Block Lea & Febiger. 1991.
21. Fundamentals of Anaerobic Bacteriology as Related to the Clinical Laboratory. ASM 1980.
22. Handbook of Laboratory Safety. 2<sup>nd</sup> Edition. Norman V. Steere. CRC Press Inc. 1971.
23. Identification of Enterobacteriaceae. 4<sup>th</sup> Edition. William H. Ewing. Elsevier Science Publishing Co. Inc. New York. 1986.
24. Identification of Medical Bacteria. 2<sup>nd</sup> Edition. Cowan and Steel. Cambridge University Press. 1979.
25. Interpretive Medical Microbiology. Harry P. Dalton et al. Churchill Livingstone. New York. 1986.
26. Laboratory Acquired Infections. C.H. Collins. Butterworth and Co. Ltd. 1983.
27. Laboratory Exercises in Microbiology. Robert A. Pollack et. al. John Wiley & Sons Ltd., New York, U.S.A. 2002.
28. Laboratory Safety CSLT Guidelines. 4<sup>th</sup> Edition. Gene Shematek. Wayne Wood. CSLT 1996.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 53 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	



29. Laboratory Safety: Principles and Practices. Brinton M. Miller et al. ASM Washington DC. 1986.
30. Laboratory Safety - Principles and Practices. 2<sup>nd</sup> Edition. Diane Fleming et al. ASM Press Washington. 1995.
31. Manual of Clinical Microbiology. 4<sup>th</sup> Edition. Edwin H. Lenette et al. ASM Wahington DC. 1985.
32. Manual of Clinical Microbiology. 5<sup>th</sup> Edition. Albert Balows et al. ASM Washington DC. 1991.
33. Manual of Clinical Microbiology. 6<sup>th</sup> Edition. Ellen Jo Baron et al. ASM Washington DC. 1995.
34. Manual of Clinical Microbiology. 7<sup>th</sup> Edition. Patrick R. Murray et al. ASM Washington DC. 1999.
35. Manual of Clinical Microbiology. 8<sup>th</sup> Edition. Patrick R. Murray et al. ASM Washington DC. 2003.
36. Manual of Clinical Microbiology. 9<sup>th</sup> Edition. Patrick R. Murray et al. ASM Washington DC. 2007.
37. Manual of Commercial Methods in Clinical Microbiology. Allan L. Truant. ASM Press. Washington, DC, U.S.A. 2002.
38. Medically Important Fungi. A Guide to Identification, 2<sup>nd</sup> Edition. Davire H. Laeone. Elsevier Science Publishing Co. Inc., N.Y., U.S.A. 1987.
39. Microscopy from the very beginning. Friedrich Karl Mollring. Carl Zeiss, Oberkochen. 1981.
40. Parasitic Diseases. Katz Michael et. al. Springer-Verlag, N.Y., 1984.
41. Practical Atlas for Bacterial Identification. D. Roy Cullimore. Lewis Publishers. 2000.
42. Practical Laboratory Mycology, 3<sup>rd</sup> Edition. Elmer W. Koneman and Glenn D. Roberts. Williams and Wilkins, Baltimore, U.S.A. 1985.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 54 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	



43. Productivity Power. Jim Temme. Shillpath Publishing, U.S.A. 1993.
44. Sexually Transmitted Diseases: An Illustrated Guide to Differential Diagnosis. Lewis H. Kaminester, M.D., Burroughs Wellcome Co., North Carolina, U.S.A., August 1991.
45. A Guide to Specimen Management in Clinical Microbiology. J. Michael Miller. ASM Press, Washington, DC, U.S.A. 1996.
46. Textbook of Diagnostic Microbiology. 3<sup>rd</sup> Edition. C.R. Mahon et al. Saunders, St. Louis Missouri. 2007.
47. The Laboratory Quality Assurance System. A Manual of Quality Procedures with Related Forms. Thomas a. Ratliff, JR., Van Nostrand Reinhold, N.Y., U.S.A. 1990.
48. The Compound Light Microscope. Dale Sinclair R.T. TIMT. 1986.
49. The Shorter Bergey's Manual of Determinate Bacteriology. 8<sup>th</sup> Edition. Williams and Wilkins. 1977.
50. Wadsworth Anaerobic Bacteriology Manual. 2<sup>nd</sup> Edition. Vera L. Sutton et al. The Regents of the University of California. 1975.
51. Wadsworth Anaerobic Bacteriology Manual. 3<sup>rd</sup> Edition. Vera L. Sutton et al. C.V. Mosby Co. St. Louis, Missouri. 1980.
52. Wadsworth Anaerobic Bacteriology Manual 4<sup>th</sup> Edition. Vera L. Sutton et al. Star Publishing Co. 1986.
53. Workplace Hazardous Materials Information System. (WHIMS): A Guide to the Legislation. Queens Printer for Ontario. May 1989.
54. Worthwhile facts about fluorescence microscopy. 3<sup>rd</sup> Edition. H.M. Holz. Carl Zeiss. D-7082 Oberkochen.
55. Zinsser Microbiology. Wolfgang K. Joklib et al. Appleton-Century Crofts N.Y. 1980.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 55 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

### Record of Edited Revisions

#### Manual Section Name: Education Manual



Page Number / Item	Date of Revision	Signature of Approval
Annual Review	May 26, 2004	Dr. T. Mazzulli
Annual Review	May 12, 2005	Dr. T. Mazzulli
Annual Review	July 5, 2006	Dr. T. Mazzulli
Annual Review	August 13, 2007	Dr. T. Mazzulli
Annual Review	August 12, 2008	Dr. T. Mazzulli
Annual Review	August 12, 2009	Dr. T. Mazzulli
Annual Review – Micro Rotation, CSMLS, Append.1, Books in Micro.	September 13, 2010	Dr. T. Mazzulli
Annual Review	March 30, 2011	Dr. T. Mazzulli
Pg.3 – Training - deleted “Mycology”	June 22, 2012	Dr. T. Mazzulli
Pg.4 – Revised Performance Assessments	June 22, 2012	Dr. T. Mazzulli
Pg.9 – Training – deleted “Mycology”	June 22, 2012	Dr. T. Mazzulli
Pg.10 – Revised Performance Assessments	June 22, 2012	Dr. T. Mazzulli
Pg.21 – Antimicrobial Sus. Testing – added CRES (7.0)	June 22, 2012	Dr. T. Mazzulli
Pg.21 – Serology changed AXSYM to Architect i2000	June 22, 2012	Dr. T. Mazzulli
Pg. 23 – Virology – Revised Observation	June 22, 2012	Dr. T. Mazzulli
Pg. 29 – General Micro Review – Added #6 – GBS	June 22, 2012	Dr. T. Mazzulli
Pg.33 – Appendix II(Media) - Revisions	June 22, 2012	Dr. T. Mazzulli
Annual Review	June 22, 2012	Dr. T. Mazzulli
Annual Review	August 19, 2013	Dr. T. Mazzulli
Pg. 3 & 9 – Remove Dr. Low & MacDonald	June 6, 2014	Dr. T. Mazzulli
Annual Review	June 6, 2014	Dr. T. Mazzulli
Update Headers	June 12, 2014	Dr. T. Mazzulli
Update Serology/Virology Tests & Competencies		
p.2 – Added perform AST to end of Welcome message	April 30, 2015	Dr. T. Mazzulli
p.5 – caveat: classical methods due to use of vitek MS - correlate grams read with culture results.	April 30, 2015	Dr. T. Mazzulli
p. 17 Residents: added to organism to know list: <i>Corynebacterium pseudodiphtheriae</i> <i>Aerococcus urinae</i> , <i>Vibrio</i> spp.	April 30, 2015	Dr. T. Mazzulli
p. 21 – changed to architrect i4000	April 30, 2015	Dr. T. Mazzulli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_EDU	Page 56 of 56
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Education Manual</b>	

p.30 - change ocbl to B.cepacia agar – selective medium	April 30, 2015	Dr. T. Mazzulli
p.37 – remove prepare a batch of medium	April 30, 2015	Dr. T. Mazzulli
p. 38/38 – add vibrio to stool, p.aer and B. cep to resp, GAS and pneumo to misc. pneumo to BC	April 30, 2015	Dr. T. Mazzulli
Removed RPR VD in virology sections	April 30, 2015	Dr. T. Mazzulli
Updated bench evaluations forms added	April 30, 2015	Dr. T. Mazzulli
Updates to: -Evaluation Form for Specimen Processing -Bench Evaluation for Student -Bench Evaluation Forms for Technologist	January 14, 2016	Dr. T. Mazzulli
Annual Review Updated MSH logo in header	April 30, 2016	Dr. T. Mazzulli
Update Simplexa to “Influenza” testing	September 7, 2016	Dr. T. Mazzulli
Updated - Trainee Evaluation Form For Virology/Serology Section (p.51)	November 2, 2016	Dr. T. Mazzulli
Annual Review	April 18, 2018	Dr. T. Mazzulli
Minor formatting change	September 14, 2018	Dr. T. Mazzulli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\