



<div><div>Toronto General Toronto Western Princess Margaret Toronto St. Luke's</div><div>Mount Sinai Hospital Sinai Health System Hospital of Mount Sinai Sinai Samaritan Sinai St. Elizabeth</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 1 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	
Prepared by QA Committee		
Issued by: Laboratory Manager	Revision Date: 1/27/2025	
Approved by Laboratory Director: Microbiologist-in-Chief	Next Review Date:	

**Uncontrolled When Printed**

## TABLE OF CONTENTS


BACTERIOLOGY SEND OUTS AND PROCEDURES .....	5
SoftMic Test Send out to PHOL .....	5
Special Instruction Specimens: .....	5
<i>Streptococcus pneumoniae</i> .....	5
TB Send out: .....	6
QuantiFERON® TB (QFT) Gold Plus: .....	6
<i>H. pylori</i> Send-out:.....	8
Parasitology Send out: .....	9
Legionella & Mycoplasma Send out.....	9
Mycology Send out .....	10
BACTERIOLOGY MOLECULAR DETECTION SENDOUTS: .....	11
Mycobacterium Molecular Detection (includes AMTD, TB PCR, Mycobacterium PCR) .....	11
Bacteria Molecular Detection (includes 16sPCR, S.aureus PCR, etc.) .....	12
Fungus Molecular Detection (includes Fungal PCR) .....	12
MetaGenomic Analysis.....	13
SEROLOGY SENDOUTS: .....	14
SoftLab Test Send outs to PHOL .....	14
CHECKING ORDERS FOR IGM PHL SENDOUTS ( SENIORS/ DESIGNATE).....	14
Prenatal Screening (MSH): .....	16
Prenatal Screening (UHN): Without Requisitions .....	16
NICU HIV Testing .....	17
Antistreptolysin O Titer (ASOT) Testing.....	17

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY " in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 2 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

Fertility and Reproductive Health Requisitions.....	18
Zika Virus Serology to PHOL.....	18
SIMIAN B VIRUS PROCESSING .....	19
Strongyloides Serology .....	21
VIROLOGY MOLECULAR DETECTION SEND OUTS .....	21
Hepatitis B DNA and Hepatitis C RNA Send out.....	21
<b>HIV VIRAL LOAD FOR MSH patients only</b> .....	22
TWH Liver Clinics .....	24
Measles/Mumps Diagnostic IgM/IgG Serology or PCR Send out .....	24
Monkeypox Sendout .....	25
Middle Eastern Respiratory Syndrome Coronavirus – MERS-CoV / AVIAN FLU PCR Send Out.....	25
WNV Molecular Testing to PHOL .....	26
Zika Virus Molecular Testing to PHOL .....	26
OTHER MISCELLANEOUS NATIONAL (NON-GTA) AND INTERNATIONAL SEND OUT SPECIMENS .....	27
COURIERS.....	28
Life Labs Courier Specimen Transportation.....	28
PHOL Courier Specimen Transportation.....	29
FEDEX Courier Specimen Send Out.....	29
KJV Courier Specimen Transportation.....	33
TAXI Specimen Transportation .....	35
Related Documents: .....	36
APPENDIX I: SEROLOGY TESTS ORDERING SCREENS .....	37
Table 1: serology in house from 1-R.....	37
Table 2: serology-prenatal.....	38

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\




 Department of Microbiology Quality Manual	Policy # MI_SO	Page 3 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

Table3: serology-donor .....	39
Tables 4, 5 & 6: serology-refer out .....	40
APPENDIX II: VIROLOGY/SEROLOGY PHOL PENDING LIST .....	42
APPENDIX III: UHN MANAGED MICROBIOLOGY SEND OUT TESTS .....	54
Record of Edited Revisions .....	55

For handling procedure for Creutzfeldt-Jakob Disease (CJD), see [Suspect RG 3 4 manual](#)

<div></div> <div>Department of Microbiology</div> <div>Quality Manual</div>	Policy # MI_SO	Page 4 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	


## REFERRAL LABORATORY TEST DICTIONARY

**Specimens can be referred out to Reference Laboratories for tests not processed in-house. See Referral laboratory Test Guide website for test availability and information.**

Referral Laboratory	Test Guide Information
Public Health Ontario Laboratory	<a href="https://www.publichealthontario.ca/en/ServicesAndTools/LaboratoryServices/Pages/Index.aspx">https://www.publichealthontario.ca/en/ServicesAndTools/LaboratoryServices/Pages/Index.aspx</a>
National Microbiology Laboratory	<a href="https://www.cnphi-rcrsp.ca/gts/faces/public/index.xhtml?lang=en">https://www.cnphi-rcrsp.ca/gts/faces/public/index.xhtml?lang=en</a>
MiraVista	<a href="http://miravistalabs.com/clinical-resources/forms/">http://miravistalabs.com/clinical-resources/forms/</a>
Center for Disease Control	<a href="http://www.cdc.gov/laboratory/specimen-submission/list.html">http://www.cdc.gov/laboratory/specimen-submission/list.html</a>
National Reference Centre for Parasitology	<a href="https://www.mcgill.ca/tropmed/nrcp">https://www.mcgill.ca/tropmed/nrcp</a>
Mitogen Advance Diagnostics	<a href="http://mitogen.ca/drug-testing/">http://mitogen.ca/drug-testing/</a> NMDA(NR1) Receptor Antibodies
Maisonneuve-rosemont Hospital	Voriconazole <a href="https://est.omni-assistant.net/hmr-labo/AnalysisRegistry/AnalysisSimpleSearch.aspx">https://est.omni-assistant.net/hmr-labo/AnalysisRegistry/AnalysisSimpleSearch.aspx</a> testing <a href="https://est.omni-assistant.net/hmr-labo/AnalysisRegistry/AnalysisSimpleSearch.aspx">https://est.omni-assistant.net/hmr-labo/AnalysisRegistry/AnalysisSimpleSearch.aspx</a>
Laboratory of Parasitic Diseases National Institute of Health	Filaria Serology
National Hansen's Disease (Leprosy) Clinical Center	Slit Skin Smear for Leprosy
Mayo Clinic Medical Laboratories	Gastric Biopsies for <i>H. pylori</i> <a href="https://www.mayocliniclabs.com/">https://www.mayocliniclabs.com/</a>

See [Fedex Send out Table](#) for a summary of key instructions for referred out tests with special directives.

See Appendix III for Microbiology samples sent out directly by UHN for UHN samples.

 <div>Department of Microbiology</div> <div>Quality Manual</div>	Policy # MI_SO	Page 5 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

## **BACTERIOLOGY SEND OUTS AND PROCEDURES**

### **SoftMic Test Send out to PHOL:**

1. Open **'Results Entry'** and scan order #.  
 -For Virology Send outs: select **PHLV** from the keypad. Within the comment type the test the sample is being sent to PHL for.  
 -For Bacteriology Isolate Send outs: in the isolate comment field for the corresponding isolate being send out, from the keypad select ">PRELI" and the comment corresponding to what you are sending to PHOL (eg. ID/AST to PHL).
2. **ctrl 'I'** for interim
3. Under Media "PHL" document the date the specimen was sent to PHOL.
4. To print PHOL report press instant report (from the top toolbar) , choose the correct printer, and format should be **MIC10A**, press ok
5. Refer to [PHOL Courier](#) section for specimen packaging instruction.

***Note:** When possible (when sufficient quantity is received), especially when only referred out test is requested, a portion of sample should be maintained in house for tissue, sterile sites etc.*

### **Special Instruction Specimens:**

#### *Streptococcus pneumoniae*

*Streptococcus pneumoniae* sent to PHOL for Penicillin confirmation requires a specific protocol be followed for successful testing to be completed.


1. Monday - Friday, subculture must be ready/packaged to send **before noon**
2. Subculture requirements for same day set up at PHOL:
  - Subculture must be **18-24 hours old** and **pure**
  - Plate must indicate **date of subculture**
3. Send out package must indicate:
  - For Microbiology Ref ID and Susceptibility section ,  
 661 University Avenue, Toronto, On M5G 1M1  
 Attention Deirdre Soares

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 6 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	


4. Bench technologist to phone PHOL next day for result @11:00
  - Tel # 416-235-6556
  - Bench to call Saturday morning. PHOL will have a technologist Saturday morning to read isolates sent Friday.
5. After 12:00 (noon) pick up on Friday or on the weekend:
  - PHOL does not **set up** *S. pneumoniae* susceptibility on the weekend.
  - Keep subbing isolate for send out on Monday morning.

### **TB Send out:**

1. Load the printer with the White Paper marked [Mycobacterium culture PHOL TB Form](#)
  2. Logon to Softmic
  3. Click '**Results**' tab from the left screen side , open **QC/Send out List** and choose **TB Culture New**
  4. Select **Mark (F8)**, scan specimen barcode, should have red check mark, **Print the PHOL form:** press **Instant Report** on the upper right side of the screen. choose the correct printer, and format should be **MIC12**, press ok
  5. Match names & LIS number, and pack sample with the PHOL report (no need to attach the LIS barcode). Follow the same for all samples
  6. **MACRO Reporting:** All samples scanned should have a red check mark on the left side.
  7. Choose Macro for **TB culture new**. Macro will result all specimen as “**@PHOL....**” and interim both **FA, TBS/TBA/TB**. All samples should have a **green!** after Macro is done.  
Print out worklist (CTRL+P) & file the work list in the BACTI send out binder
  8. Pack samples into blue PHOL bag and put into walk in fridge
- \*\*Keep blood sample for TB at room temp\*\***
- \*\* Put the type of specimen and site if not indicated\*\***

### **QuantiFERON® TB (QFT) Gold Plus:**

- MSH-out patient, testing not offered - the patient should directly go to Dynacare lab for sample collection and testing.
- MSH-in patient, we refer sample to [UHN Core Lab for 24/7](#)
- UHN-in and out patient, to TGH core lab, contact TGH specimen processing manager.

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 7 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

**Collection Tubes:** Collection tubes are available at Microbiology Laboratory (which are currently obtained from UHN Core Lab Inventory room)

There are 4 QFT Gold Plus blood collection tubes (0.8 – 1.2 mL in each tube):

1. **Grey Cap: NIL Tube**
2. **Green Cap: TB1 Antigen Tube**
3. **Yellow Cap: TB2 Antigen Tube**
4. **Purple Cap: Mitogen Tube**

#### **Specimen Collection Instructions all sites:**

1. Blood collection tubes should be at room temperature at time of collection.
2. For each patient, collect 1 mL of blood by venipuncture directly into each of the 4 QFT-Plus Blood Collection Tubes.

#### **Please ensure the following:**


- a. As 1 mL tubes draw blood relatively slowly, keep the tube on the needle for 2-3 seconds once the tube appears to have completed filling.
- b. The black mark on the side of the tubes indicates the validated range of 0.8 – 1.2 mL.  
**Important:** Collection volume should be within the black mark on the side of the tube. **If the level of blood in any tube is outside of the indicator mark, a new blood sample should be obtained.** Under or over-filling of the tubes outside of black mark range of 0.8 – 1.2 mL may lead to erroneous results.
- c. If a "butterfly needle" is being used, a "purge" tube should be used to ensure that the tubing is filled with blood prior to the QFT-Plus tubes being used.

3. Immediately after filling the tubes, shake them 10 times just firmly enough to ensure the entire surface of the tube is coated with blood, to solubilize antigens on tube walls.

**Note:** Overly vigorous shaking may cause gel disruption and could lead to aberrant results.

#### **Sample Transport**

- Send the collected samples at room temperature as soon as possible to Microbiology Laboratory, which will refer/send them out to TGH Core Lab.
- Incubation at TGH core-lab must start within 16 hours of collection.

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 8 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

### At TGH Core Lab:

1. Re-mix the tubes by inverting 10 times immediately prior to incubation at  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ . The specimen rack should sit off the floor of the incubator and the tubes should be upright.

#### Note:

- a. Each of the four tubes must be transferred to a  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$  incubator as soon as possible and within 16 hours of collection, and incubate for 16-24 hours.
  - b. If samples are not incubated for the set time or if samples are incubated longer than recommended, then aberrant results can be obtained.
2. Following Incubation: Centrifuge each of the 4 incubated collection tubes for 15 minutes at 2000-3000 RCF (g).
  3. Store samples in Refrigerator #22 Biochemistry Sample Drop-off fridge. Samples will be stable for 28 days at  $2-8^{\circ}\text{C}$  (refrigerated).

**Specimen Stability is 28 days for incubated, centrifuged & refrigerated samples. Specimen stability is up to 3 days following incubation at  $4-27^{\circ}\text{C}$ .**

### Reject Criteria

- Delayed arrival of specimen to TGH core lab (i.e. outside of 16hr from collection)
- Specimens not collected in the QuantiFERON®-TB Gold Plus blood Collection tubes
- Over filled or Under filled tubes (outside the black marking on the QFT collection container)
- Specimen received frozen or on ice

### H. pylori Send-out:



- Log in: MayoMedicalLaboratories.com
  - Click sign in (right top corner)
  - Email: MicrobiologySpecialQueries-MBX@siniahealthsystem.ca
  - Password: Microbiology\_\_ (change every 180 days)
- Under Test Catalog (left top corner) enter ZMMLS, click GO
- Click ZMMLS-Antimicrobial Susceptibility, Aerobic Bacteria, MIC
  - Click Order This Test (green box, right top corner)
  - Enter all the information under:

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is **Uncontrolled When Printed**.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 9 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

- Patient Information
- Ordering Physician Information
- Specimen Information
- Test-Specific Information: Gastric Biopsy-Helicobacter pylori
- Billing Information: Payer choose client
  - Term and Conditions: check mark box “I have reviewed and order and accept these terms”
  - Click Submit Now
- Print Forms
- Print Labels
- Sign out

### **Parasitology Send out:**

1. For Parasitology, load the printer with the White Paper [Stool Parasites PHOL Parasitology Form](#) Logon to Softmic
2. Click ‘Results’ tab from the left screen side, open **QC/Send out List** and choose **Parasitology New**
3. Follow 4-6 steps from above, please stick the diagnosis label on the requisition.
4. Choose Macro for **O&P**. Macro will result all specimen as “@PHOL....” and interim all parasites **O&P** send out. All samples should have a **green!** after Macro is done Print out & file the work list in the BACTI send out binder
5. Pack samples into blue PHOL bag and put into walk in fridge
6. **UHN samples:** To print the PHOL forms go EPR under Microbiology then press PHOL form, scan or type Auxiliary order, type Y for yes or A for accept to print out the PHOL form(has all comments & diagnosis). Please stick a small label for LIS order number.

**\*\* Put the type of specimen and SAF if not printed \*\***

**\*\*Filarial (Tissue Parasites Form) send out to PHOL is EDTA BLOOD \*\***

**\*\*Acanthamoeba corneal scraping (Arthropods & Miscellaneous Parasites Form) send out to PHOL\*\***

**\*\*For Leishmania PHOL form print from EPR if UHN patient.**

### **Legionella & Mycoplasma Send out:**



1. For **Legionella**, load the printer with the **White Paper** [Legionella Detection for PHOL Legionella Form](#)/ **White Paper** [Urogenital Mycoplasma Culture for PHOL Form](#).
2. Logon to Softmic
3. Click ‘Results’ tab from the left hand side, open **QC/Send out List** and choose **Send Out List**
4. Follow 4-6 steps from above.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*


Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>University Health Network Sciences General Sciences Women Sciences Children Sciences Mental Health</div><div>Mount Sinai Hospital Hôpital du Mont-Sinai Hospitall de Mont-Sinai Hospitall de Mont-Sinai</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 10 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

5. Choose Macro for **Legionella/Mycoplasma...** Macro will result all specimen as “**@PHOL....**” and interim all samples. All samples should have a **green!** after Macro is done Print out & file the work list in the BACTI send out binder
  6. Pack samples into blue PHOL bag and put into walk in fridge
- \*\* Put the type of specimen and site if not written \*\***
- \*\* Use Urogenital Mycoplasma Culture for all respiratory samples from newborn \*\***
- \*\* On each week Friday night pack all Urogenital Mycoplasma specimens in one brawn paper bag and store it in -70 freezers and on Monday or Tuesday morning the serology Technician 1 will take the bag out and pack it with DRY ICE then send to PHOL.**

#### **Mycology Send out:**

1. For nail/hair/skin scraping for fungus, load the printer with the **White Paper**  
[Mycology Form](#)
2. Logon to Softmic
3. Open ‘**Results Entry**’ and scan order #. Under CALCO and FUD choose }**PHOL** from keypad it will result specimen as “**@PHOL....**” and ctrl ‘**I**’ for interim
4. To print PHOL report press instant report (from the top toolbar) , choose the correct printer, and format should be **MIC12**, press ok
5. Match names & LIS number, and pack sample with the PHOL report(no need to attach the LIS barcode) follow the same for all samples
6. Pack samples into blue PHOL bag and keep in planting area at **room temperature**  
**\*\* Put the type of specimen and site if not printed\*\***

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 11 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

### **BACTERIOLOGY MOLECULAR DETECTION SENDOUTS:**

Molecular detection may be requested from a variety of specimens – all require permission prior to testing.

- 1) Physician is to contact the microbiology senior technologists to ensure there is enough sample remaining for testing.
- 2) If sufficient sample is available, a [UHN/MSH Request For Molecular Tests Form](#) will be emailed to the requesting physician to begin the approval process as outlined on the Physician Procedure and request form for Ordering Molecular.
- 3) The physician will then proceed to contact appropriate individuals for specific test permission as listed below to complete the Molecular Test Form.
- 4) Once completed laboratory staff will proceed to prepare specimen for testing as per instructions below.

#### **Mycobacterium Molecular Detection (includes AMTD, TB PCR, Mycobacterium PCR)**


- **Contact for permission:** Requesting physician to call PHOL Customer Service for approval. Request will be documented on [UHN/MSH Request For Molecular Tests Form](#) and forwarded to a senior or charge tech through email to [MicrobiologySpecialQueries-MBX@mtsinai.on.ca](mailto:MicrobiologySpecialQueries-MBX@mtsinai.on.ca)
- Lab orderables: MSH –use TBPCR (G), UHN –cannot order in EPR but we will add TBPCR (I) and enter results here when they come back.
- PHOL Requisition: Print PHOL requisition –Press instant report, choose appropriate printer, format should be MIC10. Match names and pack sample with PHOL requisition to send out. Make sure “Test requested” area is filled in with test “TB PCR”.
- Frozen tissue aliquots are acceptable if there is no more tissue sample left

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 12 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

### **Bacteria Molecular Detection (includes 16sPCR, S.aureus PCR, etc.)**

- **Contact for permission:** Requesting physician to call PHOL Customer Service for approval. Request will be documented on [UHN/MSH Request For Molecular Tests Form](#) and forwarded to a senior or charge tech through email to [MicrobiologySpecialQueries-MBX@mtsinai.on.ca](mailto:MicrobiologySpecialQueries-MBX@mtsinai.on.ca)
  - Lab orderables: MSH –use MPCRX (G), UHN –cannot order in EPR but we will add MPCR (I) and enter results here when they come back.
  - PHOL Requisition: Print PHOL requisition –Press instant report, choose appropriate printer, format should be MIC10. Match names and pack sample with PHOL requisition to send out. Make sure “Test requested” area is filled in with test.
- For approved Molecular Microbiology tests performed at Sick Kids hospital, complete the [Sick Kids Molecular Microbiology PCR order form](#). Pack and send sample through [KJV Courier](#).

### **Fungus Molecular Detection (includes Fungal PCR)**


- **Contact for permission:** Requesting physician should contact the on-call UHN/MSH Microbiologist. Request will be documented on [UHN/MSH Request For Molecular Tests Form](#) and forwarded to a senior or charge tech through email to [MicrobiologySpecialQueries-MBX@mtsinai.on.ca](mailto:MicrobiologySpecialQueries-MBX@mtsinai.on.ca)
- Lab orderables: MSH –use FPCR (G), UHN –cannot order in EPR but we will add FPCR (I) and enter results here when they come back.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

---

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\



 Department of Microbiology Quality Manual	Policy # MI_SO	Page 13 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

- Sick Kids requisition: Complete the [Sick Kids Molecular Microbiology PCR order form](#). Pack and send sample through [KJV Courier](#).

### MetaGenomic Analysis

- **Contact for permission:** Requesting physician should email Dr. Samir Patel ([Samir.patel@oahpp.ca](mailto:Samir.patel@oahpp.ca)) for permission. Request will be on [UHN/MSH Request For Molecular Tests Form](#) and forwarded to a senior or charge tech through email to [MicrobiologySpecialQueries-MBX@mtsinai.on.ca](mailto:MicrobiologySpecialQueries-MBX@mtsinai.on.ca)
- Lab orderables: MSH –use MGASX (G), UHN –cannot order in EPR but we will add MGAS (I) and enter results here when they come back.
- Send out with comment:  

"This specimen has been sent to Public Health Ontario Laboratories (PHOL) for metagenomic analysis (a research use only assay). Results will be communicated verbally directly from PHOL to the requesting physician, **Dr. [REDACTED]**. If you have any questions or would like to receive a verbal result, please contact Dr. Samir Patel from PHL at 647-792-3172."
- Requisition: PHOL Requisition: Print PHOL requisition –Press instant report, choose appropriate printer, format should be MIC10. Match names and pack sample with PHOL requisition to send out. Make sure “Test requested” area is filled in with test.

  <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 14 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

## **SEROLOGY SENDOUTS:**

### **SoftLab Test Send outs to PHOL**

1. Ensure sendout test is ordered, received and verified in the LIS.
2. Go to Print Scheduler. Choose the ALT RL10A REP = PHOL/OR HIV REPORT.
3. Page down to the Order# field and scan in or type the order number, choose the correct printer, press run, print out PHOL form
4. No need to put LIS Label on the PHOL form only exception is the UHN PHOL form attach LIS bar code
5. Write down the test requested (e.g. Hepatitis B surface antigen, Parvovirus B19 IgG/IgM, CMV IgG/IgM, VZ IgG.....) **under section-3 Test(s) Requested**  
**Do not use PREVIOUS PHOL CODE(S) V05, S17...**
6. Scan all samples into existing Task list – ‘9SERO’. Save & Print task list when all send outs have been entered, choose Layout as ‘TL-LIST’ from drop down menu, and file in “Serology PHOL Send out Binder”
  - a. If tasklist is already printed or if sending out a sample from another bench, place an LIS sticker on printed tasklist for current date.
7. Refer to [PHOL Courier](#) for sample preparation for send out.

### **CHECKING ORDERS FOR IGM PHL SENDOUTS ( SENIORS/ DESIGNATE)**

- 1) Order 8COM
- 2) Check Order Comments in Original Requisition Form / Order Comment Window for clinical information, UHN pts check diagnosis , history of ordering , last resort call doctor/ward if needed.

ACUTE FETAL INFECTION or? ACUTE FETAL INFECTION  
US (Ultrasound Findings)  
ACUTE, RECENT INFECTION or EXPOSURE


- 3) Document in 8COM  
Result Tab –IgM sent to PHL  
  \ M: Enter Clinical Information  
Also In- Patient: ordered in EPR or MSH or no Req. Form submitted

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 15 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

- 4) Match patient information , match Orders in SOFTLAB with PHL printed form ,
- 5) For PARVO Ig M, remove auto-comment “ Not sent out ” and “ Not tested @nci” replace with “\*NOTE\*@ MOH” , Technician will verify once sent out.
- 6) Add coloured label “Please do not cancel tests requested”
- 7) Scan all documents including PHOL printed form in SOFTMEDIA under the Order.
- 8) Give all forms back to Technicians for processing.


GROUP ORDERS Group Tests	Tests Included	Testing Lab
MMR Titers /Immune Status	Measles IgG Antibody	PHOL Serology
	Mumps IgG Antibody	PHOL Serology
	Rubella IgG Antibody	MSH Micro
MMRV Titers /Immune Status	Measles IgG Antibody	PHOL Serology
	Mumps IgG Antibody	PHOL Serology
	Rubella IgG Antibody	MSH Micro
	Varicella IgG Antibody	MSH Micro
TORCH TITER IgG PANEL TORCH Screening (Immune Status)	Toxoplasma gondii IgG Antibody	PHOL Serology
	Rubella Virus IgG Antibody	MSH Micro
	Cytomegalovirus (CMV) IgG Antibody	MSH Micro
	Herpes Simplex Virus (HSV) IgG Type 1 and Type 2 IgG Antibody	PHOL Serology
TORCH Screening IgG & IgM (Immune Status and Active Infection) Based on Clinical Information	Toxoplasma gondii IgG Antibody	PHOL Serology
	Rubella Virus IgG Antibody	MSH Micro
	Cytomegalovirus (CMV) IgG Antibody	MSH Micro
	Herpes Simplex Virus (HSV) IgG Type 1 and Type 2 IgG Antibody	PHOL Serology
	Toxoplasma gondii IgM Antibody	PHOL Serology
	Rubella Virus IgM Antibody	PHOL Serology
	Cytomegalovirus (CMV) IgM Antibody	MSH Micro
	Herpes Simplex Virus (HSV) IgM - TEST NOT AVAILABLE	NOT AVAILABLE

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


 Department of Microbiology Quality Manual	Policy # MI_SO Version: 4.2 CURRENT	Page 16 of 59
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: Send-out Information Manual</b>	

### **Prenatal Screening (MSH):**

1. Samples are stored in the Virology specimen fridge in a basket labeled Serology
2. Separate prenatal from fertility & Reproductive clinic (MIREF) samples
3. Logon to Soft Lab
4. Scan specimen barcode in “Order Entry” and compare test codes ordered in LIS with requisition (**9HIV: Prenatal HIV, 9HAG: Prenatal Hepatitis B surface Antigen, 9RUB: Prenatal Rubella IgG AB, 9VDP: Prenatal syphilis screening**)
5. Receive all tests under ‘Specimen’ tab but do not verify
6. If a HIV requisition is attached with the regular prenatal requisitions, remove HIV requisition and discard. **Do not separately send a sample for HIV**
7. **Check to make sure Dr. Mazzulli's name is stamped in (section-1 Clinician/Referring Laboratory) on Prenatal screening requisition**
8. **Write referring Dr's name from LIS/requisition under Clinician ID and “cc” ie: cc: Dr. Smith.**
9. Place LIS barcode label on requisition if not present and check the name matches on requisition and blood tubes
10. Once finished processing all samples, create a new task list – ‘9SERO’ and scan each sample
11. Save & Print task list when all send outs have been entered, choose Layout as ‘TL-LIST’ from drop down menu, and file in “PHOL Send out Binder”
12. **Arrange all samples in order in sample foam** and keep PHOL forms in order as well; pack both in a clear plastic zip lock. Put everything inside blue PHOL bags and store in walk-in fridge until morning pick-up

### **Prenatal Screening (UHN): Without Requisitions**

1. Logon to Soft Lab, scan specimen barcode in ‘AUX ORDER’
  2. Receive all tests under ‘Specimen’ tab but do not verify
  3. Attach LIS barcode label on each tube match names.
  4. Log in EPR, select Microbiology then select Prenatal report
  5. Scan or type the EPR Auxiliary number, type ‘A’ to accept and ‘Y’ to print out the PHOL form
  6. Attach LIS label to the prenatal form and match the names
  7. Please follow steps 7-12 above
- \*If additional send out tests are requested (other than the 4 specified on the prenatal requisition), print another PHOL form from LIS and add the extra tests**

 Department of Microbiology Quality Manual	Policy # MI_SO Version: 4.2 CURRENT	Page 17 of 59
<b>Section:</b> Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

## **NICU HIV Testing**

For Blood samples (whole blood in EDTA) received from the NICU for **HIV PCR**. These are **STAT** tests. Please process immediately upon receipt.

1. Test will be already ordered in Powerchart. If not, please add HIV PCR (9HIVP) and verify the send out. This is important, an auto email is programmed to send out once verified.
  - Serology is not necessary. If ordered, do not cancel. Send to PHOL with HIV PCR ordered as well. PHOL is aware PCR is priority if specimen quantity is an issue.
    - i. Please email Susan Guest (MSH) if the incorrect test has been ordered by the nurse
  - If CD4 is requested and two blood tubes are sent, you may forward one to core lab. If CD4 is requested and only 1 tube of blood was sent, please forward to PHOL for HIV PCR only.
2. Print Requisition using Print Scheduler “HIV PHL form” and write Baby HIV PCR on the requisition.
  - Bring sample promptly to a senior. The senior must fax a copy of the HIV requisition form to PHOL HIV/HTLV contact at 416-235-6194. The senior technologist will email [SupplementalVirologyLab@oahpp.ca](mailto:SupplementalVirologyLab@oahpp.ca) that a sample is coming.
3. Package specimen in its own brown paper bag labelled as “Baby HIV PCR”
  - Send **EDTA Whole Blood**
  - **Do NOT centrifuge** the blood
4. Specimen are to be sent with the morning PHOL shipment (10:30) Mon-Friday.
  - Specimen which misses the routine morning pick-up can be sent the next morning. From specimen collection, test must be performed within 5 days.

## **Antistreptolysin O Titer (ASOT) Testing**


1. Testing **REQUIRES** microbiologist approval. Ensure a microbiologist has approved before sending out.
2. Serum is required for ASOT testing.
3. Accession/Receive sample in LIS.
4. Prepare a [Hospital for Sick Kids serology requisition](#)
5. Send sample according to TGD regulation to Sick Kids by [KJV Courier](#).

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 18 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

### **Fertility and Reproductive Health Requisitions**

1. Samples are stored in Virology specimen fridge in a basket labeled Serology
2. Separate prenatal from fertility & Reproductive clinic (MIREF) samples
3. Logon to Soft Lab
4. Scan specimen barcode in 'Order Entry'
5. Re-order tests requested on requisition as 9 Serology codes in LIS (e.g. Parvovirus B19 IgG = 9PARG...)
6. Receive all tests under specimen tab
7. All tests should include the comment saying "The sample has been sent to the Ministry of Health Lab (PHOL)....." Add comment if not present: '\ 'F5' 'choose #1'
8. Verify all tests
9. **Print PHOL Form:** See [Printing Serology PHOL and HIV forms from Soft Lab](#)

### **Zika Virus Serology to PHOL**

For full information visit [PHOL Test Information Sheet website](#) for Zika Virus

Routine: Specimens: 2 Blood tubes required:  
1 serum tube (+ 1 extra Serum recommended)

Other Specimens: Test must be pre-approved by Dr. Tony Mazzulli  
5mL Urine  
400uL Amniotic fluid, tissue or CSF

#### **Required information needed for Zika Virus:**

The required information is mandatory upon ordering in EPR/CERNER and can be found in the LIS order comment when accessioning specimen.

The required information is as below:


- Relevant travel history (countries travelled to in the preceding 3 months)
- Dates of travel
- Symptom onset date
- Symptoms (e.g. fever, rash, arthralgia, conjunctivitis)
- Date of specimen collection
- If female, indicate if the patient is pregnant

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 19 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

- History of receiving any flavivirus vaccine (e.g. Japanese encephalitis vaccine, yellow fever vaccine)

#### Order and Shipment instructions to PHOL:

##### EPR (UHN)

- LIS order: 9ZIKX (Zika Virus IgM Antibody)  
ZIPCR (PCR for Zika Virus)
- PHOL Form: Print PHOL form from EPR; All clinical information will be printed.  
Send form with 2 serum specimens.

##### Cerner (MSH)

- LIS order: 9ZIKX (Zika Virus IgM Antibody)  
ZIPCX (PCR for Zika Virus)
- PHOL Form Print PHOL for using standard PHOL serology format.

##### Other:

- Blood and Urine specimens: store at 2-8°C, ship to PHOL on ice packs.
- Other specimens: store at 2-8°C, ship PHOL on ice packs
  - if delivery of "other specimens" to PHOL will take more than 72 hours, specimens should be frozen and shipped on dry ice.

#### Testing:

Serology TAT is 6-8 weeks

Molecular testing TAT up to 14 days

STAT testing not available

### **SIMIAN B VIRUS PROCESSING**

#### **TECHNICIAN:**

##### **SIMIAN B VIRUS BASELINE**

Sample taken from TGH New Employee in Animal House for Baseline


- Do not enter in LIS
- Enter Information in EXCEL File  
T:\microbiology\Virology\Simian B Virus Storage

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 20 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

<file:///T:/microbiology/Virology/Simian%20B%20Virus%20Storage/Simian%20B%20Virus%20Storage%202021.xls>

Enter Last Name, First Name, Date Received and Box #

3. Store serum in Simian Freezer MIT7 according to next available Freezer Box
4. Scan Requisition in V Drive /Microbiology /ScannedRequest/Planting  
File Name: example SimianBVirus20210428

### TECHNOLOGIST:

Simian B Virus Antibody Testing

1. Samples and requisition are prepared by OHS at TGH.
2. Ship samples on dry ice to

Dr. Julia K Hillard  
National B Virus Resource Center,  
Georgia State University  
161 Jesse Hill Jr. Drive, Atlanta GA 30303

3. With FEDEX Waybill , include CDC Permit to Import or Transfer Etiological Agents or Vectors of Human Disease .

[http://biotech.gsu.edu/virology/PDFs/DHS\\_Import\\_Documents\\_2017.pdf](http://biotech.gsu.edu/virology/PDFs/DHS_Import_Documents_2017.pdf)

4. On FEDEX on –line form, change the Bill transportation to “third party” enter UHN  
OHS Cost Centre  
OHS Cost Centre #711209015
5. Email/ phone National B Virus Resource to inform prior to shipping for them to schedule sample for testing  
Email: [bvirus@gsu.edu](mailto:bvirus@gsu.edu)  
Phone: 404.413.6550  
Fax: 404.413.6556
6. Examples of requisition, CDC form and FEDEX form and how to fill up the documentations can be found in send out book

Requisition Form:



<http://biotech.gsu.edu/virology/PDFs/2015%20Submission%20Form.pdf>

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 21 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

### Strongyloides Serology

1. Notify PHO before sending a priority specimen by contacting our Customer Service Centre at 416-235-6556 / 1-877-604-4567
2. Ship priority specimens separately from routine specimens in a clearly marked package (in accordance with the Transportation of Dangerous Goods Regulations) indicating "Priority" for easier retrieval.
3. We will use this bag to transport the samples with special stickers attached in each requisition form.

## VIROLOGY MOLECULAR DETECTION SEND OUTS

### Hepatitis B DNA and Hepatitis C RNA Send out



1. **Print out PHOL form:** Log on in EPR, press Microbiology then choose PHOL form, scan EPR#, press OK, type A for accept or Y for Yes
2. Log on to Soft Lab
3. Scan Auxiliary number into LIS Order entry
4. Click on 'Specimen' Tab, 'Receive all' print and label the blood purple tubes one for 9HBVD OR 9HCVQ and the other tube as 8ALI, attach a small barcode label to PHOL(EPR) form beside patient's name
5. Pack the tubes in the foam and send them to 6<sup>th</sup> floor
6. Pick up samples from 6<sup>th</sup> floor
7. Check if the sample is from **Liver Clinic** place the 2<sup>nd</sup> aliquot in designated rack in the fridge to be sent to TWH Liver Clinic the following day
8. If **not** from Liver Clinic place 2<sup>nd</sup> aliquot in freezer box for HBVD or HCVQ
9. Both Quantitative HCV RNA & Genotype can be written on the same PHOL Form send one tube to PHOL(double the volume)
10. Place 9HBVD OR 9HCVQ/9HCVG aliquots in designated rack in the freezer for the next morning to send to PHOL
11. The next morning Technician 2 Scan ALL 9HBVD/9HCVQ Samples on the Task list (create a new 9SERO & save it) separately,
12. Go to Batch Resulting Work list, choose task list, scan Task list ID-put 9HBVD or 9HCVQ on test box, press OK, the list will show up, mark all, press add result, click the small box verified all then ok, save your work
13. Pack all 9HBVD , 9HCVQ/9HCVG aliquots in a white box sort them all in order with the PHOL form, pack both in a clear plastic zip lock
14. All PHOL samples are packed in blue PHOL bags to be ready for pick-up Mon-Fri at

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>University Health Network Sciences General Sciences Women Sciences Children Sciences Mental Health</div><div>Mount Sinai Hospital Hôpital du Mont-Sinai Hospital de Mont-Sinai Hospital de Gracia</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 22 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

11:00am

## **HIV VIRAL LOAD FOR MSH patients only**

**\*\*Do not enter in Softlab\*\***

Specimen Type Required:

EDTA or PPT (inert gel and spray-dried K2EDTA anticoagulant)

1. Match sample label information with patient's information in requisition form.
2. Complete Section 5: B and D of HIV Viral Load Test Requisition Form; C (if sample will be frozen)

<input type="checkbox"/> Lamivudine (3TC)	3TC	<input type="checkbox"/> Reyataz (Atazanavir)	ATV	<input type="checkbox"/> Vocabria (Cabotegravir)	CAB
<b>4 - Comments:</b>					
<b>5 - Collection Information</b> Sections A, B and C must be completed by the submitter accordingly with each sample submitted.					
A. Collected (yyyy-mm-dd):		Time (24hr - HH:MM)	Initials:	C. Frozen (< -20°C) (yyyy-mm-dd): Time (24hr - HH:MM) Initials:	
B. Plasma separated (yyyy-mm-dd):		Time (24hr - HH:MM)	Initials:	D. Received (yyyy-mm-dd): Time (24hr - HH:MM) Initials:	


### CONFIDENTIAL WHEN COMPLETED

The personal health information is collected under the authority of the Personal Health Information Protection Act, 2004, s.36 (1)(c)(iii) for the purposes specified in the Ontario Agency for Health Protection and Promotion Act, 2007, s.1 and will be used for surveillance and other public health purposes. If you have questions about the collection of this personal health information please contact the PHOL Manager of Customer Service at 416-235-6556 or toll free 1-877-604-4567. Form No. F-C-HV-139-008 (2023-04-18)



3. Label one Greiner tube (provided by PHOL) with Patients Last Name, First Name and MRN and Date of Birth.

Please refer to Greiner tube aliquoting instructions provided by PHOL below.

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 23 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	




#### SERUM/PLASMA

Centrifuge unspun specimen tubes at 3,000 rpm for 10 minutes. Separate the serum/plasma into the screw-capped cryotubes. Follow the aliquot volume guide below. **Do not use parafilm to seal tubes.**

Tube type	Maximum aliquot
Greiner external thread cryovial (4.5mL) or equivalent	3mL

Please refer to the diagram below on how to place barcodes on the test tubes. The barcodes must be placed accordingly to prevent instrument barcode reading errors.

|




- Aliquot minimum 2.5 ml of plasma maximum 3 ml to 3.5 ml. Store original tube in Done Alinity Rack for scanning

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is **Uncontrolled When Printed**.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 24 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

5. Sample received after PHL pick up time 11AM (Monday –Friday)

As per PHOL Submission and Collection Notes

“Plasma can be stored and/or transported at 2°C - 8°C for up to 6 days after separation or at  $\leq -18^{\circ}\text{C}$  for up to 12 weeks. For extended storage of plasma must be frozen at  $\leq -60^{\circ}\text{C}$ .”

6. Document by writing patient information and test requested in PHL Send out logbook
7. Send sample with to PHL with completed Viral Load Requisition Form.

### **TWH Liver Clinics:**

Samples are sent Tuesday -Friday at 09:00 am Place a sticker on bag (blue bags used for UHN samples) samples on Friday after 9:00am should be frozen at -20C. The TWH liver clinic will send someone on **Mondays** to pick up the **Frozen Samples**.

#### **LIVER CLINIC (aliquots)**

1. Save all Liver Clinic aliquots in fridge – MIRT13 – do not put in freezer
2. Send the aliquots the following day. Place sample inside the Blue Bag with ice pack and Sticker / Label at Plating Area- Deliver to UHN rack
3. Prepare a sticker or label package as:

**To: TWH Core Lab**

**PLEASE  
CALL 6232 FOR PICK-UP**

### **Measles/Mumps Diagnostic IgM/IgG Serology or PCR Send out**


- Nasopharyngeal, throat swab, urine or CSF can be tested for measles PCR
  - Buccal swab, throat swab, saliva, urine or CSF can be tested for mumps PCR
  - Serum is used for Mumps/Measles IgM testing.
1. Receive, print and complete a MIC10 format PHOL report from result entry in the LIS.
  2. Pack the sample according to TDG guidelines.
  3. Email ICP of submitting hospital specimen was received and is being to PHOL for Mumps / Measles testing with specimen LIS#.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 25 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

4. Phone PHOL customer service first to let them know a sample is on the way for the next available run. If unable to reach customer service, contact **Duty officer 416 605 3113**
5. [Send samples by KJV](#). **Send specimens to PHOL once per day at 6am during weekdays.**

Note: If STAT testing is requested, IPAC or requesting physician can contact and arrange STAT testing directly with PHOL.

### **Monkeypox Sendout**

See [Monkeypox Procedure](#)

### **Middle Eastern Respiratory Syndrome Coronavirus – MERS-CoV / AVIAN FLU PCR Send Out**

Specimens for MERS-CoV testing will be processed same day if received by PHOL by noon.

All primary specimens can be handled in a containment level 2 lab.

Public Health Agency of Ontario instructions regarding specimen collection for MERS-CoV that includes appropriate specimen types, biosafety procedures:

<http://www.health.gov.on.ca/en/pro/programs/publichealth/coronavirus/guidance.aspx>

Instructional video regarding NP swab collection

<http://microbiology.mtsinai.on.ca/research/icuInfluenza/noseSwab.asp>

1. Requests for MERS-CoV or Avian flu requires STAT notification to:
  - ICP for the relevant hospital/ward to be paged;
  - Prior to specimen submission, contact PHO's Laboratory Customer Service Centre at 416-235-6556/1-877-604-4567 or the After-Hours Emergency Duty Officer at 416-605-3113.

#### **NOTE:**


- Ensure the clinical information is available before calling PHOL
  - Page microbiologist-on-call for any missing information as identified by PHOL Customer Service/After-Hours Emergency Duty Officer
2. Clinical information is **required** by PHOL
  3. Specimens requested can include:
    - Upper respiratory tract: NP swab AND throat swab

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is **Uncontrolled When Printed**.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 26 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

- Lower respiratory tract: sputum or BAL
- EDTA blood (purple top tube)
- If the patient has gastrointestinal symptoms, submit a stool sample in a dry sterile container - do not use other transportation media containers such as Cary-Blair.
- acute (when patient first seen with acute respiratory illness) and convalescent (21 to 28 days after illness onset) serology

**Note for suspected Avian Flu respiratory specimens:**


- Heat inactive primary specimens
  - Perform Flu/RSV testing in house on Aries platform
  - Send specimen to PHOL for AVIAN Flu PCR whether in-house flu is positive or negative
4. Receive, print and complete a MIC10 format PHOL report from result entry in the LIS with clinical information.
  5. Pack the sample according to TDG guidelines. Mark package as “STAT” testing.
  6. Notify PHOL customer service if not done by IPAC and [Send samples by KJV](#) to PHOL immediately. Contact PHOL **Duty officer 416 605 3113** if unable to reach customer service.

**WNV Molecular Testing to PHOL**

Non-donor/non-recipient blood and CSF for WNV molecular testing needs approval of the Microbiologist-on-call

**Zika Virus Molecular Testing to PHOL**

Please submit Mandatory Information Intake Form for Zika Virus Testing together with PHOL Requisition Form, refer to PHOL website for form if required.



 Department of Microbiology Quality Manual	Policy # MI_SO Version: 4.2 CURRENT	Page 27 of 59
<b>Section:</b> Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

## **OTHER MISCELLANEOUS NATIONAL (NON-GTA) AND INTERNATIONAL SEND OUT SPECIMENS**

Other miscellaneous routine specimens require testing to reference laboratories outside the GTA. These specimens are sent daily through Fedex.

1. Using the [Sendout Shipping Information Table](#) search for the specified test requested. Information on sample type, requisition, packaging and shipping addresses among other information can be found in this table.
2. Once the specimen is accessioned in the LIS, parafilm the container lid and put the specimen in a plastic biohazard bag containing and absorption pads (to absorb any possible spills). (You can use the plastic SAF T Pak bag). Place requisition in the outside sleeve of the biohazard bag.
3. Package bagged specimen(s) as per [TDG guidelines](#). Ensure to affix sender and recipient address labels.
4. See [Specimen packaging](#) section for Fedex specific packaging instructions
5. See [Courier](#) section for requesting Fedex.

Note: all Chlamydia positive rectal swabs will be reflexed for LGV testing and sent to NML.

  Department of Microbiology Quality Manual	Policy # MI_SO Version: 4.2 CURRENT	Page 28 of 59
<b>Section:</b> Bacteriology Procedures	<b>Subject Title:</b> <b>Send-out Information Manual</b>	

## **COURIERS**

	<b>Courier</b>	<b>Monday - Friday</b>	<b>After hours/ Weekends/ Holidays</b>
<b>Routine</b>	Life Labs Courier (UHN, CAMH, Bridgepoint)	Delivery: 8:00am-9:00pm & 1 :00am, 4:00am	Delivery: 8:00am-5:30pm
		Pick up: 8:00am - 8:00pm	Pick-up: 8:00am-5:00pm
	PHOL Courier (PHOL specimens)	Pick up only: 10:30am	N/A
<b>Special</b>	Fedex Courier (Non-GTA routine specimens)	Pick up only: 4:00pm	N/A
	<a href="#">KJV Courier</a> Transportation not mentioned above	9:00am- 5:00pm (online)	Anytime (by phone request)
	<a href="#">Taxi</a> (If all above unavailable)	Anytime	Anytime

### **Life Labs Courier Specimen Transportation**

Life Labs courier is used for transportation of specimens between UHN and UHN client hospitals.

To send out a sample using Life Labs Courier:


1. Place specimens in a biohazard specimen bag. If there is a requisition, place it in the outer sleeve of the bag.
2. Ensure final packaging is a Blue Life Labs Courier Bag within the hard plastic container.
3. Ensure the insert at the top of the bag, indicates receiver hospital.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 29 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

4. Place Life Labs Courier Bag on the shelf beside the microbiology specimen receiving window.

### **PHOL Courier Specimen Transportation**

PHOL courier is used for transportation of routine specimens to PHOL.

To send out a sample using PHOL Courier:

1. Place specimens in a biohazard specimen bag, with the [PHOL requisition](#) in the outer sleeve.
2. If applicable, arrange all samples in order in a foam carrier. Keep requisitions in the same order as the samples and place both in a clear plastic zip lock bag
3. Ensure final packaging is a Blue PHOL Courier Bag within the hard plastic container.
4. Place PHOL Courier Bag in a fridge if courier is not coming imminently, if not place the Courier Bag on the trolley by the serology specimen management bench.

For PHOL directions and location for all deliveries (including after-hours, weekends and STAT specimens) see link: [PHOL Directions and Receiving Dock Instructions](#)


Note: For STAT deliveries, contact PHOL **Duty officer 416 605 3113** if unable to reach customer service.

### **FEDEX Courier Specimen Send Out**

Fedex courier is used for International and National non-GTA routine miscellaneous specimens.

To send out a sample using Fedex:


1. Package specimen(s) as per [TDG guidelines](#).
2. For Courier service go to : [http://www.fedex.com/ca\\_english/](http://www.fedex.com/ca_english/)  
Log in using current ID and password.

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 30 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title: Send-out Information Manual</b></p>	

3. Complete a FedEx Express form online.

a. After logging in, the following window pops-up.

In the first section-My Shipment Profile, click drop down arrow and choose the the name of the person to whose attention shipment is being sent.



The screenshot shows the FedEx Ship Manager interface. A red box labeled "Drop Down Arrow" points to the dropdown menu in the "My Shipment Profiles" section. The page includes sections for shipping information, billing details, pickup/drop-off, and contact information.

b. Scroll down to the Package & Shipment Details, fill in No. of packages and Weight.


Scroll Bar


UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is *Uncontrolled When Printed*.


Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 31 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> Send-out Information Manual</p>	



c. Click “SHIP” in section 6



d. Click “Ship” in the section of Confirm your shipment details

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 32 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

**Confirm your shipment details**

**Outbound Shipment** [Help](#)

<b>From</b> Micro Lab Mount Sinai Hospital 600 University Avenue Microbiology Room 1470 Toronto, ON Canada M5G1X5 4165864432	<b>Ship date</b> 04/01/2014 <b>Number of packages</b> 1 <b>Total weight</b> 1.00 kgs <b>Declared/Carriage Value</b> 0.0 CAD <b>Pricing option</b> <b>Service type</b> Priority Overnight <b>Package type</b> FedEx Box <b>Dimensions</b> <b>Bill transportation to</b> MicroLab-821 <b>Your reference</b> <b>Pickup/Drop-off</b> Use an already scheduled pickup at my location  <b>P.O. no.</b> <b>Invoice no.</b> <b>Department no.</b> <b>Special Services</b>
<b>To</b> Dr. Michael Mulvey David Boyd National Microbiology Laboratory Nosocomial Infections 1015 Arlington Street, Suite T2380 WINNIPEG, MB Canada R3E3R2 204-789-5000	


[Edit](#) [Ship](#)


- e. A screen will appear displaying your waybill. Print **TWO** copies of the waybill.

[Print](#) [Receipt](#) [View](#)


[Create return shipment](#) [Edit shipment](#) [Cancel shipment](#) [Repeat last](#) [Create new](#)

[Schedule a Pickup](#)

From: (416) 586-4432      Origin ID: YOOA            Ship Date: 01APR14  
Micro Lab      Mount Sinai Hospital      Act/Ngt: 1.0 KG  
600 University Avenue      Microbiology Room 1470      CAD: 103260135/NCA3480  
Toronto, ON M5G1X5

Delivery Address Bar Code  


4. Retrieve a FedEx International Shipment Documentation holder. Open the purple zip-lock holder and place the completed [FedEx sheet](#) face up (so that information is showing on the clear plastic side with all barcodes showing).

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 33 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	

5. Affix the prepared FedEx International Shipment Documentation holder outside the prepared package. Ensure Sender and Receiver addresses are also on the outside of the package. Ensure final specimen packaging is labelled UN3373.
6. Staple the sender's portion of the FedEx Express online form to the photocopied send out requisition and file in the **Shipping Documents Binder** on the send out bench.
7. Leave the completed FedEx Clinical Pak package for pickup by FedEx on the ledge at the reception desk outside the Lab. near Secretary's desk. The courier makes pickup at 16:00 hours Monday to Friday.

### **KJV Courier Specimen Transportation**

KJV Courier service is used for any GTA deliveries not serviced by routine couriers during regular hours or for deliveries outside of regular courier hours see [Courier Table](#).

To send out a sample using KJV:

1. Package specimen(s) as per TDG guidelines.
2. Ensure final specimen packaging is a SAF T PAK labelled UN3373.
3. For Courier service:
  - a. If KJV's service is required on Monday-Friday between the hours of 9am to 5pm use their website online to request their services by following to Step 4.
  - b. If KJV's service is required after hours (outside 9am - 5pm), on weekends or holidays, call 647-212-2395 to request their service. Follow instructions given as per phone conversation.



**If wait time is longer than expected, call for delivery by taxi using [STAT Transportation by Taxi - Procedure QEPMI5007i](#)**
4. Click KJV website link to open website: <http://www.kjvcourier.com>  
Beside "Existing Accounts" click "Online Services"
5. Enter Account No. & Password and click "Yes" circle for saving account number on this computer. Then click on "Login".

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 34 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p>Subject Title: <b>Send-out Information Manual</b></p>	



6. On the Main Page under Main Menu: Click on “Enter New Order”. A page will pop up to create a “New Order” page.
  7. Under Main Menu:
    - a. Click the down arrow ▼ to access the drop down bar at the top of the form, select an address for **pick up** and click  to the **left**.
    - b. again, click the down arrow ▼ to access the drop down bar at the top of the form. This time select an address **for delivery** and click  to the **right**.  
Addresses will automatically populate.
  8. Complete remaining form filling in all sections with \*asterisks and any “Special Instructions”.  
For the sections below always select the following:
    - a. Person Placing Order → Use your own name or “Microbiology Tech”
    - b. Service → Same day VIP
    - c. Vehicle → Car
    - d. Number of Pieces → Type in number
    - e. Description → Parcel
  9. Click on “Continue” at the bottom of the page to proceed.
  10. After you confirm order information is correct, select
  11. After you process the order, a Waybill Page will pop up.
  12. Print Waybill page and keep it open. Secure the bottom half of the page; “KJV’s Copy” onto the package. The top half of the page is not needed.
  13. Go back to the Waybill and save a copy in the file “KJV Courier Waybills” found here:  
**T:\Microbiology\KJV COURIER - Waybills**. Save the file as the date. (eg. 2015.04.21)
  14. Once Steps 1-9 are completed, click on “Return to Main” in the middle of the page and SIGN OUT.
  15. Leave the packed specimen on the Specimen Management Bench Sorting Area for pick up
- For PHOL directions and location for all deliveries (including after-hours, weekends and STAT specimens) see link below:

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*



Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div></div> <div>Department of Microbiology</div> <div>Quality Manual</div>	Policy # MI_SO	Page 35 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

## [PHOL Directions and Receiving Dock Process](#)

### **TAXI Specimen Transportation**

1. Package bagged specimen(s) as per [TDG guidelines](#).
2. Ensure final specimen packaging is a SAF T PAK labelled UN3373.
3. Refer to [STAT Transportation by Taxi - Procedure QEPMI5007i](#) for instructions on sending specimens by Taxi.

<div></div> <div>Department of Microbiology</div> <div>Quality Manual</div>	Policy # MI_SO	Page 36 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

**Related Documents:**


Document Title	Document Number
<a href="#">Sendout Shipping Information Table</a>	
<a href="#">Referred Out Test Results Entry</a>	
<a href="#">Post-Analytical - Refer-out Test Results Reporting Procedure</a>	QPCMI18001

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*


*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	Policy # MI_SO	Page 37 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

## APPENDIX I: SEROLOGY TESTS ORDERING SCREENS

**Table 1: serology in house from 1-R**



The screenshot shows the 'Order Entry - [New Order - Edit Mode]' window. The 'Keypad' menu is open, displaying a list of serology tests. The tests are organized into columns and rows, with checkboxes for selection. The tests include:

- 1 Hep A IgM Antibody
- 2 Hep A IgG Antibody
- 3 Hep B Surface Antigen
- 4 Hep B Surface Antibody (HBsAb)
- 5 Hep B Core Antibody
- 6 Hep B Core IgM Antibody
- 7 Hep Be Antigen
- 8 Hep Be Antibody
- 9 Hep C Antibody
- A HTLV 1&2 Antibody
- B HIV 1&2 Antibody/p24 Antigen Test
- C Hep B DNA
- D Hep C RNA Quantitative
- E Hep B Genotype
- G HBV Antiviral Resistance Report
- H CMV IgG Antibody
- I Rubella IgG Antibody
- J Syphilis Screening
- K VZV IgG Antibody
- L Infectious Mono Screen
- M EBV VCA IgG
- N EBV EBNA IgG
- O Galactomannan - Blood
- P Galactomannan - BAL
- Q Serum Crypto Antigen
- R Serum Storage
- S Prenatal PRE
- U Donor DONOR
- V Ref-out REF1
- W More Ref REF2
- X More Ref REF3
- Y Back To OEMIC
- Z Main OE DETST


The window also shows patient information on the left, including Last name, DOB, Country, and HCN. The bottom status bar indicates the user is MISAA, Env: LIVE2, Term: M4705, NUM: EDIT, 2013.05.14, 12:32.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*


*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


<div><p>Department of Microbiology</p></div> <div>Quality Manual</div>	Policy # MI_SO	Page 38 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	


**Table 2: serology-prenatal**




 <p>Department of Microbiology</p> <p>Quality Manual</p>	Policy # MI_SO	Page 39 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	


**Table3: serology-donor**




 <p>Department of Microbiology</p> <p>Quality Manual</p>	Policy # MI_SO	Page 40 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

### Tables 4, 5 & 6: serology-refer out



 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 41 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> Send-out Information Manual</p>	




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*


Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 42 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	


## **APPENDIX II: VIROLOGY/SEROLOGY PHOL PENDING LIST**


Procedure in LIS:

1. In SoftLab, open “Resulting Worklist”.
2. Use Template “**9REF**”; specify date range (i.e. From 4 months ago until 2 weeks prior); status use “Pending+Non-verified”; include both received and unreceived samples (uncheck “Received only”)



3. Once pending list appears, right click mouse anywhere on this list to customize what columns to display.

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 43 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	





Note: the column parameters in Red box are what appear at default.

- Add to show columns C Date (Collection Date) and R Date (Receive Date) by selecting each of those on the *left*, then press the “->” Key.

REMOVE the C Time (Collection Time) and R Time (Receive Time) by selecting each of those on the *right*, then press “<-” Key.

Press OK.



 Department of Microbiology Quality Manual	Policy # MI_SO	Page 44 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

- Once list appears again, filter first by Receive Date by pressing on the “R Date” column tab until all the unreceived samples appear on the top (the “R Date” would be BLANK). Cancel ALL orders that were not received by lab for  $\geq 7$  days as follows:

**Result field:** CANCELLED


**Comment field:** This specimen was documented as collected but not received



by the Microbiology/Virology Laboratory after 7 days. If this test is still required, please place a new order.


- Sort the Pending List by Ward. Result all tests ordered under any Occupational Health and Safety wards (i.e. OHS, GOCC, PMOCC, WOCC) with:

“Report is sent directly to Occupational Health and Safety.”


Verify results!


- Refresh the Pending List by pressing “”
- Sort the Pending list by Receive Date again so that oldest specimens received are at the top.
- Go to SoftLab Desktop → Reports → SQL Reports → EQ Print SQL Reports

<div><div><div><div>UHN</div></div><div><div>Toronto General Toronto Western Princess Margaret SickKids</div></div></div><div><div><div>Mount Sinai Hospital</div></div><div><div>Mount Sinai St. Michael's Hospital</div></div></div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 45 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	



10. Once screen bridges to SQL Report program, use arrow keys to scroll to and highlight “PENDING\_PHOL”. Press “V” for View (Red Circle).




 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 46 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> Send-out Information Manual</p>	

11. Enter same date range as the pending list pulled in LIS (i.e. 20171102 [Tab] 20180313)



Press “Enter” and wait until list appears. (May take up to a minute)




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 47 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> Send-out Information Manual</p>	

12. Press “P” for Print. Another list will appear on the top left. Use ↑↓ or PageUp keys to scroll to the very top of that list until cursor highlights “Save Save to File”.



Press “Enter”.

13. At Save to box, Enter “pend+today’s date” (i.e. pend20180327).




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*


Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	Policy # MI_SO	Page 48 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	


Press “Enter” → Select “Y”. However, this file will ONLY be saved for **6 days**.

Now you can close this window.

14. Go back to SoftLab Desktop → [SHIFT]+[F2] → Type in File name of “Save to” file (i.e.pend20180327) → Click “Find”. When file appears, double click on file to open.



15. When list appears on screen; Open File tab → Export to Excel




16. Uncheck “Comma” box → Check “Other” Box → Type in “|” ([SHIFT]+[|]) in box beside “Other” → Click “Export”

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 49 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title: Send-out Information Manual</b></p>	



17. When exported excel file appears, expand all columns to show all contents.

18. FORMAT the excel file as follows:

- a. **BOLD FONT** and Add ***FILTER*** (Black arrow) to Row 1.

Order	ColDate	ColTime	RcvDate	RcvTime	LastName	FirstName	HCN	DOB	Test	
2	S9161795	01/12/2018	12:54	01/16/2018	12:32	WILLEMSMA	JUDY	19610817	9HSG1	Herpes Simplex IgG
3	T0062355	02/06/2018	13:39	02/06/2018	16:01	LEUNG	YIU	3997747161YX	19550102	Hepatitis B DNA
4	T0061235	02/06/2018	10:40	02/06/2018	14:35	HUSSAIN	SHAHABUB	6169284970LA	19660505	Hepatitis B DNA
5	T0061715	02/06/2018	11:36	02/06/2018	16:21	KWOK	FANNY	8943755143GE	19750517	Hepatitis B DNA
6	T0061762	02/06/2018	11:36	02/06/2018	16:07	CHOW	HONG	7873200485DI	19670718	Hepatitis B DNA

- b. Insert Header by:

- i. Go to Insert Tab and select “Header & Footer” (Red Arrow below)

Order	ColDate	ColTime	RcvDate	RcvTime	LastName	FirstName	HCN	DOB	Test	
2	S9161795	01/12/2018	12:54	01/16/2018	12:32	WILLEMSMA	JUDY	19610817	9HSG1	Herpes Simplex IgG
3	T0062355	02/06/2018	13:39	02/06/2018	16:01	LEUNG	YIU	3997747161YX	19550102	Hepatitis B DNA
4	T0061235	02/06/2018	10:40	02/06/2018	14:35	HUSSAIN	SHAHABUB	6169284970LA	19660505	Hepatitis B DNA
5	T0061715	02/06/2018	11:36	02/06/2018	16:21	KWOK	FANNY	8943755143GE	19750517	Hepatitis B DNA

- ii. Type in Header (Use **BOLD** font and Font size 16):


**Please fax back reports to (416) 619 5522-VERIFIED BY PHL. Thank you.**


UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*


Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 50 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> Send-out Information Manual</p>	




- iii. Click left mouse key anywhere OUTSIDE of the header box → Go to View Tab → Select “Normal” view (Red arrow below)


This will bring the excel worksheet back to the normal view of list of orders and tests as before.



- c. **FREEZE** First Row by:


Go to View Tab → Under “Freeze Panes” (Red Arrow) → Select “Freeze Top Row”

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 51 of 59</p>
<p><b>Section: Bacteriology Procedures</b></p>	<p><b>Subject Title: Send-out Information Manual</b></p>	



Order	ColDate	ColTime	RcvDate	RcvTime	LastName	FirstName	HCN	Test	Result
2	S9161795	01/12/2018	12:54	01/16/2018	12:32	WILLEMSMA	JUDY	19610817	9HSG1 Herpes Simplex IgG
3	T0062355	02/06/2018	13:39	02/06/2018	16:01	LEUNG	YIU	19550102	9HBD1 Hepatitis B DNA
4	T0061235	02/06/2018	10:40	02/06/2018	14:35	HUSSAIN	SHAHABUB	19660505	9HBD1 Hepatitis B DNA
5	T0061715	02/06/2018	11:36	02/06/2018	16:21	KWOK	FANNY	19750517	9HBD1 Hepatitis B DNA
6	T0061762	02/06/2018	11:29	02/06/2018	16:07	CHOW	HON	19620218	9HBD1 Hepatitis B DNA
7	T0062352	02/06/2018	13:41	02/06/2018	16:05	ZHENG	ENHAO	19810510	9HBD1 Hepatitis B DNA
8	T0061963	02/06/2018	12:44	02/06/2018	16:06	KOLA	KOSTANDIN	19410110	9HBD1 Hepatitis B DNA
9	T0071949	02/07/2018	10:15	02/07/2018	15:40	HEDDERSON	SEAN	19610322	9HBD1 Hepatitis B DNA
10	T0080628	02/08/2018	06:40	02/08/2018	10:51	RUMBLE	GAIL	19481003	9HIS1 Histoplasma Serology
11	T0080277	02/08/2018	07:02	02/08/2018	10:52	SAVAGLIO	RITA	19520802	9CMM1 CMV IgM
12	T0080277	02/08/2018	07:02	02/08/2018	10:52	SAVAGLIO	RITA	19520802	9HSG1 Herpes Simplex IgG
13	T0110533	02/11/2018	07:41	02/12/2018	12:55	ZAROD	ADAM	19480706	9HPY1 Helicobacter pylori Serology
14	T0191601	02/19/2018	19:57	02/20/2018	18:45	BARRETT	DEAN	19910325	9TSC1 HIV
15	T0200816	02/20/2018	13:36	02/20/2018	13:36	CHAGANY	ABDULREHMAN	19540719	9CMM1 CMV IgM
16	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9HPY1 Helicobacter pylori Serology
17	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9PST1 Strongyloides Serology
18	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9SCH1 Shistosoma Serology

19. Sort by Receive Date (Click where black arrow points and select Sort A to Z)



Order	ColDate	ColTime	RcvDate	RcvTime	LastName	FirstName	HCN	Test	Result
2	S9161795	01/12/2018	12:54	01/16/2018	12:32	WILLEMSMA	JUDY	19610817	9HSG1 Herpes Simplex IgG
3	T0062355	02/06/2018	13:39	02/06/2018	16:01	LEUNG	YIU	19550102	9HBD1 Hepatitis B DNA
4	T0061235	02/06/2018	10:40	02/06/2018	14:35	HUSSAIN	SHAHABUB	19660505	9HBD1 Hepatitis B DNA
5	T0061715	02/06/2018	11:36	02/06/2018	16:21	KWOK	FANNY	19750517	9HBD1 Hepatitis B DNA
6	T0061762	02/06/2018	11:29	02/06/2018	16:07	CHOW	HON	19620218	9HBD1 Hepatitis B DNA
7	T0062352	02/06/2018	13:41	02/06/2018	16:05	ZHENG	ENHAO	19810510	9HBD1 Hepatitis B DNA
8	T0061963	02/06/2018	12:44	02/06/2018	16:06	KOLA	KOSTANDIN	19410110	9HBD1 Hepatitis B DNA
9	T0071949	02/07/2018	10:15	02/07/2018	15:40	HEDDERSON	SEAN	19610322	9HBD1 Hepatitis B DNA
10	T0080628	02/08/2018	06:40	02/08/2018	10:51	RUMBLE	GAIL	19481003	9HIS1 Histoplasma Serology
11	T0080277	02/08/2018	07:02	02/08/2018	10:52	SAVAGLIO	RITA	19520802	9CMM1 CMV IgM
12	T0080277	02/08/2018	07:02	02/08/2018	10:52	SAVAGLIO	RITA	19520802	9HSG1 Herpes Simplex IgG
13	T0110533	02/11/2018	07:41	02/12/2018	12:55	ZAROD	ADAM	19480706	9HPY1 Helicobacter pylori Serology
14	T0191601	02/19/2018	19:57	02/20/2018	18:45	BARRETT	DEAN	19910325	9TSC1 HIV
15	T0200816	02/20/2018	13:36	02/20/2018	13:36	CHAGANY	ABDULREHMAN	19540719	9CMM1 CMV IgM
16	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9HPY1 Helicobacter pylori Serology
17	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9PST1 Strongyloides Serology
18	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9SCH1 Shistosoma Serology
19	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9TSC1 HIV
20	T0202803	02/20/2018	14:23	02/21/2018	12:24	SHAW	BIBI	19521031	9BBI1
21	T0262496	02/26/2018	14:29	02/26/2018	18:08	JAO	YUNG-WO	19390203	9HBD1

**\*\*This excel list should now essentially contain identical patient orders as the LIS 9REF pending worklis**

Save file type As *Excel Workbook* with the same name Today's Date + Initials (i.e. 20180327jkc) under: **T:\Microbiology\Virology\Faxed PHL pending Lists\2018**


**NOTE:** If date file saved is not the same as date faxed, save using the DATE FAXED!

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>University Health Network Toronto General Toronto Western St. Michael's Toronto Mount Sinai Hospital Division of Health Sciences Department of Microbiology</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 52 of 59
	Version: 4.2 CURRENT	
Quality Manual		
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

20. Go through LIS pending worklist by bridging to Order Entry to see which orders need to be faxed to PHOL according to priority.

If an order needs to go on a PHL pending list and be faxed, order 8PFAX and document the date faxed (F6) in the comment field. DO NOT VERIFY YET! Order 8PFAX in **each** order that was faxed.

\*\*\*Verify 8PFAX **only** when all results from PHOL have been entered and verified for that order.

i.e. If an order is already on a previous pending list and have been faxed but still waiting for results, it should have 8PFAX ordered but unverified.

**Note:** For 2<sup>nd</sup> time and subsequent faxes → Add subsequent dates those faxes were faxed again in the comment field underneath first time it was faxed

Sort by:		
Test	Receive Date	Ward
To deal with STAT tests first (i.e. HIV, prenatal tests etc...)	To get rid of orders not received by the lab for >7 days	To result and verify tests ordered under any Occ. Health & Safety wards
To verify any unverified test results or sent out comments (i.e. tests with result but not verified; any tests ending in "2" such as 9TSC2, 9HBD2 etc...)	To deal with oldest test orders received first	To deal with critical or certain wards first upon request
To deal with tests based on their TAT from PHOL		


21. On exported SQL pending list, DELETE all orders that DO NOT require fax. Retain **ONLY** the ones you will fax now, in batches of no more than 10 specimen order numbers per fax. (Same order number can appear on different rows for a different test)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_SO</p> <p>Version: 4.2 CURRENT</p>	<p>Page 53 of 59</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> <b>Send-out Information Manual</b></p>	



If need to fax >10 orders, break into batches of 10 orders per fax in multiple excel sheets by cutting and pasting into a new excel tab within the same file (ONLY if these batches are faxed on the same day).

22. Create a different excel file name if a second batch of orders are faxed on a different day! Save using the DATE FAXED.
23. Translate LIS test codes to corresponding test names and enter in a separate column to the right of the test codes (i.e. 9TSCX → HIV See diagram in Step #19)
24. Print this list by:

Highlighting the cells you wish you print → File → Print → Under Settings, select “Print Selection” (Red Circle) and choose “Fit All Columns on One Page” (Red Arrow); choose Landscape Orientation (Blue Circle) → Print (Blue Arrow)

Fax printed sheet to PHOL @ 416-265-6552.



<div><div>University Health Network Toronto General Toronto Western Toronto St. Michael's Toronto St. Paul's Toronto St. Joseph's Toronto St. Elizabeth's Toronto St. Michael's Toronto St. Paul's Toronto St. Joseph's Toronto St. Elizabeth's</div><div>Mount Sinai Hospital Department of Microbiology Division of Infectious Diseases Division of Clinical Microbiology Division of Laboratory Medicine</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 54 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

### **APPENDIX III: UHN MANAGED MICROBIOLOGY SEND OUT TESTS**

For UHN patient samples only, the follow list of send out tests are managed by the UHN specimen processing team directly and sent directly to PHOL.

If Received at Sinai:

- DO NOT “add on” test to orders within our department for the list below.
- Direct samples back to UHN laboratory, phone them for awareness.

If encountered:

- Request test be ordered in EPIC by the requestor if needed
- If an add-on request or circumstances require sample retrieval or submission between UHN/MSH, work collaboratively with the UHN core laboratory to make arrangements.


AMOEBIC SEROLOGY	DENGUE IgM SEROLOGY	HEPATITIS C (HCV) DRUG RESISTANCE ANALYSIS	MUMPS ACUTE/ACTIVE DISEASE (IGM)	PINWORM PREP	TOXOPLASMA TOTAL ANTIBODY
ANAPLASMA SEROLOGY	EBV EA IGG	Hepatitis C Ab PHOL Confirmation	MUMPS IMMUNITY (IGG/TOTAL ANTIBODY)	Q FEVER SEROLOGY	TRICHINELLA SEROLOGY
ARBOVIRUS SEROLOGY	EBV VCA IGM	HEPATITIS E ANTIBODY	MUMPS MOLECULAR DETECTION, OTHER	RABIES SEROLOGY	TRYPANOSOMA SEROLOGY
BARTONELLA HENSELAE SEROLOGY	ECHINOCOCCUS SEROLOGY	HERPES SIMPLEX IGG ANTIBODY	MUMPS MOLECULAR DETECTION, URINE/CSF	RICKETTSIA SEROLOGY	TULAREMIA SEROLOGY
BLASTOMYCES SEROLOGY	FILARIA (SKIN SNIP FOR MICROFILARIAE)	HISTOPLASMA SEROLOGY	MURINE TYPHUS SEROLOGY	ROCKY MOUNTAIN SPOTTED FEVER SEROLOGY	VARICELLA ZOSTER IGM
BRUCELLA SEROLOGY	FILARIA SEROLOGY	HISTOPLASMA SERUM ANTIGEN	MYCOPLASMA DETECTION	RUBELLA, IGM	WEST NILE VIRUS ANTIBODY, IGG
CHIKUNGUNYA IGG ANTIBODY	FUNGUS CULTURE - NAIL/HAIR/SKIN	LEPTOSPIRA SEROLOGY	OVA AND PARASITE, OTHER	SCHISTOSOMA SEROLOGY	WEST NILE VIRUS ANTIBODY, IGM
CHIKUNGUNYA IGM ANTIBODY	H.PYLORI SEROLOGY	LYME DISEASE SEROLOGY	OVA AND PARASITE, STOOL	SPOROTHRIX ANTIBODY	WORM/INSECT IDENTIFICATION
CHLAMYDIA PNEUMONIAE SEROLOGY	HBV DNA	MEASLES ACUTE/ACTIVE DISEASE (IgM)	OVA AND PARASITE, URINE	STRONGYLOIDES SEROLOGY	YELLOW FEVER IgG ANTIBODY
COCCIDIOIDOMYCES SEROLOGY	HCV GENOTYPING	MEASLES IMMUNITY (IGG/TOTAL AB)	PARVOVIRUS B19 ACUTE/ACTIVE DISEASE (IGM)	TETANUS SEROLOGY	YELLOW FEVER IgM ANTIBODY
CYSTICERCOSIS SEROLOGY	HCV RNA	MEASLES MOLECULAR DETECTION, CSF	PARVOVIRUS B19 IMMUNITY (IGG)	TOXOCARA SEROLOGY	YERSINIA SEROLOGY
DELTA VIRUS ANTIBODY	HCV RNA (TRANSPLANT)	MEASLES MOLECULAR DETECTION, URINE	PHL Blood Parasite Confirmation and Speciation	TOXOPLASMA ANTIBODY, IGG	ZIKA VIRUS IgG SEROLOGY
DENGUE IgG SEROLOGY	HEPATITIS B SURFACE ANTIGEN (PHL)			TOXOPLASMA ANTIBODY, IGM	ZIKA VIRUS IgM SEROLOGY

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is *Uncontrolled When Printed*.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 55 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

### **Record of Edited Revisions**

#### **Manual Section Name: Send-out Information Manual**



<b>Page Number / Item</b>	<b>Date of Revision</b>	<b>Signature of Approval</b>
Annual Review	May 1, 2002	Dr. T. Mazzulli
Annual Review	May 12, 2003	Dr. T. Mazzulli
Annual Review	May 26, 2004	Dr. T. Mazzulli
Multiple resistant organism multidrug bactericidal testing at CHEO – added	September 11, 2004	Dr. T. Mazzulli
CHEO instruction for Transplant team when requesting testing.	December 16, 2004	Dr. T. Mazzulli
Annual Review	May 12, 2005	Dr. T. Mazzulli
Send-out Bench duties	April 01, 2006	Dr. T. Mazzulli
Malaria send-out protocol moved from Parasitology Manual	April 01, 2006	Dr. T. Mazzulli
SEND OUT to PHOL procedure for after-hours, weekends and holidays added	April 01, 2006	Dr. T. Mazzulli
Annual Review	April 1, 2006	Dr. T. Mazzulli
Histoplasma antigen send-out procedure added	August 26, 2006	Dr. T. Mazzulli
Slit Skin send-out procedure added	August 26, 2006	Dr. T. Mazzulli
Removed Malaria coverage section	March 14, 2007	Dr. T. Mazzulli
Added Parasite reporting procedure	March 14, 2007	Dr. T. Mazzulli
Modified printing PHOL requisition (from LIS)	March 14, 2007	Dr. T. Mazzulli
Removed UHN printing of PHOL requisition	March 14, 2007	Dr. T. Mazzulli
Added - Printing PHOL Requisition from Softlab for Serology Tests	March 14, 2007	Dr. T. Mazzulli
Revised Histoplasma Ag, Blastomysis Ag and Slit Skin Send out	March 14, 2007	Dr. T. Mazzulli
Annual Review	March 14, 2007	Dr. T. Mazzulli
Send-out Bench duties #24, revised – remove ureaplasma, mycoplasma	May 01, 2007	Dr. T. Mazzulli
New PHOL forms implemented	November 15, 2007	Dr. T. Mazzulli
Annual Review	April 15 2008	Dr. T. Mazzulli
Slit Skin for Leprosy specimen collection link added	January 25, 2009	Dr. T. Mazzulli
Annual Review	April 15 2009	Dr. T. Mazzulli
Annual Review	June 10, 2010	Dr. T. Mazzulli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div></div> <div>Department of Microbiology</div> <div>Quality Manual</div>	Policy # MI_SO	Page 56 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	


Page Number / Item	Date of Revision	Signature of Approval
Annual Review	June 10, 2011	Dr. T. Mazzulli
<a href="#">Send-out Shipping Information Table</a> added	June 05, 2012	Dr. T. Mazzulli
Annual Review	June 05, 2012	Dr. T. Mazzulli
Annual Review	June 05, 2013	Dr. T. Mazzulli
Updated Headers/footers, UHN log, Miravista req Review of Requisitions printed and procedures	August 23, 2014	Dr. T. Mazzulli
Added CJD procedure link	August 23, 2014	Dr. T. Mazzulli
Annual Review	August 23, 2014	Dr. T. Mazzulli
General PHOL Delivery Information	January 20, 2015	Dr. T. Mazzulli
Addition of MERSCoV and BACTERIOLOGY MOLECULAR DETECTION SENDOUTS	February 17, 2015	Dr. T. Mazzulli
Updated after hours/stat send to PHOL, merged in General PHOL inquiries with this. Updated link to sick kids fungus pct	April 1, 2015	Dr. T. Mazzulli
Added KJV courier protocol Added Courier table Moved Serology Workflow to Workflow manual	May 1, 2015	Dr. T. Mazzulli
Added FEDEX courier protocol Annual Review	May 5 <sup>th</sup> , 2015	Dr. T. Mazzulli
Updated Molecular bacte test requisition with July 14 <sup>th</sup> version.	August 1, 2015	Dr. T. Mazzulli
Added comment when sending out for metagenomic assay: "This specimen has been sent to Public Health Ontario Laboratories (PHOL) for metagenomic analysis (a research use only assay). Results will be communicated verbally directly from PHOL to the requesting physician, <b>Dr. [REDACTED]</b> . If you have any questions or would like to receive a verbal result, please contact Dr. Samir Patel from PHL at 647-792-3172."	October 1, 2015	Dr. T. Mazzulli
Added MERS under VIROLOGY MOLECULAR DETECTION SEND OUTS section Created instructions for Measles PCR send out.	October 20, 2015	Dr. T. Mazzulli
Added link to PHL sendout for Zika Virus	January 27, 2016	Dr. T. Mazzulli
<i>Streptococcus pneumoniae</i> send out instructions added.	February 25, 2016	Dr. T. Mazzulli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is *Uncontrolled When Printed*.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_SO	Page 57 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	



Page Number / Item	Date of Revision	Signature of Approval
NICU HIV Testing moved to serology section. HIV testing procedure changed from sick kids to PHOL send out.	March 24, 2016	Dr. T. Mazzulli
Zika Virus changes: <ul style="list-style-type: none"> <li>Changed from a serum &amp; plasma to 2 serum samples</li> <li>Updated orderable specific for Cerner and Citrix</li> </ul> Added for required information: "The required information is mandatory upon ordering in EPR/CERNER and can be found in the LIS order comment when accessioning specimen. "	April 5, 2016	Dr. T. Mazzulli
Annual Review Updated MSH logo in header Addition of Referral laboratory test dictionary section Removed Resulting report section to Refer-out Test Results Reporting Procedure QPCMI17001.	May 5, 2016	Dr. T. Mazzulli
Under Bacteria Molecular testing section of BACTERIOLOGY MOLECULAR DETECTION SENDOUTS added: "For approved Molecular Microbiology tests performed at Sick Kids hospital, complete the <a href="#">Sick Kids Molecular Microbiology PCR order form</a> . Pack and send sample through <a href="#">KJV Courier</a> . "	July 14, 2016	Dr. T. Mazzulli
Updated NICU HIV procedure to include updated requisition format and step to highlight both tests requested on requisition.	July 22, 2016	Dr. T. Mazzulli
ASOT serology to HSC added	July 24, 2016	Dr. T. Mazzulli
Updated NICU HIV procedure to: <ul style="list-style-type: none"> <li>Process STAT</li> <li>Always order HIV PCR</li> <li>Notify if there was problem ordering</li> <li>Use print scheduler for requisition</li> <li>Show specimen to senior before sending</li> <li>Package separately: HIV PCR attn. Jocelyn</li> <li>Send by Courier if misses routine PHL pick up in am</li> <li>Instructions to manage samples if serology/CD4/1 or 2 samples are received</li> </ul>	December 9, 2016	Dr. T. Mazzulli

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is **Uncontrolled When Printed**.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>University Health Network Sciences General Sciences Women Sciences Children Sciences Mental Health</div><div>Mount Sinai Hospital Hôpital du Mont-Sinai Hospital de Mont-Sinai Hospital de Mont-Sinai</div></div> <div>Department of Microbiology</div>	Policy # MI_SO	Page 58 of 59
	Quality Manual	
Section: Bacteriology Procedures	Subject Title: <b>Send-out Information Manual</b>	

Page Number / Item	Date of Revision	Signature of Approval
Annual Review	May 05, 2017	
Annual Review Updated NICU HIV procedure to: -Use HIV form not neonatal HIV PCR form -updated email to fax requisition for PHOL notification. Inserted links for PHOL after hours specimen dropoff	March 05, 2018	Dr. T. Mazzulli
Page 34-47 Appendix II Virology PHOL Pending List added	May 31, 2018	Dr. T. Mazzulli
Annual Review pg 3 added Mayo Clinic to Referral table pg 5-6 updated send-out procedure for QFTB ; added send-out procedure for H. pylori	October 08, 2019	Dr. T. Mazzulli
Pg 20 all CT+ rectal swabs to be sent to NML for LGV	October 30, 2019	Dr. T. Mazzulli
Addition of 16s notification form	August 20, 2020	Dr.T.Mazzulli
Annual Review Updated Measles/Mumps/MERS: <ul style="list-style-type: none"> <li>Removed notification to microbiologist</li> <li>Updated process to notify IPAC and send samples to PHOL right away (do not wait until morning)</li> </ul>	August 31, 2020	Dr.T.Mazzulli
Updated QUANTIFERON TB (QFT) procedure from SickKids requisition	September 25, 2020	Dr.T.Mazzulli

**Full document review included in all updates. Bi-annual review conducted when no revision had been made within 2 years.**


Page Number / Item	Date of Revision	Edited by:
Updated Baby HIV PCR with tube type, contact info, and test collect/receipt info	Jan 13, 2021	Dorna Zareianjahromi
Minor formatting change	April 11, 2021	Jessica Bourke
Added section on Simian B virus antibody testing	May 6, 2021	Wayne Chiu
Clarified instructions TB quantiferon from Sinai outpt	July 9, 2021	Wayne Chiu
Added PHOL Duty officer contact information to measles section and MERS section, also PHOL section	Oct 15, 2021	Wayne Chiu
Added section for strongyloides serology, PHOL update	Feb 10, 2022	Wayne Chiu
Added approval needed for non-donor/non-recipient blood and CSF for WNV PCR testing to PHL	Mar 11, 2022	Qin Liu
Added info regarding zika molecular sendout	Mar 16, 2022	Wayne Chiu

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

<div><div>Department of Microbiology</div></div> <div>Quality Manual</div>	Policy # MI_SO	Page 59 of 59
	Version: 4.2 CURRENT	
Section: Bacteriology Procedures	Subject Title: Send-out Information Manual	

Page Number / Item	Date of Revision	Edited by:
Added section “checking orders for IGM PHL sendout”	March 26, 2022	Wayne Chiu
Updated sendout frequency for measles mumps Add link for monkeypox	June 23, 2022	Wayne Chiu
Minor formatting change	November 24, 2022	Wayne Chiu
Minor formatting change	March 24, 2023	Oliver Li
Added the followings to Baby NICT HIV testing <ul style="list-style-type: none"> <li>Send <b>EDTA Whole Blood</b></li> <li><b><u>Do NOT centrifuge</u></b> the blood</li> </ul>	April 03, 2023	Qin Liu
Addition: APPENDIX III: UHN MANAGED MICROBIOLOGY SEND OUT TESTS FOR UHN PATIENTS Updated Measles/Mumps STAT procedure. For bacteriology send-outs added: <i><b>Note:</b> When possible (when sufficient quantity is received), especially when only referred out test is requested, a portion of sample should be maintained in house for tissue, sterile sites etc</i>	May 4, 2023	Jessica Bourke
Updated send-out for MERS-CoV or Avian Flu (pg 23): Added notes for Avian Flu PCR send-out	June 21	Qin Liu
Added instructions for HIV Viral Load MSH patients only in section of VIROLOGY MOLECULAR DETECTION SEND OUTS	September 19, 2023	Oliver Li
Added hyperlink from Sickkids for QuantiFERON TB testing	August 1, 2024	Oliver Li
Updated the instructions for sample collection and send out to UHN Core Lab for QuantiFERON® TB (QFT) Gold Plus testing	November 12, 2024	Oliver Li
Updated bulletin point 3 and 4 in the section of HIV VIRAL LOAD FOR MSH patients only; Added Greiner tube aliquoting instructions provided by PHOL (page 22-23).	January 24, 2025	Oliver Li

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\