

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # M\MYC\10\03\v01	Page 1 of 2
Section: Mycology Bench Manual	Subject Title: Appendix III - To Determine Cycloheximide Resistance of an Isolate	
Issued by: LABORATORY MANAGER	Original Date: March 22, 2001	
Approved by: Laboratory Director	Revision Date:	

APPENDIX III

DETERMINING CYCLOHEXIMIDE RESISTANCE OF AN ISOLATE

I. Purpose

To rule out or help confirm the presence of a possible dimorphic fungus or dermatophyte. Normally used in the identification of white moulds.

The following pathogenic fungi are resistant to Cycloheximide:

Blastomyces dermatitidis
Histoplasma capsulatum
Coccidioides immitis
Sporothrix schenckii
Paracoccidioides brasiliensis
Trichophyton sp.
Microsporum sp.
Epidermophyton floccosum

Cycloheximide inhibition rules out the above fungi. Resistance to cycloheximide **may** indicate one of the above pathogens.

II. Procedure

1. Subculture the isolate to Potato Dextrose Agar (PDA) and Mycosel Agar (MYC) and incubate at 28⁰C (or RT). MYC contains cycloheximide.
2. Observe periodically for 7-10 days or until good growth on one or both media.

III. Interpretation

- i) PDA+/ MYC + : Fungus **MAY** be one of the above listed pathogens.
- ii) PDA +/- MYC - : Fungus is **NOT** one of the ABOVE listed pathogens.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\MYC\10\03\v01	Page 2 of 2
Mycology Bench Manual		

iii) PDA -/ MYC - : Test invalid, repeat.

NB: *Penicillium marneffei* (dimorphic fungus) is inhibited by cycloheximide.

IV. Reference

M.R. McGinnis. Laboratory Handbook of Medical Mycology 1980, Academic Press.