

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy #MI\TECH\42\01\v01</b>	Page 1 of 2
Section: <b>Technical Manual</b>	Subject Title: <b>Acid Fast Stain for Mycobacteria (Kinyoun)</b>	
Issued by: <b>LABORATORY MANAGER</b>	Original Date: July 31, 2000	
Approved by: Laboratory Director	Revision Date: February 15, 2002	

## ACID FAST STAIN FOR MYCOBACTERIA (KINYOUN)

### Principle

To stain Mycobacteria present in specimens and cultures.

Mycobacteria are different to stain with common aniline dyes. However, they will stain with basic fuchsin. Once stained, they retain the dye despite treatment with mineral acids i.e. HCl H<sub>2</sub>SO<sub>4</sub>. This property of acid fastness may be due to a lipid fraction called mycolic acid. Mycobacteria also exhibit degrees of resistance to decolourization with alcohol.

### Materials

Kinyoun Carbol fuchsin  
3% HCl in 95% ethanol  
Brilliant green

### Procedure

1. Prepare smear over an area of 2-3 sq. cm.
2. Heat fix smear on heating block (56<sup>0</sup>C/1 hr). Then hold to incinerator for 10 secs.
3. Place slide on stain rack and allow to cool. Flood with Kinyoun Carbol fuchsin for 5 min.
4. Rinse off stain with water.
5. Decolourize with 3% acid alcohol for 3 mins.
6. Rinse with water.
7. Repeat decolourization for 1-2 mins. or until no red appears.
8. Rinse with water.
9. Flood slide 3-4 mins. with Brilliant green.

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy #MI\TECH\42\01\v01</b>	Page 2 of 2
<b>Technical Manual</b>		

10. Rinse with water.

11. Air dry. DO NOT BLOT.

### **Microscopy**

Place a drop of oil between the specimen and coverslip and oil again on top. Smears are examined with oil immersion lens. The coverslip prevents cross contamination. An average of 15 mins. is spent on each slide. The total area of the specimen must be examined.

### **References**

1. Baker, F.J., Breach, M.R. 1980. Medical Microbiological Technique, p. 15