

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI04\03\03	Page 1 of 5
Section: Antimicrobial Susceptibility Testing Manual	Subject Title: Appendix III - Double Disk Test for ESBL	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 2005	

APPENDIX III - DOUBLE DISK TEST FOR ESBL

I. Introduction

Class A or Bush Group 1 extended spectrum beta-lactamases (ESBLs) are inhibited by clavulanic acid. This may be detected by testing the suspected organism to a 3rd generation cephalosporin alone and in combination with clavulanic acid. If the combination results in an expanded zone of inhibition compared to that of the 3rd generation cephalosporin alone, it is indicative of the presence of an ESBL.

II. Materials

Mueller-Hinton (MH) agar (150) mm
 20/10 mg amoxicillin-clavulanate disc
 30 mg ceftazidime disc
 30 mg ceftriaxone or cefotaxime disc
 30 mg aztreonam disc
 10 mg cefpodoxime disc (optional)
 30 mg cefoxitin disc
 Quality control strain: *E. coli* ATCC 51446


III. Procedure

1. Prepare a bacterial suspension of the organism to be tested that has a turbidity equivalent to that of a 0.5 McFarland standard.
2. Inoculate a Mueller-Hinton agar plate with this suspension in accordance with NCCLS M100-S10 (M2) guidelines for disc diffusion testing.
3. Place the amoxicillin-clavulanic acid disc towards the centre of the plate.
4. Carefully measure 15 mm out from the edge of that disc at 90° angles marking the plate.
5. Place a ceftazidime disk on the plate so that its inner edge is 15 mm (the mark) from the amoxicillin-clavulanic acid disc (See Figure 1 KB-ESBL Template).
6. Do the same with cefotaxime (or ceftriaxone), aztreonam and cefpodoxime discs so that they are spaced 90° apart and 15mm from the centre disc.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI04\03v03	Page 2 of 5
Antimicrobial Susceptibility Testing Manual		

7. Place a cefoxitin disc in any available space remaining on the plate.
8. Incubate 35°C, in O₂ x 18-24 hours and record the zone diameters for the all cephalosporins as per NCCLS guidelines.

Figure 1. KB-ESBL Template



TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI04\03v03	Page 3 of 5
Antimicrobial Susceptibility Testing Manual		


IV. Interpretation

Note: The following applies to cefpodoxime-nonsusceptible *E. coli*, *Klebsiella* species and *Proteus* species only.

1. Document zone size for all antibiotics.
2. Observe for **potentiation** of the inhibition zone (i.e. **increase** in the inhibition zone) of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid (enter **Yes** or **No** to the “drug” named ESBL Inhibitor in the LIS).
3. If a reduction of zone of inhibition of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid is observed, recheck the identification of the isolate and repeat testing. Notify the charge technologist if result remains unchanged.

Class A ESBL present:

- i) Potentiation of the inhibition zone of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid (see below for examples of different patterns of potentiation that can be seen with organisms that contain Class A ESBLs)
- ii) Susceptibility to ceftazidime
- iii) Susceptibility or resistance to any one of ceftazidime, ceftriaxone or aztreonam



Class A and Class C ESBL present:

- i) Potentiation of the inhibition zone of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid
- ii) Resistant or Intermediate to ceftazidime.
- iii) Susceptibility or resistance to any one of ceftazidime, ceftriaxone or aztreonam

Class C-ESBL present:

- i) No potentiation with clavulanic acid
- ii) Resistance or Intermediate to ceftazidime
- iii) Resistance to any one of ceftazidime, ceftriaxone or aztreonam.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI04\03v03	Page 4 of 5
Antimicrobial Susceptibility Testing Manual		

ESBL not Class A or Class C present:

- i) No potentiation with clavulanic acid
- ii) Susceptibility to cefoxitin
- iii) Resistance to any one of ceftazidime, ceftriaxone or aztreonam

ESBL absent:

- i) No potentiation with clavulanic acid
- ii) Susceptibility or resistance to cefoxitin
- iii) Susceptibility to all of ceftazidime, ceftriaxone or aztreonam

V. Reporting

Reporting Comment	Potentiation of the inhibition zone of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid (enter Y or N to the “drug” named ESBL Inhibitor in the LIS)	Cefoxitin	Ceftazidime, ceftriaxone or aztreonam
The susceptibility pattern suggests that this organism contains a class A extended spectrum beta-lactamase (ESBL).	Yes	S	S/R
The susceptibility pattern suggests that this organism contains class A and C extended spectrum beta-lactamases (ESBL).	Yes	I/R	S/R
The susceptibility pattern suggests that this organism contains a class C extended spectrum beta-lactamase (ESBL).	No	I/R	R
The susceptibility pattern suggests that this organism contains an extended spectrum beta-lactamase (ESBL) other than class A or C.	No	S	R
Not ESBL – no reporting comment	No	S/R	S

VII. References

1. Cormican MG, Marshall SA, Jones RN. 1996. Detection of extended-spectrum beta-lactamase (ESBL)-producing strains by the Etest ESBL screen. J. Clin. Microbiol. 8:1880-1884.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI04\03v03	Page 5 of 5
Antimicrobial Susceptibility Testing Manual		

2. Jacoby GA, Han P. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. J. Clin. Microbiol. 34:908-911.
3. Jacoby GA. 1994. Genetics of extended-spectrum beta-lactamases. Eur. J. Clin. Microbiol. Infect. Dis. Suppl 1,13:2-11.
4. Livermore DM. 1995. Beta-lactamases in laboratory and clinical resistance. Clin. Microbiol. Rev.8:557-584.
5. Moland ES, Sanders CC, Thomson KS. 1998. Can results obtained with commercially available MicroScan microdilution panels serve as an indicator of beta-lactamase (ESBL)-production among *Escherichia coli* and *Klebsiella* isolates with hidden resistance to expanded spectrum cephalosporins and aztreonam? J. Clin. Microbiol. 36:2575-2579.
6. NCCLS. Performance standards for antimicrobial susceptibility testing; Eighth informational supplement. NCCLS document M100-SB [ISBN 1-26238-337-X]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1998)
7. Sanders CC, Barry AL, Washington JA, Shubert C, Moland ES, Traczewski MM, Knapp C, Mulder R. 1996. Detection of extended-spectrum beta-lactamase-producing members of the family Enterobacteriaceae with the Vitek ESBL test. J. Clin. Microbiol. 34:2997-3001.
8. Swenson JM, Ferraro MJ, Jorgensen JH, Tenover FC. 1996. Can different breakpoints be used to detect extended-spectrum beta-lactamase (ESBL)-producing organisms? Abstr D-22, ICAAC meeting, New Orleans, Louisiana, Sept 15 -18, 1996.
9. Vercauteren E, Descheemaeker P, Ieven M, Sanders CC, Goossens H. 1997. Comparison of screening methods for detection of extended-spectrum beta-lactamases and their prevalence among blood isolates of *Escherichia coli* and *Klebsiella* spp. in a Belgian teaching hospital. J. Clin. Microbiol. 35:2191-2197.