

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI\STER\06\v01	Page 1 of 1
Section: Sterility Testing Manual	Subject Title: Dog (Cardiovascular Lab) Specimens	
Issued by: LABORATORY MANAGER	Original Date: July 17, 2001	
Approved by: Laboratory Director	Revision Date: February 15, 2002	

DOG (CARDIOVASCULAR LAB) SPECIMENS

I. Introduction

These specimens are collected from the research laboratory. Dr. Wilson is the contact person (ext. 4795).

II. Specimen Collection and Transport

Specimens are sent in a clean sterile container or in their original container.

III. Processing of Specimens

iii) Direct examination: Gram stain

iv) Culture:

Media	Incubation
Blood Agar (BA)	O ₂ , 35 ⁰ C x 48 hours
MacConkey Agar (MAC)	O ₂ , 35 ⁰ C x 48 hours
Fastidious Anaerobic Broth (THIO)	AnO ₂ , 35 ⁰ C x 7 days

IV. Isolation and Identification

All isolates are to be identified. Gram stain and subculture all turbid FAB.

V. Sensitivity Testing

Not required.

VI. Reporting

Telephone all positive reports to ward / physician.

Negative Report: "No Growth"
Positive Report: Report all isolates