

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/VIR/16/08/v01	Page 1 of 3
Section: Virology Manual	Subject Title: Appendix VIII Media	
Issued by: LABORATORY MANAGER	Original Date: March 14, 2001	
Approved by: Laboratory Director	Revision Date:	

Appendix VIII

MEDIA

GROWTH MEDIUM (Not in routine use)

500 mL Eagles' MEM with Hanks' salt/without glutamine
 10 mL Vitamins
 5 mL L-glutamine (200 mM)
 50 mL Inactivated fetal calf serum

Store at 4°C. Stable for 2 weeks.

MAINTENANCE MEDIUM

500 mL Eagles' MEM with Hanks' salt/without glutamine (stored at 4°C)
 5 mL L-glutamine (200 mM, stored frozen at -20°C in clean room freezer)
 10 mL Inactivated fetal calf serum (stored frozen at -20°C in clean room freezer)
 5 mL Fungizone (250 µg/mL, stored frozen at -20°C in clean room freezer)
 *5 mL Gentamicin (1 mg/mL, stored frozen at -20°C in clean room freezer)
 **5 mL Vancomycin (10 mg/mL, stored frozen at -20°C in clean room freezer)

Store at 4°C. Stable for 2 weeks.

- * **Gentamicin Solution (1 mg/mL):**
 Dilute entire contents (10 mL) of 10 mg/mL Gentamicin Sulphate vial into 90 mL of d. H₂O to achieve 1 mg/mL. Filter sterilize, then dispense into 5 mL aliquots (stored frozen at -20°C in clean room freezer).

- ** **Vancomycin Solution (10 mg/mL):**
 Add contents of 1 gram vial of Vancomycin to 100 mL of distilled H₂O. Filter sterilize, then dispense into 5 mL aliquots (stored frozen at -20°C in clean room freezer).

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/VIR/16/08/v01	Page 2 of 3
Virology Manual		

Media Quality Control:

- a. A 2 mL aliquot of maintenance media is to be placed into each of CMK, HFF, HEp-2 and RD tubes for a check on sterility/toxicity before the media is used.
- b. Maintenance Media is to be registered into the LIS (micqc) using the date made as the lot number. The individual components and their lot numbers are recorded at this time.
- c. QC results are all entered into the LIS when reading tube culture QC on Mondays, Wednesdays and Fridays.

Failed Media QC:

- a. Any cell lines showing microbial contamination or toxicity is considered abnormal.
- b. Do not release media lot for use pending resolution of the QC failure.
- c. Inform charge/senior technologist to discuss further actions including the following.
- d. Record in Incident Report (if necessary).
- e. Repeat QC/re-make media to evaluate if it is the cell media, the supplement(s) or other materials causing the problem.

CELL WASHING MEDIUM (Purchased)

Hanks' Balanced Salt solution without CaCL₂, MgCL₂ and MgSO₄.7H₂O

CELL FREEZING MEDIUM (SIGMA C6164) (Purchased)

MEM supplemented with a mixture of fetal bovine serum and calf serum containing DMSO.

VIRAL FREEZING MEDIUM (SIGMA C6039) (Purchased)

MEM supplemented with a mixture of fetal bovine serum and calf serum containing 10% glycerol.

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/VIR/16/08/v01	Page 3 of 3
Virology Manual		

TRANSPORT MEDIA:

BARTEL'S VIRAL TRANSPORT MEDIUM (B1029-36A) (Purchased)

Eagles' MEM with nonessential amino-acid and L-glutamine in HBSS.

2% FBS

10 µg/mL Gentamicin

50 µg/mL Streptomycin

50 µg/mL Penicillin

4 µg/mL Amphotericin B

15 mM HEPES

9.5 mM Sodium bicarbonate

CHLAMYDIA TRANSPORT MEDIUM (Not in routine use)

500 mL Eagles' MEM with Hank's salt/without glutamine.

50 mL Fetal Calf Serum (normal)

5 mL 3M Glucose

5 mL Gentamicin (1000 µg/mL)

5 mL Vancomycin (1000 µg/mL)

5 mL Fungizone (250 µg/mL)

Combine all components aseptically. Dispense 1.5 mL transport medium into sterile screw-cap vials containing 3 glass beads. Store at -20°C.