

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy #MI\VIT\15\01\v01</b>	Page 1 of 2
Section: <b>Vitek Manual</b>	Subject Title: <b>Quality Control</b>	
Issued by: <b>LABORATORY MANAGER</b>	Original Date: October 4, 1995	
Approved by: Laboratory Director	Revision Date: May 30, 2001	

## QUALITY CONTROL

QC all new lots of **identification cards** once when received in the lab. No further QC is required if all results are within limits. Perform QC's for 5 consecutive days if out-of-control results are observed (see QA technologist).

QC all new lots of **susceptibility cards** once when received and once weekly when lot is in use. Perform daily QCs when any out-of-control results are observed (see QA technologist).

### Quality Control Set Up For Vitek Bench

#### **DAY 1**

Subculture purity plates from the stock slants to BA for the next days QC requirements.

#### **DAY 2**

1. Remove fresh subculture plates from incubator and check for purity.
2. At Main Menu:
  - a. click "**VITEK**"
  - b. click "**QC**"
  - c. click "**LOADLIST**"
  - d. select card by scrolling screen by screen; hi-light appropriate card type and lot # by holding left button down
  - e. click "**PRINT**"
  - f. click "**SELECTED ITEM**"
3. Proceed to set up cards according to the load list.

#### **DAY 3**

1. Remove purity plates from incubator and check purity.
2. Upon completion of the test, results will be transferred automatically to the permanent QC data base.

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI\VIT\15\01\v01</b>	Page 2 of 2
<b>Vitek Manual</b>		

3. Go to QC file to check for errors:
  - a. click on "**VITEK**"
  - b. click on "**QC**"
  - c. click on "**VIEW RESULTS**"
  - d. enter: **test type**  
**lot #**  
**date range**  
**deviations only**  
click "**OK**"
  - e. if exceptions are present, click "**YES**" to load results
  - f. click "**PRINT**"
  - g. click "**◆**" if there is more than one exception for this card type.
  - h. click "**FILE**"
  - i. click "**QUIT**"
  - j. click "**SYSTEM**"
  - k. move cursor to the problem QC card and click
  - l. click "**LAB REPORT**" button, hold left button and move cursor to "**RAW DATA REPORT**"
  - m. click "**PRINT**"
  - n. See Pauline/Glen for corrective action.
  - o. File load list, exception list, QC exception reports and raw data printouts in QC binder.
  - p. Repeat steps b to o for all card types.

For QC Results that appear in the "REVIEW" file eg. unidentified organism, delete result as in p. 8 and repeat test.