

|  |                                   |             |
|--|-----------------------------------|-------------|
| TML\MSH Microbiology Department<br>Policy & Procedure Manual | <b>Policy # MI\QC\10\v01</b>      | Page 1 of 1 |
| Section: <b>Quality Control Manual</b>                       | Subject Title: <b>Serology QC</b> |             |
| Issued by: <b>LABORATORY MANAGER</b>                         | Original Date: April 11, 2001     |             |
| Approved by: Laboratory Director                             | Revision Date: February 15, 2002  |             |

### **SEROLOGY QC**

All new reagents lots will be tested using samples and/or commercially available reference material (where available) before being placed in service.

All Serology tests require QC run in parallel with each clinical sample test run. Refer to the appropriate tests for procedure and method.

Record all QC results into the LIS.

Inform the charge technologist or senior technologist of all out-of-range results for the appropriate corrective action.

Record all corrective action into the LIS.