

TML\MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI\VIR\17\08\v02</b>	Page 1 of 2
Section: <b>Serology Manual</b>	Subject Title: <b>Appendix VIII Autoverification for the Probetec</b>	
<b>Issued by:</b> LABORATORY MANAGER	Original Date: June 25, 2002	
Approved by: Laboratory Director	Revision Date:	

### APPENDIX VIII

#### AUTOVERIFICATION FOR THE PROBETEC

**POLICY:**

The Probetec performs the following tests for the Virology department:

1. Chlamydia trachomatis detection
2. Neisseria gonorrhoeae detection.

**Autoverification(Autofinalizing) will occur as follows:**

TEST NAME	SOFT TEST CODE	RESULT	TRANSLATION	ACTION	FINAL RESULT
Chlamydia Detection	CHLA	POSITIVE	?POSITIVE	No Autoverification	Canned Message: }CHLP
Chlamydia Detection	CHLA	NEGATIVE		Autoverifies/Autofinalizes	Canned Message: }CHLN
Neisseria gonorrhoeae	NGON	POSITIVE	?POSITIVE	No Autoverification	Canned Message: }GCPO
Neisseria gonorrhoeae	NGON	NEGATIVE	?NEGATIVE	Autoverifies/Autofinalizes	Canned Message: }GCNE

TML\MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI\VIR\17\08\v02</b>	Page 2 of 2
Section: <b>Serology Manual</b>		

**Canned Messages:**

}CHLP - POSITIVE for Chlamydia trachomatis by Nucleic Acid Amplification.

}CHLN - Negative for Chlamydia trachomatis by Nucleic Acid Amplification.

}GCPO - POSITIVE for Neisseria gonorrhoeae by Nucleic Acid Amplification.

}GCNE - Negative for Neisseria gonorrhoeae by Nucleic Acid Amplification.